

Sep 19, 2022

Midbrain organoid differentiation in spinner flasks.

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzbnr4vx1/v1>

gustavo.parfitt Parfitt¹

¹Icahn School of Medicine

Ahfeldt Lab

Organoid and Assembloid



gustavo.parfitt Parfitt

University College of London

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzbnr4vx1/v1>

Protocol Citation: gustavo.parfitt Parfitt 2022. Midbrain organoid differentiation in spinner flasks.. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vzbnr4vx1/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 16, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 70150

Keywords: ASAPCRN, midbrain organoid differentiation in spinner flask, using spinner flask, spinner flask, midbrain organoid differentiation, midbrain differentiation protocol

Funders Acknowledgements:

Lily Sarrafha

Tim Ahfeldt

Joel Blanchard

Abstract

Midbrain differentiation protocol using spinner flasks.



Materials

Media for iPSC

StemFlex medium

Media for spheres making

StemFlex medium + 10 μ M Y-27632 + 1:100 Pen/Strep

Media composition for diferenciation

D0-1 - DMEM/F-12+Glutamax + 13 B-27 minus vitamin A + 13 N-2 + 100 nM LDN193189 + 10 μ M SB431542

D2-3 - DMEM/F-12+Glutamax + 13 B-27 minus vitamin A + 13 N-2 + 100 nM LDN193189 + 10 μ M SB431542+ 2 μ M Purmorphamine + 1 μ M SAG

D4-7 - DMEM/F-12+Glutamax + 13 B-27 minus vitamin A + 13 N-2 + 100 nM LDN193189 + 10 μ M SB431542 + 2 μ M Purmorphamine + 1 μ M SAG + 3 μ M CHIR99021

D8-11 - DMEM/F-12+Glutamax + 13 B-27 minus vitamin A + 13 N-2 + 100 nM LDN193189 + 3 μ M CHIR99021

D12-21 (Terminal Media) DMEM/F-12 + Glutamax + 13 B-27 minus vitamin A + 13 N-2 + 20 ng/mL BDNF + 20 ng/mL GDNF + 0.2 mM Ascorbic Acid + 10 μ M DAPT + 0.1 μ M dcAMP

D22-34 (Terminal Media) DMEM/F-12+Glutamax + 13 B-27 minus vitamin A + 13 N-2 + 10 ng/mL BDNF + 10 ng/mL GDNF + 0.2 μ M Ascorbic Acid + 10 μ M DAPT + 0.1 mM dcAMP

D35+ (Long-Term Media) DMEM/F-12+Glutamax + 13 B-27 minus vitamin A + 13 N-2 + 10 ng/mL BDNF + 10 ng/mL GDNF + 0.2 μ M Ascorbic Acid



iPSCs aggregation in spinner flasks

8m

- 1 Passage iPSC lines to four 10-cm dishes to generate 40×10^6 cells.
- 2 Prepare 120 mL of StemFlex per flask with rock inhibitor (refer to material section) and add 115 to the spinner flask.
- 3 Place the spinner flask on the stir plate and set to 65 rpm.
- 4 Wait for the cells to reach 80% confluence and seed iPSC 40×10^6 cells into a spinner flask.
- 5 For the passaging add  3 mL of Accutase for  00:05:00 .
- 6 Add  5 mL of StemFlex , gently mix and transfer to a 15 ml conical tube.
- 7  300 x g, 25°C, 00:03:00
- 8 Resuspend the cells in  1 mL of StemFlex with rock inhibitor and Count the cells.
- 9 Put the cells in the correct cell concentration and add to  5 mL of StemFlex with rock inhibitor.
- 10 From the main opening of the spinner flask add the cells and plate back on the stir plate.
- 11 Every other day change  60 mL of media using StemFlex until the spheres reach 300 μ M.

5m

3m

Midbrain Differentiation

1w 5d



- 12 On the starting day D0, filter the spheres in a 300 μ M to 500 μ M range.
- 13 Aspirate all the media in the flask and replace by D0-D1 media (refer to materials). Add the filtered spheres and place the spinner flask back on the stir plate.
- 14 D1 change half the media with D0-1 medium.
- 15 D2 change half the media with D2-3 medium.
- 16 D3 Change half the media with D2-3 medium.
- 17 D4 change half the media with D4-7 medium.
- 18 D5 change half the media with D4-7 medium.
- 19 D6 change half the media with D4-7 medium.
- 20 D7 change half the media with D4-7 medium.
- 21 D8 change all the media with D8-11 medium.
- 22 D9 change half the media with D8-11 medium.
- 23 D10 change half the media with D8-11 medium.
- 24 D11 change half the media with D8-11 medium.



- 25 On D12 change to D12 terminal media for final differentiation. From this point on the organoids can be transfer to low attachment plates and place on shaker.
- 26 On D35 change to long-term maintaining media.
- 27 Protocol Media schedule
- <!--td {border: 1px solid #cccccc;}br {mso-data-placement:same-cell;}-->

A	B
Day x--1	StemFlex
Days 0-1	DMEM/F12+Glutamax + B27 + N2 + SB + LDN
Days 2-3	DMEM/F12+Glutamax + B27 + N2 + SB + LDN + Purmorphamine + SAG
Days 4-7	DMEM/F12+Glutamax + B27 + N2 + SB + LDN + Purmorphamine + SAG + CHIR
Days 8-11	DMEM/F12+Glutamax + B27 + N2 + LDN + CHIR
Days 12-21	DMEM/F12+Glutamax + B27 + N2 + BDNF + GDNF + DAPT+ Ascorbic Acid + dCAMP
Days 22-35	DMEM/F12+Glutamax + B27 + N2 + 0.5(BDNF + GDNF) + DAPT + Ascorbic Acid + dCAMP
Days 35-end	DMEM/F12+Glutamax + B27 + N2 + 0.5(BDNF + GDNF) + Ascorbic Acid