

Nov 17, 2025

Midbrain organoid dissociation for Single-cell transcriptome library preparation and sequencing

 In 2 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.261gek7odg47/v1>

Hariam Raji¹, michela.deleidi¹

¹Institute Imagine Paris (24 Bd du Montparnasse, 75015 Paris, France)

Hariam Raji: Mechanisms and Therapy of Genetic Brain Diseases

michela.deleidi : Mechanisms and Therapy of Genetic Brain Diseases

Team Deleidi



Hariam Raji

Institut Imagine Paris - Inserm U1163, Aligning Science Acro...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261gek7odg47/v1>



Protocol Citation: Hariam Raji, michela.deleidi 2025. Midbrain organoid dissociation for Single-cell transcriptome library preparation and sequencing. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261gek7odg47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 05, 2025

Last Modified: November 17, 2025

Protocol Integer ID: 231606

Keywords: enzymatic dissociation of midbrain organoid, midbrain organoid dissociation, midbrain organoid, cell transcriptome library preparation, 10x genomics scrna, preparation of 10x genomics scrna, organoid, enzymatic dissociation, illumina novaseq platform, cell suspension, midbrain, seq library

Abstract

This protocol describes the enzymatic dissociation of midbrain organoids into single-cell suspensions, preparation of 10x Genomics scRNA-seq libraries, and high-depth sequencing on an Illumina NovaSeq platform.

Guidelines

iPSC-derived midbrain organoids were dissociated into single-cell suspensions using a Worthington Papain Dissociation System kit according to manufacturer's instructions. Organoids were minced and incubated with 2.5 mL papain/DNAase solution on an orbital shaker for 30 min at 37 °C. Organoid suspensions were triturated using scalpels and 1 mL low-attachment pipette tips and then incubated on an orbital shaker for 20 min at 37 °C twice. After enzyme inactivation, the cell suspensions were filtered with a 40-µm cell strainer. The final cell suspension was centrifuged at 300×g for 7 min, and the cell pellet was resuspended in PBS-0.04% BSA at a final concentration of 1000 cells/µL, ensuring that more than 95% of the cells were viable. scRNA-Seq libraries were generated with the 10X Chromium Next GEM Single Cell 3' Reagent Kit v3.1 according to the manufacturer's instructions. Libraries were pooled and subjected to paired-end sequencing on the Illumina NovaSeq X platform with a sequencing depth of 300 million reads per library.

Materials

Worthington Papain Dissociation System kit (Worthington Biochemical), 2.5 mL papain/DNAase solution, scalpel, 1 mL low-attachment pipette tips, 40-µm cell strainer, PBS-0.04% BSA, 10X Chromium Next GEM Single Cell 3' Reagent Kit v3.1 (Cat. number 1000128, 10x Genomics), Illumina NovaSeq X platform (SP Flow Cell, Illumina)

Midbrain organoid dissociation for Single-cell transcriptome library preparation and sequencing

- 1 Dissociate iPSC-derived midbrain organoids into single-cell suspensions using a Worthington Papain Dissociation System kit according to manufacturer's instructions.
- 2 Mince individual organoids using scalpels and incubate with 2.5 mL papain/DNAase solution on an orbital shaker for 30 min at 37 °C.
- 3 After 30 min at 37 °C, triturate organoid suspensions using 1 mL low-attachment pipette tips, then incubate on an orbital shaker for another 20 min at 37 °C.
- 4 After enzyme inactivation using the inhibition solution, filter the cell suspensions with a 40- μ m cell strainer. Centrifuge the final cell suspension at 300 $\times$ g for 7 min, and resuspend the cell pellet in PBS-0.04% BSA at a final concentration of 1000 cells/ μ L, ensuring that more than 95% of the cells are viable.
- 5 scRNA-Seq libraries were generated with the 10X Chromium Next GEM Single Cell 3' Reagent Kit v3.1 (Cat. number 1000128, 10x Genomics, Pleasanton, CA, USA) according to the manufacturer's instructions. Libraries were pooled and subjected to paired-end sequencing on the Illumina NovaSeq X platform (SP Flow Cell, Illumina) with a sequencing depth of 300 million reads per library.