

Nov 21, 2025

Midbrain sample preparation for 10xGenomics Multiome Protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.14egnr17zl5d/v1>

Matthew Kearney^{1,2}, Zhixiang Liao^{1,3}, Idil Tuncali^{1,3}, Nathan Haywood^{1,4}, Monika Sharma^{1,2}, Madison Cline^{1,5}, Geidy E Serrano^{1,5}, Thomas G Beach^{1,5}, Joshua Z Levin^{1,4}, Clemens Scherzer^{1,2,6,7,8}

¹Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA;

²Stephen & Denise Adams Center for Parkinson's Disease Research of Yale School of Medicine, New Haven, CT 06510;

³Brigham and Women's Hospital;

⁴Stanley Center for Psychiatric Research, Broad Institute of Harvard and MIT, Cambridge, MA 02142, USA;

⁵Banner Sun Health Research, Sun City, AZ 85351, USA;

⁶APDA Center for Parkinson Precision Medicine, Yale, New Haven, CT 06510;

⁷Department of Neurology, Yale, New Haven, CT 06510; ⁸Department of Genetics, Yale, New Haven, CT 06510

Scherzer Lab



Beatrice Weykopf

Yale

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egnr17zl5d/v1>

Protocol Citation: Matthew Kearney, Zhixiang Liao, Idil Tuncali, Nathan Haywood, Monika Sharma, Madison Cline, Geidy E Serrano, Thomas G Beach, Joshua Z Levin, Clemens Scherzer 2025. Midbrain sample preparation for 10xGenomics Multiome Protocol. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egnr17zl5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 29, 2025

Last Modified: November 21, 2025

Protocol Integer ID: 231075

Keywords: midbrain sample preparation for 10xgenomics multiome protocol, 10xgenomics multiome protocol, processing of human frozen postmortem midbrain tissue, human frozen postmortem midbrain tissue, midbrain sample preparation, nuclei multiome library generation, rnaseq, isolating single nuclei, single nuclei, nuclei

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000301

Abstract

This protocol describes the processing of human frozen postmortem midbrain tissue for isolating single nuclei. The nuclei are purified and sorted and subjected to 10X single-nuclei multiome library generation and sequencing (RNAseq + ATAC-seq). Four samples can be processed simultaneously with this protocol.

Materials

1. 20X Nuclei Buffer (10x Genomics, 2000207, Supplied in Multiome Kit), -20°C
2. Sigma Protector RNase Inhibitor (Sigma,3335402001), -20°C
3. Digitonin (Thermo Fisher, BN2006), 4°C, incubate at 65 °C before use
4. 1000mM DTT (Sigma, 43816-10ml), 4°C, freshly dilute to 100mM with water
5. 10% MACS BSA Stock Solution (Miltenyi Biotec, 130-091-376), 4°C
6. 1mg/ml 7-ADD Solution (Invitrogen, A1310)
7. RNase-Free Disposable Pellet Pestles (Fisher Scientific, 12-141-368), RT
8. DPBS (Gibco, 14190-144), RT
9. Sodium Chloride (Fisher, S671), RT
10. Magnesium Chloride (Sigma, M9272), RT
11. 1M pH7.4 Trizma Hydrochloride Solution (Sigma, T2196), RT
12. 100% Nonidet P40 Substitute (Sigma, 74385), RT, dilute to 10% with water
13. 10% Tween 20 (Bio-Rad, 1662404), RT
14. 70µm Nylon Mesh (Fisher, 22363548), RT
15. 40µm Flowmi Cell Strainer (Bel-Art, 136800040), RT
16. Wide bore tips (Thermo Scientific, 2079G), RT

	A	B	C
	All Buffer Amounts Sufficient for Processing 4 Samples		
	NP40 Lysis Buffer (NPLB)		
		Stock	Amount
	pH 7.4 Tris-HCl	1M	57µl
	NaCl	1M	57µl
	MgCl ₂	1M	17.1µl
	Nonidet P40 Substitute	10% (Diluted from 100% with water before use)	57µl
	DTT	100mM (Diluted from 1M with water before use)	57µl
	RNase Inhibitor	40U/µl	142.5µl
	Nuclease-free Water		5.31ml
	Total		5.7mL

A	B	C
PBS Buffer (PBSB)		
	Stock	Amount
BSA	10%	1ml
RNase Inhibitor	40U/ μ l	250 μ l
DPBS		8.75ml
Total		10ml
1X Lysis Buffer (LB)		
	Stock	Amount
pH 7.4 Tris-HCl	1M	10 μ l
NaCl	1M	10 μ l
MgCl ₂	1M	3 μ l
Tween-20	10%	10 μ l
Nonidet P40 Substitute	10% (Diluted from 100% with water before use)	10 μ l
Digitonin (Incubate at 65°C before use)		2 μ l
BSA	10%	100 μ l
DTT	100mM (Diluted from 1M with water before use)	10 μ l
Nuclease-free Water		845 μ l
Total		1ml
Lysis Dilution Buffer (DLB)		
	Stock	Amount
pH 7.4 Tris-HCl	1M	10 μ l
NaCl	1M	10 μ l

A	B	C
MgCl ₂	1M	3μl
BSA	10%	100μl
DTT	100mM (Diluted from 1M with water before use)	10μl
Nuclease-free Water		867μl
Total		1ml
0.1X Lysis Buffer (0.1XLB)		
	Stock	Amount
LB		44μl
Rnase Inhibitor	40U/μl	11μl
DLB		385μl
Total		440μl
Wash Buffer (WB)		
	Stock	Amount
pH 7.4 Tris-HCl	1M	22μl
NaCl	1M	22μl
MgCl ₂	1M	6.6μl
BSA	10%	220μl
Tween-20	10%	22μl
DTT	100mM (Diluted from 1M with water before use)	22μl
Rnase Inhibitor	40U/μl	55μl
Nuclease-free Water		1.83ml
Total		2.2ml

	A	B	C
	Diluted Nuclei Buffer (DNB)		
		Stock	Amount
	20X Nuclei Buffer		5µl
	DTT	100mM (Diluted from 1M with water before use)	1µl
	Rnase Inhibitor	40U/µl	2.5µl
	Nuclease-free Water		91.5µl
	Total		100µl

Protocol materials

 7-AAD (7-Aminoactinomycin D) Thermo Fisher Catalog #A1310

Troubleshooting



Cryosectioning

- 1 Frozen Brain Samples were received from the brain bank - Four 50µm sections were collected in a 1.5 ml microcentrifuge tube and shipped on dry ice to the Lab.

Prepare Buffers and Equipment

- 2 All buffers are prepared the day of the experiment. DTT and RNase inhibitor are added to the buffers immediately before using. Prepare & maintain all buffers  On ice . See all buffer recipes in the Materials tab.
- 3 Pre-cool the microcentrifuge to  4 °C

Nuclei Isolation and Sorting

25m

- 4 Add  650 µL **NPLB** to the tube with tissue and homogenize 15x using a pestle  On ice .
- 5 Add  650 µL **NPLB** to the tube and incubate for  00:05:00  On ice , pipette mix a few times during incubation with wide bore tips (use in step 4, 5 and 6).  
- 6 Pass the suspension through a 70 µm Nylon Mesh into a pre-chilled 50 ml Falcon tube  On ice .
- 7 Transfer the filtered lysate to a pre-chilled Lo-Bind 2 ml tube  On ice .
- 8 Centrifuge at  500 rcf, 4°C,  00:05:00 .  
- 9 Remove most of the supernatant, leaving around  50 µL .
- 10 Add  1 mL **PBSB**, DO NOT mix. Incubate for  00:05:00  On ice .  

11 Pipette mix to resuspend the pellet



12 Centrifuge at 500 rcf, 4°C, 00:05:00 .

5m



13 Remove most of the supernatant, leaving around 50 μ L .

14 Resuspend with 1.4 mL **PBSB**.

15 Add 1.4 μ L 7-AAD (7-Aminoactinomycin D) Thermo Fisher Catalog #A1310 to the 1.4 mL sample, mix by inverting.



16 Incubate for 00:05:00 On ice .

5m



17 Pass the suspension through a 40 μ m FlowMi strainer into a pre-chilled 5 ml FACS tube On ice :

17.1 Aspirate 700 μ L of the suspension into a regular 1 mL pipette tip. Add the FlowMi filter to the end of the tip. Press the plunger slowly to dispense into the new tube. Discard tip and filter.

17.2 Repeat step 17.1 one additional time for the remaining volume.

18 Pre-load 1.5 ml microcentrifuge tubes with 15 μ L **PBSB** for nuclei collection. Maintain On ice .

19 Sort nuclei using a 100 μ m nozzle on BD FACSAria or similar. Collect sorted nuclei in chilled, pre-loaded 1.5 ml microcentrifuge tubes.



Note

Gate first by SSC and FSC to remove doublets. Then, **collect 7-AAD positive events - these are nuclei.**

Nuclei Permeabilization

12m

20 Centrifuge sorted nuclei at  500 rcf, 4°C, 00:05:00

5m



21 Remove the supernatant without disrupting the nuclei pellet.

22 Resuspend the pellet in  100 μ L **0.1XLB** and pipette mix 5x.



23 Incubate for  00:02:00  On ice .

2m



24 Add  500 μ L **WB** and pipette mix 5x.



25 Centrifuge at  500 rcf, 4°C, 00:05:00 .

5m



26 Remove the supernatant without disrupting the nuclei pellet.

27 Based on the nuclei count post-sorting & targeted nuclei recovery, resuspend in chilled **DNB**.

Note

See table in Step 1.1 of the 10X Multiome User Guide (Rev E) for help with calculating **DNB** volume.



Note

Based on empirical data, actual nuclei counts are approximately 1/3rd of collected FANS events.

- 28 Proceed immediately to Chromium Next GEM Single Cell Multiome ATAC + Gene Expression User Guide (CG000338 Rev E, page 30). Library preparation and sequencing for both RNAseq + ATAC-seq is performed following strictly 10x Genomics protocol (CG000338 Rev E).