

Jan 13, 2025

🌐 Mild Immunoprecipitation with Low Background

DOI

dx.doi.org/10.17504/protocols.io.eq2lywe9pvx9/v1

Cole Sitron¹, F Ulrich Hartl¹

¹Max Planck Institute of Biochemistry



Cole Sitron

Max Planck Institute of Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lywe9pvx9/v1>

Protocol Citation: Cole Sitron, F Ulrich Hartl 2025. Mild Immunoprecipitation with Low Background. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eq2lywe9pvx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 02, 2024

Last Modified: January 13, 2025

Protocol Integer ID: 103395

Keywords: ASAPCRN, immunoprecipitation of protein, immunoprecipitation protocol, mild immunoprecipitation with low background, mild immunoprecipitation, immunoprecipitation, bait protein, bait protein of interest, antibody complex to immuno, antibody, antibody complex, immuno, protein aggregate, protein complex, protein, using gentle lysis condition, bind tube, gentle lysis condition

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Aligning Science Across Parkinson's

Grant ID: ASAP- 024268

Abstract

There is often a trade off between buffer stringency and background when performing an immunoprecipitation. One can reduce background with a stringent buffer, but this stringency could also cause a loss of interactors. This immunoprecipitation protocol circumvents this dilemma by using gentle lysis conditions to preserve protein complexes, while lowering the background with BSA blocking, Tween, and low-bind tubes. We have found that this protocol works well for immunoprecipitation of proteins from lysates that contain "sticky" protein aggregates. This protocol uses a protein G-antibody complex to immuno-precipitate the bait protein of interest.

Materials

1. 4 80% confluent 6-wells of cells expressing your protein of interest.
2. Dynabeads™ Protein G Immunoprecipitation Kit (Invitrogen Cat. no. 10007D)
 -  Dynabeads™; Protein G Immunoprecipitation Kit **Thermo Fisher Catalog #10007D**
3. An appropriate antibody (Ab) against your protein of interest
4. Low bind 1.5 ml centrifuge tubes (Eppendorf cat. no. 0030108116)
 -  Protein LoBind Tubes, 1.5 mL **Eppendorf Catalog #0030108116**
5. Microcentrifuge
6. Magnetic tube rack
7. Tube rotator
8. Pierce Rapid Gold BCA Protein Assay Kit (Thermo Fisher Scientific A53225)
 -  Pierce™; Rapid Gold BCA Protein Assay Kit **Thermo Fisher Catalog #A53225**
9. PBS pH 7.2 (Thermo Fisher Scientific cat. no. 20012068)
 -  PBS, pH 7.2 **Thermo Fisher Catalog #20012068**
10. Lysis buffer: PBS pH 7.2, 0.1% Triton X-100 (Thermo Fisher Scientific cat. no. 85111)
 -  Triton™; X-100 Surfact-Amps™; Detergent Solution **Thermo Fisher Catalog #85111** , 0.02% Tween-20 (Sigma-Aldrich cat. no. P9416-100ML)
 -  Tween 20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9416-100ML** , 1X cComplete, Mini EDTA-free Protease Inhibitor Cocktail (Roche cat. no. 11873580001)
 -  cComplete™ EDTA-free Protease Inhibitor Cocktail **Merck MilliporeSigma (Sigma-Aldrich) Catalog #11873580001**
- + 1X PhosStop (Sigma-Aldrich cat. no. 04906837001)
 -  Roche PhosSTOP™ **Merck MilliporeSigma (Sigma-Aldrich) Catalog #4906837001** + 0.75 U/ml Benzonase (Max Planck Institute of Biochemistry Core Facility)
11. 2X Blocking Buffer: PBS pH 7.2, 0.1 % Triton X-100, 0.02% Tween-20, 6% BSA (Cell Signaling Technologies cat. no. 9998S)
 -  BSA **Cell Signaling Technology Catalog #9998S**
12. Wash Buffer: PBS pH 7.2, 0.1 % Triton X-100, 0.02% Tween-20, additional 113 mM NaCl (final concentration 250 mM)
13. 5 M NaCl
14. 1 M DTT
15. 4X NuPAGE LDS Sample Buffer (Invitrogen cat. no. NP0007)
 -  NUPAGE LDS sample buffer (4x) **Thermo Fisher Scientific Catalog #NP0007** + 5% β-mercaptoethanol (Sigma Aldrich cat. no. M6250-100ml)
 -  2-mercaptoethanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M6250**
16. TrypLE Express (Gibco cat. no. 12605036)
 -  TrypLE™; Express Enzyme (1X), phenol red **Thermo Fisher Catalog #12605036**

17. 10% FBS (Thermo Fisher Scientific cat. no. 10270106) in PBS

 Fetal Bovine Serum **Gibco - Thermo Fisher Scientific Catalog #10270106**

Troubleshooting

Bead preparation:

30m

- 1 For each IP, prepare two low-bind tubes with  30 μL Protein G beads and  300 μL Ab Binding Buffer (buffer comes from the Dynabeads Protein G IP Kit; one tube will be the no antibody control to assess background).



Note

Remember to always pipette the beads with a cut (or wide-bore) pipette tip to avoid damaging the beads.

- 2 Add the antibody to one tube and an equivalent amount of PBS to the other and rotate all tubes at  10 rpm, Room temperature, 00:30:00 .
- 2.1 The amount of antibody should be appropriate for  200 μg lysate. The exact amount must be determined empirically or gleaned from the manufacturer's suggestions.
- 2.2 Check the Protein G Dynabeads manufacturer datasheets to ensure that the antibody you are using can be bound by Protein G.
- 3 Wash 1X with Ab binding buffer on the magnetic tube rack.
- 4 Block  01:00:00 with 1X Blocking Buffer (diluted in Lysis Buffer) on the tube rotator.
- 5 Wash 1x with Lysis Buffer.

30m



1h



Lysate preparation:

3m

- 6 While blocking beads, remove the medium from the wells and trypsinize with  500 μL TrypLE Express.
- 7 Quench the TrypLE Express with  500 μL 10% FBS and transfer cells into a centrifuge tube.



8 Pellet cells at  1500 x g, 4°C, 00:03:00 .

3m



9 Wash cells with  1 mL PBS and pellet again.



10 Lyse pellets in  150 µL Lysis Buffer by gentle pipetting.



11 Incubate  On ice for  00:20:00 .

20m



11.1 To reduce nonspecific binding of lysate proteins to the beads, the lysate can be centrifuged for  1000 x g, 4°C, 00:05:00 . The supernatant can then be further processed and prepared for immunoprecipitation. While this will reduce background, it can come at the expense of losing some signal if the protein of interest pellets after centrifugation.

5m



12 Quantify total protein concentration after diluting a small sample 1:10 in PBS (BCA Gold).

13 Set up samples  480 µg in  300 µL . To each sample, add:  287 µL Lysis Buffer +  13.56 µL  5 Molarity (M) NaCl.



IP:

1h

14 Add  250 µL lysate ( 200 µg) to each tube of beads.



Note

The remaining amount of sample left can be reserved to run as an input fraction.

15 Rotate at  4 °C for  01:00:00 .

1h



16 Take off flow through and set aside to run as a flow-through fraction.

17 Wash beads 3x with Wash Buffer.



18 Allow the tubes to rotate briefly on the first wash and transfer to a new tube for the last wash.



19 Elute proteins by adding  30 μ L 4X NuPAGE LDS Sample Buffer + 5% β -mercaptoethanol and boiling the sample at  95 °C for  00:05:00 .

5m



20 Optional: dilute the eluate 1:2 with Wash Buffer +  10 millimolar (mM) DTT and boil again.



20.1 This will ensure total reduction of the antibody to reduce its detection by secondary antibodies upon immunoblotting. We use Veriblot (Abcam cat. no. ab131366) as a secondary antibody to more specifically detect the primary antibody used in the immunoblot rather than the denatured antibody used in the immunoprecipitation.