


Oct 24, 2019

Miracle Prep for Plasmid Isolation

DOI

dx.doi.org/10.17504/protocols.io.8p2hvvqe



Blake Flood¹

¹University of Chicago



Blake Flood

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.8p2hvvqe

Protocol Citation: Blake Flood 2019. Miracle Prep for Plasmid Isolation. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.8p2hvvqe>

Manuscript citation:

Pronobis MI, Deutch N, Peifer M (2016) The Miraprep: A Protocol that Uses a Miniprep Kit and Provides Maxiprep Yields. PLoS ONE 11(8): e0160509. <https://doi.org/10.1371/journal.pone.0160509>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 24, 2019

Last Modified: October 24, 2019

Protocol Integer ID: 29146

Keywords: DNA isolation, cloning, plasmid prep



Abstract

Based on:

<https://journals.plos.org/plosone/article?id=10.1371/journal.pone.0160509>

Materials

MATERIALS

- ☒ Sterile deionized H₂O
- ☒ LB Broth
- ☒ Buffer P1 **Qiagen Catalog #19051**
- ☒ Buffer P2 **Qiagen Catalog #19052**
- ☒ Buffer N3 **Qiagen Catalog #19064**
- ☒ 96% EtOH
- ☒ Econospin Mini Spin Column 250ct **Epoch Life Science Catalog #1920-250**
- ☒ Buffer PB **Qiagen Catalog #19066**
- ☒ Buffer PE **Qiagen Catalog #19065**

STEP MATERIALS

- ☒ Buffer P1 **Qiagen Catalog #19051**
- ☒ Buffer P2 **Qiagen Catalog #19052**
- ☒ Buffer N3 **Qiagen Catalog #19064**
- ☒ Buffer PB **Qiagen Catalog #19066**



Protocol materials

⊗ 96% EtOH

⊗ Econospin Mini Spin Column 250ct **Epoch Life Science Catalog #1920-250**

⊗ Buffer PB **Qiagen Catalog #19066**

⊗ Buffer P1 **Qiagen Catalog #19051**

⊗ Sterile deionized H₂O

⊗ LB Broth

⊗ Buffer P2 **Qiagen Catalog #19052**

⊗ Buffer P2 **Qiagen Catalog #19052**

⊗ Buffer PE **Qiagen Catalog #19065**

⊗ Buffer N3 **Qiagen Catalog #19064**

⊗ Buffer PB **Qiagen Catalog #19066**

⊗ Buffer P1 **Qiagen Catalog #19051**

⊗ Buffer N3 **Qiagen Catalog #19064**

⊗ Buffer P1 **Qiagen Catalog #19051**

⊗ Buffer P2 **Qiagen Catalog #19052**

⊗ Buffer N3 **Qiagen Catalog #19064**

⊗ Buffer PB **Qiagen Catalog #19066**



Grow bacteria containing plasmid

- 1 Pick colony from plate. Inoculate 50 mL LB broth per vector.
- 2 Allow to grow overnight at 37 °C .

Isolate Plasmid

- 3 Spin culture 00:10:00 4000 x g 4 °C

Isolate Plasmid

- 4 Remove supt. Resuspend cells in 2 mL Buffer P1.
- Buffer P1 **Qiagen Catalog #19051**
- .

Isolate Plasmid

- 5 2 mL Buffer P2 (invert 3-4 times)















Buffer P2 **Qiagen Catalog #19052**

- 6 Incubate 00:03:00 @ Room temperature



- 7 2 mL Buffer N3 (invert 3-4 times)

Buffer N3 **Qiagen Catalog #19064**



- 8 Distribute lysate into  4 x 1.5mL tubes
- 9 Centrifuge  00:10:00 @  13200 x g
- 10 Collect supt, add 1 volume 96% EtOH, mix
- 11 Load on  4 x "QIAprep 2.0 spin columns" or any generic silica membrane mini spin column (Econospin).
- 12 Spin  00:01:00 @  17900 x g
- 13 Re-load and re-spin until everything is loaded.
- 14  500 μ L Buffer PB to each column
 Buffer PB **Qiagen Catalog #19066**
- 15 Spin  00:01:00 @  17900 x g
- 16 750 uL Buffer PE
- 17 Spin  00:01:00 @  17900 x g
- 18 Spin again to dry  00:01:00 @  17900 x g
- 19 Elute with 30-50 uL H₂O



- 20 (Optional) Re-load H₂O on columns and re-spin increase yield.
-  00:01:00 @  17900 x g to