



Mar 27, 2023

Mitochondrial complex activity assays

DOI

dx.doi.org/10.17504/protocols.io.eq2ly753elx9/v1

michela.deleidi¹, Maria Jose Perez J.¹, Hariam Raji¹, Pascale Baden¹, Federico Bertoli¹

¹German Center for Neurodegenerative Diseases (DZNE), Tübingen, 72076 Germany



Federico Bertoli

ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly753elx9/v1>

Protocol Citation: michela.deleidi, Maria Jose Perez J., Hariam Raji, Pascale Baden, Federico Bertoli 2023. Mitochondrial complex activity assays. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.eq2ly753elx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 17, 2023



Last Modified: March 27, 2023

Protocol Integer ID: 78971

Keywords: Mitochondrial complex activity assays, mitochondrial electron transport chain , mitochondria complex activity assay, mitochondrial complex activity, mitochondrial electron transport chain, mitochondria, complex activity assay, activity levels of the different complex, activity level, different complex

Abstract

Mitochondria complex activity assays measure the activity levels of the different complexes of the mitochondrial electron transport chain (ETC).

Attachments



[676- 1425.docx](#)

13KB




Materials

Materials

- pyruvate
- malate
- ADP
- Succinate
- rotenone
- antimycin A
- TMPD (N,N,N',N'-tetramethyl-p-phenylenediamine dihydrochloride, Santa Cruz Biotechnology)
- ascorbic acid
- azide
- Qproteome Mitochondrial isolation kit

 MitoCheck® Complex I Activity Assay Kit **Cayman Chemical Company Catalog #700930**

 Qproteome Mitochondria Isolation Kit **Qiagen Catalog #37612**

MAS buffer

	A	B
	Sucrose	70 mM
	Mannitol	220 mM
	KH ₂ PO ₄	5 mM
	MgCl ₂	5 mM
	EGTA	1 mM
	HEPES pH 7.4	2 mM

Troubleshooting

Mitochondrial complex activity assays

35m

- 1 Isolate mitochondria from HEK cells, iPSC-derived neurons, or midbrain organoids using the Qproteome Mitochondrial isolation kit (QIAGEN, Cat. No. / ID: 37612) according to manufacturer's instructions.
- 2 Measure Complex I (NADH oxidase/coenzyme Q reductase) using the MitoCheck Complex I Activity Assay kit (Cayman Chemical, cat# 700930).
- 3 Determine the rate of NADH oxidation, which is proportional to CI activity, by a decrease in absorbance at 340 nm over 00:15:00 in the presence of ubiquinone and potassium cyanide to inhibit complex IV and prevent oxidation of ubiquinone.

15m

Note

To assess CI, CII, and CIV function, we used a respirometry approach based on XFp Extracellular Flux Analysis and then proceed with steps 4-9.

- 4 To this end, resuspend 3 mg of purified fresh mitochondria in 200 μ L of MAS buffer (70 millimolar (mM) sucrose, 220 millimolar (mM) mannitol, 5 millimolar (mM) KH_2PO_4 , 5 millimolar (mM) MgCl_2 , 1 millimolar (mM) EGTA, 2 millimolar (mM) HEPES 7.4) and seed in XFpSeahorse microplates.
- 5 Centrifuge the plate at 2000 x g, 4°C, 00:05:00 .

5m

- 6 Measure the OCR before and after the serial addition of pyruvate + malate (5 millimolar (mM) each) + ADP 3,5 mM or 1 mM Succinate + 4 micromolar (μ M) rotenone, 4 micromolar (μ M) rotenone + 8 micromolar (μ M) antimycin A, 0,5 mM TMPD (N,N,N',N'-tetramethyl-p-phenylenediamine dihydrochloride, Santa Cruz Biotechnology) + 1 millimolar (mM) ascorbic acid, and 50 millimolar (mM) azide.
- 7 Following each injection, record three measurements for a total period of 00:15:00 .
- 8 Calculate Complex I-, II-, and IV-dependent respiration by subtracting OCR values from the substrates (Pyruvate + malate + ADP for CI, Succinate + rotenone for CII and TMPD +

15m



ascorbic acid for CIV) subtracted from the ones from the inhibitors (rotenone for CI, antimycin A + rotenone for CII and azide for CIV).

- 9 Normalize the experimental values to the protein content per well via a BCA assay.