



Aug 21, 2025

Mitochondrial isolation from cultured fibroblasts

DOI

<https://dx.doi.org/10.17504/protocols.io.81wgbwxr1gpk/v1>

Ashley Seifert¹, Vekaria Hemendra², Ebenezer Aryee²

¹University of Kentucky; ²Medical University of South Carolina



Ashley Seifert

University of Kentucky

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbwxr1gpk/v1>

Protocol Citation: Ashley Seifert, Vekaria Hemendra, Ebenezer Aryee 2025. Mitochondrial isolation from cultured fibroblasts. protocols.io <https://dx.doi.org/10.17504/protocols.io.81wgbwxr1gpk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 18, 2025



Last Modified: August 21, 2025

Protocol Integer ID: 224833

Keywords: ASAPCRN, mitochondria, fibroblasts, mitochondrial isolation from cultured fibroblast, mitochondria from cultured fibroblast, mitochondrial isolation, mitochondria, cultured fibroblast, seifert lab

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020495

NIH

Grant ID: R01AR070313

Abstract

This protocol describes the procedure used by the Seifert Lab to isolate mitochondria from cultured fibroblasts.

Protocol materials

⊗ EDTA Merck MilliporeSigma (Sigma-Aldrich) Catalog #ED-1KG

⊗ Deionized water

⊗ 0.5 M EDTA Fisher Scientific Catalog #50-983-251



- 1 Dislodge cultured fibroblasts with trypsin
- 2 Spin cell down at  1500 rpm, 00:05:00
- 3 Discard supernatant and keep pellet  On ice
- 4 Resuspend pellet in  1 mL Mitochondria Isolation Buffer (MIB;
 215 millimolar (mM) mannitol,  75 millimolar (mM) sucrose, 0.1% BSA,
 20 millimolar (mM) HEPES, and  1 millimolar (mM) EGTA, adjusted to pH 7.2
with KOH)
- 5 Spin suspension at  400 rcf, 00:05:00
- 6 Discard supernatant
- 7 Resuspend pellet in  450 μ L
 EDTA Merck MilliporeSigma (Sigma-Aldrich) Catalog #ED-1KG
- 8 Break cells with a Nitrogen cell disruptor at  1250 Bar for  00:10:00
- 9 Pass lysates through a 26-gauge syringe topped with  1.5 mL MIB by spinning at
 1300 x g, 4°C, 00:03:00
- 10 Repeat Step 9
- 11 Transfer supernatant to a new tube and spin  1300 x g, 4°C, 00:10:00



- 12 Apply supernatant to Ficoll gradient ( 25 mL STE [ 8.55 g sucrose,  32 μ L  0.5 M EDTA Fisher Scientific Catalog #50-983-251 ,  80 μ L 1 M Tris, up to  25 mL with  Deionized water] and Ficoll stock [ 10 mL STE ,  17.5 mL 20% Ficoll])
- 13 Spin at  3200 rpm, 00:30:00
- 14 Discard supernatant
- 15 Resuspend pellet in  1 mL MIB
- 16 Spin at  13000 rcf, 00:10:00
- 17 Discard supernatant
- 18 Resuspend pellet (purified mitochondria) in MIB