


Mar 13, 2025

Mitochondrial Sub-fractionation Using Differential Centrifugation - Cells

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzkwpvx1/v1>



Julia Fitzgerald¹

¹Hertie Institute for Clinical Brain Research



Julia Fitzgerald

Hertie Institute for Clinical Brain Research

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzkwpvx1/v1>

Protocol Citation: Julia Fitzgerald 2025. Mitochondrial Sub-fractionation Using Differential Centrifugation - Cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzkwpvx1/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: March 13, 2025

Last Modified: March 13, 2025

Protocol Integer ID: 124326

Keywords: mitochondrial fractionation, crude mitochondria, enrichment, isolated, cells crude fractionation, cells crude fractionation of cell, centrifugation, using differential centrifugation, cell, mitochondria

Disclaimer

NA

Abstract

Crude fractionation of cells to yield crude mitochondria and subfractions using differential centrifugation

Guidelines

GM cells should be treated according to local safety guidelines

Materials

Eppendorf tubes

Dounce tissue grinder (loose and tight fitting)

Tris

EDTA

Sucrose

Centrifuge capable of 10,000g

Ultracentrifuge (optional)

Protocol materials

 1X PBS (Phosphate-buffered saline)

Safety warnings

 No warnings

Ethics statement

NA

Before start

Prepare and cool the buffers

Prepare ice

Process collected cells and prepare crude mitochondria

2h 5m

- 1 Start with a cell pellet (for crude mitochondrial preparation start with a minimum of a confluent 10cm dish, for sub-fractionation start with minimum of 3-6 confluent T75cm flasks). The cell pellet is always kept cold  On ice
- 2 Wash the cell pellet once with  1X PBS (Phosphate-buffered saline) and  1200 rpm, 21°C, 00:05:00 to generate a washed cell pellet.  On ice
- 3 Dislodge the cell pellet with approximately  1000 μ L of TES buffer (70 mM Tris base, 1mM EDTA, 0.25 M sucrose, pH 7.4) and transfer to a loose fit hand held homogenizer (such as Eppi-pestles +1.5mL tube). With strong force by hand, turn the pestles ten times. Transfer the homogenate to a tight fitting glass homogenizer (such as dounce tissue grinder Sigma-Aldrich P1110-1EA) and with strong force by hand dounce the homogenate by turning the glass pestle as you move up and down in the glass mortar tube 10 times.  On ice

Equipment

Dounce Tissue Grinder Set with Two Glass Pestles

NAME

Wheaton

BRAND

cat# 357538

SKU

<https://us.vwr.com/store/product/4832641/dounce-tissue-grinders-wheaton> ^{LINK}



Equipment

Eppi-pestles

NAME

Pestle for 1.5mL tube

TYPE

Schuett biotec

BRAND

3200512

SKU

<https://www.schuett-biotec.de/>

LINK

- 4 Transfer the homogenate to a 2mL round bottomed "Eppendorf" tube (Eppi). Optional: aliquot a small amount of the homogenate in a new tube labelled "TE" total extract or homogenate. 🧊 On ice

- 5 🌀 1000 x g, 4°C, 00:10:00

10m

- 6 Carefully remove the supernatant and transfer to a fresh 1.5mL Eppi and label "M" 🧊 On ice . Add 500µL TES buffer to the "nuc" pellet and homogenize again (see above). 🌀 1000 x g, 4°C, 00:10:00 . Transfer the supernatant to the "M" tube on ice. Label the nuc pellet tube "N" and store if needed. 🧊 On ice

10m

- 7 Tube "M" > 🌀 9800 x g, 4°C, 00:15:00

15m

- 8 Transfer the supernatant to a new Eppi labelled "C" for cytosolic fraction. 🧊 On ice

- 9 Look carefully at the "M" pellet which is crude mitochondria. The pellet should be red/brown/yellow in colour and flaky/hydrophobic (see image in tab "description" of this protocol). If the mitochondrial pellet is surrounded by contamination (white fluffy pellet),



the contaminating white fluffy pellet can be carefully removed with a pasteur pippette pump or a pipette.

- 10 Optional: wash the mitochondrial pellet in TES buffer and centrifuge again

9800 x g, 4°C, 00:10:00

10m

- 11 Resuspend the mitochondrial pellet in TMG buffer (10 mM Tris, 5 mM MgCl₂, 20 % (v/v) glycerol) and perform protein estimation. Working range: minimum 0.7mg/mL mitochondrial protein concentration for most applications and assays.

On ice

- 12 Optional: aliquot a small amount of the crude mitochondrial extract.

On ice

- 13 Optional: to store the crude mitochondrial fraction for further use;

9800 x g, 4°C, 00:15:00 and store as a dry pellet at -80 °C

15m

Mitochondrial sub-fractions

10m

- 14 Optional: To immunoprecipitate from crude mitochondria, first resuspend the mitochondrial pellet in TES buffer and sonicate for 1 min at 5s intervals (14uA amplitude)

On ice

- 15 Incubate the crude mitochondrial fraction in ice cold digitonin (0.2 mg/mg mitochondrial protein) rolling or turning at 4 °C for 15 minutes.

- 16 9800 x g, 4°C, 00:10:00

10m

- 17 Pellet = mitoplast (IMM and matrix). Supernatant = OMM/IMS.

- 18 Optional: store supernatant or pellet resuspended in TES at -80 °C

Further fractionation

1h



- 19 Disrupt mitoplast by sonicating for 1 min at 5s intervals (14uA amplitude)  On ice . Do the same for the OMM/IMS fraction  On ice .
- 20 Centrifuge all disrupted fractions: mitoplast, OMM/IMS and if needed cytosolic fraction at high speed.  100000 x g, 4°C, 01:00:00 1h
- 21 Mitoplast: Supernatant = matrix, pellet = IMM
OMM/IMS: Supernatant = IMS, pellet = OMM
Cytosol: Supernatant=purified cytosol, pellet = membrane contaminants.
- 22 Resuspend pellets in 100-200 µL TES and store all fractions at -80 C or use straight away.

Acknowledgements

Modified after a protocol developed for tissue by Prof. Ellen Billett (retired). Nottingham Trent University, UK