



Apr 05, 2023

[Modified] DNeasy PowerSoil Pro Kit_Increased Sediment Volume & Optional Concentration

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzb9x4vx1/v1>

Grayson Huston¹

¹Maine Center for Genetics in the Environment, University of Maine



Grayson Huston

University of Maine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzb9x4vx1/v1>

Protocol Citation: Grayson Huston 2023. [Modified] DNeasy PowerSoil Pro Kit_Increased Sediment Volume & Optional Concentration. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzb9x4vx1/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Other

Protocol successful at detecting fish sedDNA collected in a stream during a fish migration Protocol unsuccessful at detecting fish sedDNA from Maine lakes

Created: March 10, 2023

Last Modified: April 05, 2023

Protocol Integer ID: 78502

Keywords: detecting fish sedDNA, fish sedDNA from lake, fish sedDNA, dneasy powersoil pro kit-increased sediment volume, streams during anadromous fish sea, anadromous fish sea, lake, increased sediment amount, sediment amount

Abstract

Protocol (both increased sediment amount up to 2.0g as well as concentrating DNA post-extraction) unsuccessful at detecting fish sedDNA from lakes in Maine, USA.

Both protocols successful at detecting fish sedDNA collected in streams during anadromous fish sea-run migrations



Modified PowerSoil Pro extraction - sample preparation & cell lysis

31m

1 **CENTRIFUGE** sediment samples briefly to separate pore water

DISCARD pore water to retain only sediment samples

2 **SPIN** the PowerBead Pro Tube briefly to ensure that the beads have settled at the bottom

ADD up to  2.0 g of wet sediment to the PowerBead Pro Tube

Note

Qiagen recommends  .25 g

In this test,  0.25 g ,  0.5 g ,  1.0 g , and  2.0 g of sediment were tested

ADD  800 μ L Solution CD1

VORTEX briefly to mix

3 **SECURE** PowerBead Pro Tubes horizontally to a 1.5mL-2.0mL Vortex Adapter

20m

VORTEX for  00:10:00

ROTATE tubes so caps are oriented in opposite direction

VORTEX for another  00:10:00

4 **CENTRIFUGE** PowerBead Pro Tube at  15000 x g for  00:01:00

1m

TRANSFER all supernatant to a clean 2 mL Microcentrifuge Tube

Modified PowerSoil Pro extraction - inhibitor removal

31m

5 **ADD**  200 μ L of Solution CD2

VORTEX briefly to mix

6 **CENTRIFUGE** at  15000 x g for  00:01:00

1m

AVOIDING the pellet, transfer **all** supernatant to a clean 2 mL Microcentrifuge Tube



Modified PowerSoil Pro extraction - bind DNA

31m

7 **ADD**  600 μL of Solution CD3

VORTEX briefly to mix

8 **LOAD**  650 μL of the lysate onto a MB Spin Column

1m

CENTRIFUGE at  15000 x g for  00:01:00

DISCARD the liquid flow-through

9 **REPEAT** step 8 to ensure all of the lysate has passed through the MB Spin Column

CAREFULLY place the MB Spin Column into a clean 2 mL Collection Tube

Modified PowerSoil Pro extraction - wash spin column

31m

10 **ADD**  500 μL of Solution EA to the MB Spin column

1m

CENTRIFUGE at  15000 x g for  00:01:00

DISCARD the liquid flow-through and place the MB Spin Column into same 2 mL Collection Tube

11 **ADD**  500 μL of Solution C5 to the MB Spin Column

1m

CENTRIFUGE at  15000 x g for  00:01:00

DISCARD the liquid flow-through and place the MB Spin Column into a **new** 2 mL Collection Tube

12 **CENTRIFUGE** at  16000 x g for  00:02:00

2m

CAREFULLY place the MB Spin Column into a new 1.5mL Elution Tube

Modified PowerSoil Pro extraction - elute the DNA

31m

13 **ADD**  100 μL of Solution C6 to the center of the white membrane in the MB Spin Column

2m



INCUBATE at Room temperature for 00:01:00

CENTRIFUGE at 15000 x g for 00:01:00

- 14 **PIPETTE** the liquid flow-through and re-add it to the center of the white membrane in the MB Spin Column

2m

INCUBATE at Room temperature for 00:01:00

CENTRIFUGE at 15000 x g for 00:01:00

DISCARD the MB Spin Column

DNA is now ready for downstream applications

(Optional) DNA Concentration with PALL Nanosep 30K Centrifugal Devices

2m

- 15 **ENSURE** that the sample reservoir is firmly placed into the filtrate receiver

PIPETTE 50-100 μ L of DNA extract into the sample reservoir

CAP the Nanosep device

- 16 **CENTRIFUGE** at 5000 x g for 00:02:00

2m

RECOVER concentrated sample from the sample reservoir with a micropipette

TRANSFER concentrated sample to a new 1.5 mL Microcentrifuge Tube

Concentrated DNA is now ready for downstream applications

Note

If the sample appears to have "spun dry", recover the sample by pipetting ~ 20uL of elution buffer onto the membrane and recovering with a micropipette