

Apr 07, 2023

## Modified Organic Extraction Protocol (Pedersen et al., 2016)

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1o3oklr2/v1>

Mikkel Pedersen<sup>1</sup>

<sup>1</sup>Centre for GeoGenetics, Natural History Museum, University of Copenhagen



Mark Louie Lopez

University of Victoria

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1o3oklr2/v1>

**Protocol Citation:** Mikkel Pedersen 2023. Modified Organic Extraction Protocol (Pedersen et al., 2016). **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.ewov1o3oklr2/v1>

**Manuscript citation:**

Pedersen, M., Ruter, A., Schweger, C. et al. Postglacial viability and colonization in North America's ice-free corridor. Nature 537, 45–49 (2016). <https://doi.org/10.1038/nature19085>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** January 20, 2023

**Last Modified:** April 19, 2023

**Protocol Integer ID:** 75593

**Keywords:** Sedimentary DNA, SedDNA, version of the organic extraction protocol, organic extraction protocol, modified organic extraction protocol, kya sediment core, extraction, biodiversity

## Abstract

A modified version of the organic extraction protocol used in Wales et al., 2014.

Protocol successful at reconstructing the biodiversity of 15 kya sediment core, including a few fish species

## Guidelines

Phenol should be handled with care as it is corrosive and toxic.

Chloroform is toxic if swallowed or inhaled. It can cause severe and irreversible health effects, including death

## Materials

Lysis buffer

Proteinase K

MO BIO's PowerClean DNA Clean Up Kit

Phenol: Chloroform: Isoamyl alcohol

Ethanol

TE Buffer

30 kDa Amicon Ultra-4 centrifugal filters

## Safety warnings

 Skin exposure to large amounts of phenol or chloroform can cause sores. Skin exposure to lesser amounts of phenol or chloroform can cause irritation. Phenol and chloroform can irritate the eyes that encounter it.

## Before start

Review the manufacturer's Safety Data Sheet and additional chemical information. Ensure that a written experimental protocol including safety information is available.



## Sample preparation & cell lysis

53m 40s

1 **ADD**  2 g of sediment sample to a sterile spin tube

40s

**ADD**  170 µg of proteinase K to the sample

**ADD**  3 mL of lysis buffer to the sample

**VORTEX** sample at max speed for  00:00:40

### Note

Lysis buffer:

-  68 millimolar (mM) *N*-lauroylsarcosine sodium salt
-  50 millimolar (mM) Tris-HCl (pH 8.0)
-  150 millimolar (mM) NaCl
-  20 millimolar (mM) EDTA (pH 8.0) and,

Immediately before extraction, for each  30 mL of lysis buffer add:

-  1.5 mL 2-mercaptoethanol
-  1 mL 1M DTT

2 **ADD** an additional  170 µg of proteinase K to each sample

**INCUBATE**  Overnight at  37 °C

## Inhibitor removal

53m 40s

3 **REMOVE** PCR inhibitor using MOBIO (MO BIO Laboratories, Carlsbad, CA) C2 and C3 buffers following the manufacturer's protocol

3m

### Note

Link to MO BIO's PowerClean DNA Clean Up Kit Handbook:

<https://www.qiagen.com/us/resources/resourcedetail?id=a757e687-7b7a-4801-96bb-6267620414de&lang=en>



**CENTRIFUGE** samples at  10000 x g for  00:03:00

**TRANSFER** supernatant to a sterile microcentrifuge tube

## Phenol:Chloroform extraction

53m 40s

4 **ADD** equal volume of Phenol: Chloroform: Isoamyl alcohol (25: 24: 1) solution to samples

15m

**ROTATE** samples gently at  Room temperature for  00:10:00

**CENTRIFUGE** at  3200 x g for  00:05:00

**DISCARD** supernatant

5 **REPEAT** step 4

**TRANSFER** supernatant to a Millipore 15 mL Amicon Ultra filter for concentration and purification

## DNA concentration and purification

40m

6 **CENTRIFUGE** filter at  4000 x g for  00:10:00

10m

**DISCARD** liquid flow-through

7 **ADD**  500 µL EB Buffer (Qiagen) to the filter

10m

**CENTRIFUGE** at  4000 x g for  00:10:00

**DISCARD** liquid flow-through

8 **REPEAT** step 7

9 **ADD**  50 µL molecular grade water to the filter

20m

**INCUBATE** at  Room temperature for  00:10:00

**CENTRIFUGE** at  4000 x g for  00:10:00



**DISCARD** filter

**DNA is now ready for downstream applications**

## Protocol references

Pedersen, M., Ruter, A., Schweger, C. et al. (2016). Postglacial viability and colonization in North America's ice-free corridor. *Nature* 537, 45–49. <https://doi.org/10.1038/nature19085>

Wales, N., Andersen, K., Cappellini, E., Avila-Arcos, M. C. & Gilbert, M. T. Optimization of DNA recovery and amplification from non-carbonized archaeobotanical remains. *PLoS One* **9**, e86827 (2014)