



Sep 26, 2023

Modified salting out method for high molecular weight gDNA extraction (oribatid mites)

DOI

dx.doi.org/10.17504/protocols.io.yxmvm3yybl3p/v1

Hüsna Öztoprak¹, Jens Bast¹

¹Institute for Zoology, University of Cologne, Köln, Germany



h.oeztoprak

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm3yybl3p/v1>

Protocol Citation: Hüsna Öztoprak, Jens Bast 2023. Modified salting out method for high molecular weight gDNA extraction (oribatid mites). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvm3yybl3p/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working. It enables the generation of HMW gDNA of ~20 ng, even up to 200 ng of a single oribatid mite, depending on the genus. It also works for collembolans and nematodes.

Created: September 25, 2023

Last Modified: September 26, 2023

Protocol Integer ID: 88360

Keywords: dna extractions from oribatid mite, weight gdna extraction protocol for mite, method for high molecular weight gdna extraction, molecular dna from oribatid mite, high molecular weight gdna extraction, weight gdna extraction protocol, dna extraction with chitinase, weight genomic dna extraction method, molecular gdna extraction, molecular gdna extraction ii, dna extraction, dna purification, oribatid mite, more gdna, genomic dna, mite, dna, additional loss of dna, extraction, molecular dna, chitinase, single minuscule specimen, phased genome, chitinous exoskeleton

Funders Acknowledgements:

German Research Foundation Emmy Noether

Grant ID: grant BA 5800/4-1

Abstract

This protocol describes a low-cost, high-molecular-weight genomic DNA extraction method for a single minuscule specimen (modified from Miller et al 1988). DNA extractions from oribatid mites are typically challenging, because of the small body size of 150-1400 µm. Their chitinous exoskeleton does not dissolve during DNA extraction, impedes DNA purification, and leads to additional loss of DNA. Therefore, high-molecular DNA from oribatid mites has been thus far unattainable, especially from single individuals. We established a high-molecular-weight gDNA extraction protocol for mites that enables the generation of high-quality phased genomes for small non-model organisms. There are three options to utilize this protocol: i) for high-molecular gDNA extraction ii) for high-molecular gDNA extraction, while preserving the exoskeleton for morphological analysis, and iii) DNA extraction with Chitinase to yield more gDNA.

As specimens are collected from natural populations and are not cultured in the lab. They are cleansed prior to DNA extraction to minimize external contamination. Cleansing includes brushing the specimen in distilled water and in distilled water with detergent (fit GmbH, Zittau, Germany). Once the external residue is not visible anymore specimen is incubated in NaClO 0.05% (DonKlorix; CP GABA GmbH, Hamburg, Germany) and ethanol 70% for 30 seconds each and rinsed in distilled water again.

Attachments



[nn7eb8kcf.pdf](#)

111KB

Materials

Material and Regents

 Proteinase K, 2mL **Qiagen Catalog #19131**

 Yeast tRNA (10 mg/mL) **Thermo Fisher Scientific Catalog #AM7119**

 RNase Cocktail[®]; Enzyme Mix **Thermo Fisher Catalog #AM2286**

Solutions

1. TNES buffer (see Recipes)

Recipes

1. Final concentration of TNES buffer freshly made before each extraction, ddH₂O used to dilute.:


	A	B
	NaCl	400 mM
	EDTA	20 mM
	Tris pH 8.0	50 mM
	SDS	0.5%

Troubleshooting





Version i) High-molecular gDNA extraction: DNA Extraction



2h 5m 3s

- 1 Submerge cleansed specimen in  195 μL TNES buffer and flash freeze by holding tube in liquid nitrogen.
- 2 Using a sterile pestle, homogenize by applying pressure to grind the specimen between pestle head and the walls of the tube.

Note






If low DNA yield is expected. Leave pestle in tube to ensure maximum digestion of material. Consider including multiple (3x) freeze and thaw cycles and vortexing to disrupt tissue fully.

- 3 Add  5 μL proteinase K.
- 4 Vortex for  00:00:03 . Centrifuge briefly.



- 5 Incubate at  55 $^{\circ}\text{C}$ for ~  01:00:00 .

Note

Completely dissolve the specimen.

- 5.1 If there is indigestible debris left over centrifuge  18000 rcf, 00:05:00 and transfer the supernatant to a fresh tube.
- 6 Add  1.5 μL yeast tRNA, flick to mix briefly then spin down.
- 7 Add  65 μL  5 Molarity (M) NaCl and  290 μL 96% EtOH, mix by inversion.



8 This should clarify the solution. Store at  -20 °C for  01:00:00 .

1h

Note

Solution can be stored  Overnight at this step.

Version i) High-molecular gDNA extraction: DNA Purification

2h 20m

9 **Optional:** add  1 µL Pellet Paint Co-Precipitant (for colorful pellet).



10 Spin down at  18000 rcf, 4°C, 00:15:00 .

15m

Note

Know the expected position of the DNA pellet, as it can be difficult to see.



11 Ghostly pellet should be visible.

12 Remove supernatant.

13 Add  0.5 mL chilled 70% EtOH (make fresh).



14 Spin at  18000 rcf, 00:05:00 .

5m



15 Repeat ethanol rinse.





16 Carefully remove supernatant.

17 Leave tube open to air dry. Pellet should have a glassy appearance.






18 Using a wide-bore pipette tip, add  21 μL TE Buffer (elution buffer) and gently resuspend the DNA pellet with pipette mixing.



19 Let DNA resuspend at  4 °C  Overnight .

1h



20 Add  2 μL RNase Cocktail. Incubate at  37 °C for approx.  01:00:00 .

1h

Note

Note: For femto pulse systems elution in TE Buffer, more specifically EDTA is not recommended. Alternatively use 0.1 mM EDTA or EB, Tris-HCl (pH 8-8.5).




Version ii) High-molecular gDNA extraction preserving exoskeleton: DNA Extraction


2h

21 To observe the specimen under a microscope, cleansed specimens are placed on a sterile slide and submerged with TNES buffer until fully covered.



22 Remove one genital plate cautiously with a sharp needle and stir tissue without destroying the exoskeleton.

23 Transfer specimen in  195 μL TNES buffer.



24 Add  5 μL proteinase K.




25 Incubate at  37 °C  Overnight .

1h











26 Transfer specimen with sterile needle to a tube containing 70% EtOH and store it for morphological analysis.

27 Add  1.5 μL yeast tRNA, flick to mix briefly then spin down.





28 Add  65 μL  5 Molarity (M) NaCl and  290 μL 96% EtOH, mix by inversion.  


29 This should clarify the solution. Store at  -20 °C for  01:00:00 . 

Note


DNA Purification follows the same procedure as above.



Version iii) gDNA extraction with chitinase digestion: DNA Extraction

3h 5m 10s



30 Submerge specimen in  195 μL TNES buffer and flash freeze by holding tube in liquid N.

31 Using a sterile pestle, homogenize by applying pressure to grind the specimen between pestle head and the walls of the tube.



32 Add  2 μL chitinase (1 U/ml). 



33 Vortex for  00:00:05 . Centrifuge briefly. 






34 Incubate at  55 °C for ~  01:00:00 . 




35 Add  5 μL proteinase K. 

36 Vortex for  00:00:05 . Centrifuge briefly. 

37 Incubate at  55 °C for ~  01:00:00 (completely dissolve the specimen). 









37.1 If there is indigestible debris left over centrifuge  18000 rcf, 00:05:00 and transfer the supernatant to a fresh tube.

5m





38 Add  1.5 μL yeast tRNA, flick to mix briefly then spin down.



39 Add  65 μL  5 Molarity (M) NaCl and  290 μL 96% EtOH, mix by inversion.



40 This should clarify the solution. Store at  -20 °C for  01:00:00 .

1h

Note

DNA Purification follows the same procedure as above.

Protocol references

References

Laumer, C. 2023. "Picogram input multimodal sequencing (PiMmS)." protocols.io <https://dx.doi.org/10.17504/protocols.io.rm7vzywy5lx1/v1>

Miller, S. A., D. D. Dykes, and H. F. Polesky. 1988. "A Simple Salting out Procedure for Extracting DNA from Human Nucleated Cells." Nucleic Acids Research 16 (3): 1215