


Mar 12, 2020

# Monitoring cell-surface expression of GPCR by ELISA

DOI

[dx.doi.org/10.17504/protocols.io.zfef3je](https://dx.doi.org/10.17504/protocols.io.zfef3je)



Elie Besserer-Offroy<sup>1</sup>, Rebecca L Brouillette<sup>2</sup>, Jean-Michel Longpré<sup>2</sup>, Philippe Sarret<sup>2</sup>

<sup>1</sup>Inserm U1086 - Anticipo;

<sup>2</sup>Department of Pharmacology-Physiology, Faculty of Medicine and Health Sciences, Université de Sherbrooke

The Besserer Lab

Tech. support email: [protocol\\_support@bessererlab.org](mailto:protocol_support@bessererlab.org)



**Elie Besserer-Offroy**

Inserm U1086 - Anticipo, Centre de lutte contre le cancer Fr...

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.zfef3je>

**Protocol Citation:** Elie Besserer-Offroy, Rebecca L Brouillette, Jean-Michel Longpré, Philippe Sarret 2020. Monitoring cell-surface expression of GPCR by ELISA. [protocols.io https://dx.doi.org/10.17504/protocols.io.zfef3je](https://dx.doi.org/10.17504/protocols.io.zfef3je)

**Manuscript citation:**

Gusach A, Luginina A, Marin E, Brouillette RL, Besserer-Offroy É, Longpré JM, Ishchenko A, Popov P, Patel N, Fujimoto T, Maruyama T, B Stauch, Ergasheva M, Romanovskaia D, Stepko A, Kovalev K, Shevtsov M, Gordeliy V, Han GW, Katritch V, Borshchevskiy V, Sarret P, Mishin A, Cherezov V. Structural basis of ligand selectivity and disease mutations in cysteinyl leukotriene receptors. Nature Communications. 2019; (10):5573. PMID: [31811124](https://pubmed.ncbi.nlm.nih.gov/31811124/).

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 22, 2019

**Last Modified:** March 12, 2020

**Protocol Integer ID:** 21702

**Keywords:** G Protein-Coupled Receptor, Seven Transmembrane Receptor, Cell Surface Expression, ELISA, HA-tag, HRP-linked anti-HA, elisa quantifying cell surface expression, surface expression of gpcr, gpcr, linked immunosorbent assay, immunosorbent assay, mutant receptor, surface expression, monitoring cell, coupled receptor, surface detection, assay, elisa,

## Abstract


Quantifying cell surface expression of G Protein-Coupled Receptors (GPCRs) can be extremely important for the expression of mutant receptors. Herein we report a useful Enzyme-Linked Immunosorbent Assay (ELISA) for the cell-surface detection of a HA-tagged version of a GPCR.

## Materials

### MATERIALS

 NESTLE CARNATION Instant Non-fat Dry Milk

 Poly-l-lysine, 0.1% (wt/vol) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P8920**

 Falcon® 24-Well Flat-Bottom Plate, Tissue Culture-Treated 50 Plates **STEMCELL Technologies Inc. Catalog #38021**


 HEK293 **ATCC Catalog #CRL-1573**

 Lipofectamine 3000 **Thermo Fisher Scientific Catalog #L3000015**


 Opti-MEM® I Reduced Serum Medium **Thermo Fisher Catalog #31985062**


 Sterile Water **Wisent Bioproducts Catalog #809-115-CL**

 PBS 1X **Wisent Bioproducts Catalog #311-011-CL**


 Trypsin 0.25% / EDTA 2.21 mM in HBSS **Wisent Bioproducts Catalog #325-043-EL**


 SafeSeal tube 5mL **Sarstedt Catalog #72.701**

 Formaldehyde Reagent Grade **Bioshop Catalog #FOR201**


 3 3'5 5'-Tetramethylbenzidine Liquid Substrate Supersensitive for ELISA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T4444-100ML**

 Anti-HA-Peroxidase High Affinity **Merck MilliporeSigma (Sigma-Aldrich) Catalog #12013819001**

 FBS (Fetal Bovine Serum) Premium Quality Endotoxin **Wisent Bioproducts Catalog #080-150**

 DMEM 4.5g/L glucose with L-glutamine sodium pyruvate and phenol red **Wisent Bioproducts Catalog #319-005-CL**

 HEPES 1M Free acid **Wisent Bioproducts Catalog #330-050-EL**






 Penicillin (5000IU) / Streptomycin (5000µg/mL) sterile filtered for cell culture **Wisent Bioproducts Catalog #450-200-EL**

Endotoxin-free purified plasmidic DNA encoding for HA-tagged GPCRs


Tris-Buffered Saline (TBS, containing 20mM Tris–HCl pH 7.5 and 150mM NaCl)

## Troubleshooting

## Day 1 - Cell Culture & Transfections

- 1 **Coat 24-well plates with Poly-L-Lysine** (this need to be done in a biological safety cabinet to ensure sterility).
  - 1.1 Add  300  $\mu$ L of [M] 0.1 mg/mL Poly-L-Lysine solution in each well of the 24-well plate and incubate  00:10:00 at  Room temperature .  
This can be done using a 300 $\mu$ L multichannel pipet fitted with 4 tips
  - 1.2 Remove the Poly-L-Lysine solution (this solution can be re-used up to 4 times to coat cell culture plasticware).  
This can be done using a 300 $\mu$ L multichannel pipet fitted with 4 tips
  - 1.3 Rinse the wells twice with  300  $\mu$ L of sterile water.
  - 1.4 Let dry the 24-well plate without lid under the biological safety cabinet for  00:20:00 before seeding cells.

### Note

Poly-L-Lysine-coated plates can be stored for several weeks at  Room temperature before use.

- 2 **Prepare transfections of plasmids encoding HA tagged-GPCRs** (this need to be done in a biological safety cabinet to ensure sterility).


### Note


In the case you want to compare expression of a mutant receptor, positive and negative controls are needed to normalize the results. The positive control should be the wild-type receptor and the negative control should be cells transfected with the empty vector (MOCK cells).




For 3 well of the 24-well plate:

- 2.1 Add  300  $\mu$ L of Opti-MEM into a sterile 5mL tube.



2.2 Add  1.5  $\mu\text{g}$  of plasmidic DNA encoding for the desired HA tagged-GPCR to the tube containing Opti-MEM and mix.

2.3 Add  3  $\mu\text{L}$  of P3000 Reagent to the tube and mix.




2.4 Add  2.25  $\mu\text{L}$  of Lipofectamine 3000 to the tube, mix, and incubate at  Room temperature for  00:15:00 .


### 3 **Prepare HEK293 cells for transfection** (this need to be done in a biological safety cabinet to ensure sterility).

#### Note

Ideally, cells were seeded at a density of 3 million cells per 10cm-petri dish 48h before transfection to ensure a high transfection rate.

3.1 Remove culture media and rinse cells with PBS.

3.2 Add  1 mL of 0.25% Trypsin to a 10-cm petri dish and incubate for  00:02:00 at  37  $^{\circ}\text{C}$  .

3.3 Add  5 mL of complete DMEM (10% FBS, 20mM HEPES, Penicilin/Streptomycin) to the petri dish and dissociate cells by pipeting up and down.

3.4 Count cells using an automated cell counter or a hemacytometer.



### Equipment

MOXI Z Mini

NAME

Automated Cell Counter

TYPE

Orflo

BRAND


MXZ001


SKU


[https://www.orflo.com/product\\_p/mxz001.htm](https://www.orflo.com/product_p/mxz001.htm)

LINK

3.5 Adjust cell concentration to 150,000 cells/mL.

3.6 Add  2 mL of the cell suspension at 150,000 cells/mL to the 5mL tube containing the plasmidic DNA and mix gently.

3.7 Dispense  450 µL of the mix of cell and plasmidic DNA to the desired wells of the Poly-L-Lysine-coated 24-well plate.

3.8 Incubate at  37 °C in humidified chamber at 5% CO<sub>2</sub> for  48:00:00 .

## Day 3 - ELISA

4 **Detection of cell surface expression by ELISA** (this part of the protocol can be done on the wet bench).



### Note

Buffer and reagents in this part of the protocol can be dispensed to the 24-well plate using combitips and a single channel repeater.

### Equipment

#### M4 Repeater

NAME

Multidispense pipet

TYPE

Eppendorf








BRAND

4982000322

SKU

<https://online-shop.eppendorf.ca/CA-en/Manual-Liquid-Handling-44563/Manual-Pipetting--Dispensing-44564/Repeater-M4-PF-44619.html>



LINK

- 4.1 Remove cell culture media and wash each well with  400  $\mu$ L of PBS.
- 4.2 Fix cells using  400  $\mu$ L of [M] 3.7 Mass Percent formaldehyde in TBS for  00:05:00 at  Room temperature .
- 4.3 Rinse cells three times with  500  $\mu$ L of TBS.
- 4.4 Block non-specific sites using  400  $\mu$ L of [M] 3 Mass Percent non-fat dry milk dissolved in TBS for  00:30:00 under gentle orbital agitation.



### Note







In case non-fat dry milk is not suitable for blocking non-specific sites, a solution of [M] 1 Mass Percent of Bovine Serum Albumine in TBS can be used. The same solution must also be used for the incubation of the antibody.



- 4.5 Remove blocking solution and add  250  $\mu\text{L}$  of 1/1000 dilution of HRP-linked anti-HA antibody diluted in  3 Mass Percent non-fat dry milk in TBS for  03:00:00 under gentle orbital agitation.

#### Note

This incubation step can also be done  Overnight at  4 °C under gentle orbital agitation.

- 4.6 Remove antibody solution and wash each well three times with  500  $\mu\text{L}$  of TBS.
- 4.7 Add  250  $\mu\text{L}$  of  Room temperature 3, 3',5 ,5'-Tetramethylbenzidine (TMB) Liquid Substrate, Supersensitive, for ELISA and incubate under gentle orbital agitation for 2 to 15 min (until the color of your positive control turn intense blue).
- 4.8 Stop the TMB reaction by adding  250  $\mu\text{L}$  of  2 Molarity (M) Hydrochloric Acid (HCl).
- 4.9 Transfer  100  $\mu\text{L}$  of the colorimetric reaction to a flat-bottom transparent 96-well plate.
- 4.10 Read the absorbance at 450nm using a multimode plate reader.



### Equipment

Mithras2 LB943	NAME
Multimode plate reader	TYPE
Berthold	BRAND
LB943	SKU

## Results analysis

- 5 As this quantification of cell surface expression is a semi-quantitative method it should not be presented as raw OD<sub>450nm</sub> values but rather as a percentage of expression compared to the positive control (or wild-type receptor).

To normalize the results, average the OD<sub>450nm</sub> of the positive control and the OD<sub>450nm</sub> of the negative control and apply the following formula:

$$y = \frac{x - \overline{x_{min}}}{\overline{x_{max}} - \overline{x_{min}}} \times 100$$

#### Normalization formula

y = normalized value

x = OD<sub>450nm</sub> value of the sample

x<sub>min</sub> = mean OD<sub>450nm</sub> value of the negative control

x<sub>max</sub> = mean OD<sub>450nm</sub> value of the positive control