



Dec 21, 2023

Mouse Brain Tissue Collection and Analysis

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4qj4zvo5/v1>

Robert Edwards¹, Shweta Jain¹

¹University of California San Francisco



Kelsey Barcomb

CU Anschutz

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4qj4zvo5/v1>

Protocol Citation: Robert Edwards, Shweta Jain 2023. Mouse Brain Tissue Collection and Analysis . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.3byl4qj4zvo5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 09, 2023



Last Modified: May 31, 2024

Protocol Integer ID: 86289

Keywords: ASAPCRN, Mouse, Brain, Catecholamine, Dopamine, mouse brain tissue collection, distribution of the vesicular monoamine transporter, vesicular monoamine transporter, measurement catecholamine level, neurochemistry core, mouse brain, striatum, tissue

Funders Acknowledgements:

ASAP-CRN

Grant ID: 020529

Abstract

This protocol describes the dissection and collection of coronal sections of the striatum and midbrain from a mouse brain. The tissue can be used in a number of applications and here we describe two: the measurement catecholamine levels using high-performance liquid chromatography (aided by a neurochemistry core) and the distribution of the vesicular monoamine transporter (VMAT2) using a radioligand binding assay.

Tissue dissection

- 1 Prepare solutions
 - 1.1 For measurement of brain catecholamine by HPLC, use cold Hank's Balanced Salt Solution (HBSS) containing 10 mM HEPES and 20 mM glucose
 - 1.2 For [³H]-Dihydrotetrabenazine binding, dissection into cold SHT buffer (320 mM sucrose, 10 mM HEPES/Tris, pH 7.4) with 0.4 mM EDTA and Complete Protease Inhibitor Cocktail (Roche)
- 2 Euthanize mice by inhalation of CO₂
- 3 Decapitate mouse and remove brain
- 4 Place brain into a rodent brain matrix (RBM-2000C, Protech International Inc.)
- 5 By inserting razor blades into the slots of the brain matrix, collect two 1 mm coronal sections containing the striatum

Note

Be sure to flip the slices and look at both sides, selecting the sections that contain most of the striatum

- 6 Remove the cortical tissue and cut out striatum, collecting tissue in cold solution (HBSS or SHT depending on application)
- 7 From the same brain, dissect the midbrain using the brain matrix
Transfer this region onto parafilm. Remove top 70% of the dorsal side and collect the ventral side (containing the VTA and substantia nigra) into cold solution (HBSS or SHT depending on application). Do not separate hemispheres; there will be one sample per mouse.

Measurement of catecholamine levels by HPLC

- 8 Transfer tissue to Eppendorf tubes and flash freeze in liquid nitrogen. Immediately transfer to  -80 °C .
- 9 Tissue is shipped on dry ice to the Vanderbilt Neurochemistry Core (<https://lab.vanderbilt.edu/vbi-core-labs/neurochemistry-core/>), where tissue catecholamine levels is measured by HPLC with coupled electrochemical detection.

[3H]-Dihydrotetrabenazine binding

- 10 Disrupt midbrain and striatal tissue with 12 strokes of a Dounce homogenizer at 500 RPM in cold SHT buffer
- 11 Sonicate tissue for 30 seconds
- 12 Sediment the debris in a centrifuge at 2000 g for 2 minutes
- 13 Collect supernatant
- 14 Measure protein content with bicinchoninic (BCA) assay
- 15 Dilute 50 µg protein into SHT buffer and add 10 nM (+)- α -dihydrotetrabenazine [9-O-methyl-³H] (ARC; 80 Ci/mmol); incubate at  30 °C for 30 minutes
Binding should be performed in triplicate for each sample; measure non-specific binding by adding 10 µM non-radioactive benazine in the assay
- 16 Stop reaction by filtration through a Supor 200 0.2 µm filter (PALL); wash 3 times in ice-cold SHT buffer with 20 mM tetrabenezine (Fluka)
- 17 Measure radioligand signal and normalize specific binding to the amount of membrane protein added to the reaction