



Aug 03, 2023

## 🌐 Mouse genetic models to manipulate enterochromaffin cell activity - Murine Organoid ELISA

DOI

[dx.doi.org/10.17504/protocols.io.8epv5jz54l1b/v1](https://dx.doi.org/10.17504/protocols.io.8epv5jz54l1b/v1)



Koki Tohara<sup>1</sup>

<sup>1</sup>University of California, San Francisco



Koki Tohara

University of California, San Francisco

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5jz54l1b/v1>

**Protocol Citation:** Koki Tohara 2023. Mouse genetic models to manipulate enterochromaffin cell activity - Murine Organoid ELISA. [protocols.io](https://dx.doi.org/10.17504/protocols.io.8epv5jz54l1b/v1) <https://dx.doi.org/10.17504/protocols.io.8epv5jz54l1b/v1>

**Manuscript citation:**

Gut Enterochromaffin Cells are Critical Drivers of Visceral Pain and Anxiety  
(<https://www.biorxiv.org/content/10.1101/2022.04.04.486775v1>)



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 12, 2022

**Last Modified:** August 03, 2023

**Protocol Integer ID:** 71277

**Keywords:** Intestinal organoid, Serotonin, 5-HT, ELISA, murine organoid elisa, enterochromaffin cell activity, murine organoid elisa this protocol, murine intestinal organoid, intestinal organoid, serotonin, immunosorbent assay, linked immunosorbent assay, elisa, murine, enzyme

**Funders Acknowledgements:**

**Damon Runyon Cancer Research Foundation**

Grant ID: 2387-30

## Abstract

This protocol describes how we maintain murine intestinal organoids and perform the enzyme-linked immunosorbent assay (ELISA) for serotonin (5-HT).



## Materials

DPBS, no calcium, no magnesium (Thermo Fisher #14190144)

DPBS, no calcium, no magnesium + 2 mM EDTA

Advanced DMEM/F-12 (Thermo Fisher #12634010)

1 M HEPES (Thermo Fisher #15630080)

GlutaMAX™ Supplement (Thermo Fisher #35050061)

Penicillin-Streptomycin (10,000 U/mL) (Thermo Fisher #15140163)

Recombinant Murine Noggin (Peprotech #250-38)

Mouse EGF Recombinant Protein (Thermo Fisher #PMG8041)

*N*-Acetyl-L-cysteine (Millipore-Sigma #A7250-5G)

B-27™ Supplement (50X), serum-free (Thermo Fisher #17504044)

R-spondin 2 (supernatant from R-spondin expressing HEK293 cells)

Corning® Matrigel® Matrix (Corning #356231)

Bovine Serum Albumin (Millipore-Sigma #A9418)

Serotonin ELISA kit (LDN #BA-E-8900)

24-well plate (Falcon® 24-Well Plate #38021)

Corning™ Falcon™ Cell Strainers 70 µm (Corning #08-771-2)

Ringer's solution (140 mM NaCl, 5 mM KCl, 2 mM CaCl<sub>2</sub>, 2 mM MgCl<sub>2</sub>, 10 mM D-glucose, and 10 mM HEPES-Na [pH 7.4]))

Following four mouse lines were used.

RC::FL-hM3Dq (Jackson Labs, Strain #026942)

Tac1Cre (Jackson Labs, strain #021877)

Pet1Flp (gift of Dr. Susan Dymecki)

RC::PFTox (gift of Dr. Susan Dymecki)

## Troubleshooting

## Preparation or basal media

10m

- 1 Add 5 mL of GlutaMAX, 1 M HEPES, and Penicillin-Streptomycin to 500 mL of Advanced DMEM/F-12 and filter through 0.22  $\mu$ m. This basal media can last up to 1 month.

## Preparation of complete organoid culture media

- 2 Add 800  $\mu$ L of B-27, Recombinant Murine Noggin (final conc. 100 ng/mL), Mouse EGF Recombinant Protein (final conc. 500 ng/mL), *N*-Acetyl-L-cysteine (final conc. 0.5 mM), and 4 mL of R-spondin 2 to 40 mL of basal media. This complete media can last up to 2 weeks.

## Mouse organoid prep

10m

- 3 Organoids were generated from male Tac1-Cre(+/-);Pet1Flp(+/-);RC::FL-hM3Dq(+/-) or Tac1-Cre(+/-);Pet1Flp(+/-);RC::OFTox(+/-) animals (5-8 weeks old).
- 4 Euthanize the animal under CO<sub>2</sub> and spray belly with 70% EtOH.
- 5 Isolate the whole small intestine and divide the tissue into three pieces. Use the middle piece (jejunum) for organoid generation.
- 6 Filet open the intestine along the mesentery.
- 7 Scrape off the villi using a glass coverslip. Scrape well so that all the villi come off.
- 8 Cut the intestine into 2-4 mm pieces and move to a 50 mL falcon tube.
- 9 Add 10 mL ice-cold DPBS. Pipette up and down with 10 mL pipettes.
- 10 Wait for a couple of minutes until the tissue settles down and discard the supernatant.
- 11 Repeat the step 9-10 until the supernatant becomes clear (usually 4-5 times).



- 12 Add 30 mL of ice-cold DPBS + 2 mM EDTA and rock the tube at 4°C for 30 minutes.
- 13 While waiting, prepare 6× 50 mL falcon tubes and label them No.1-6.
- 14 Wait for a couple of minutes until the tissue settles down and discard the supernatant.
- 15 Add 10 mL ice-cold DPBS. Pipette up and down with 10 mL pipettes.
- 16 Wait for a couple of minutes until the tissue settles down and move the supernatant to falcon tube #1 prepared at step #13.
- 17 Repeat the step 15-16 five times. Note that the supernatant should be put into separated tubes prepared at step #13 (tube #2-6).
- 18 Place 20 µL of each supernatant fraction on a glass slide and check under a microscope.
- 19 Identify fractions that contain crypts.
- 20 Combine crypts-containing fractions and filter through a 70 µm strainer.
- 21 Spin down at 900 g for 5 minutes and discard the supernatant.
- 22 Resuspend pellet with 10 mL ice-cold basal media and spin down at 900 g for 5 minutes.
- 23 Carefully remove all the supernatant.
- 24 Repeat the step #22 and #23.



- 25 Resuspend the pellet into 400  $\mu$ L ice-cold matrigel.
- 26 Pipette 50  $\mu$ L of the matrigel/organoid suspension into the center of each of four wells of a prewarmed 24-well plate to form domes in the center of each well.
- 27 Place the lid on the culture plate and quickly turn the plate upside down.
- 28 Incubate at 37°C for 10 minutes to set the matrigel.
- 29 Turn the plate back upside so that you can apply media to wells. Gently add 500  $\mu$ L complete media to each well.
- 30 Add basal media to surrounding wells to prevent the wells from drying.

## Maintenance of intestinal organoids

46m

- 31 Put an aliquot of matrigel in fridge at 4°C. Put a 24-well plate to a 37°C incubator. 1m
- 32 Put desired volume of complete media at 37°C 1m
- 33 Aspirate media from wells 2m
- 34 Add 500  $\mu$ L cold basal media to a well and let it sit for about 15 seconds 2m
- 35 Pipette up and down to dissolve matrigel, collect organoids, and move to a 15 mL falcon tube 3m
- 36 Add 500  $\mu$ L cold basal media to a well, collect residual organoids, and move to the same 15 mL falcon tube 3m
- 37 Place a 200  $\mu$ L tip inside a 1000  $\mu$ L tip and pipette up and down ~25 times to break up organoids 3m



- 38 Add 10 mL cold basal media to the falcon tube and mix well 1m
- 39 Centrifuge at 200 g for 3 min at 4°C 3m
- 40 Aspirate the supernatant (be careful not to aspirate the organoid pellet) 2m
- 41 Add 10 mL cold basal media to the falcon tube and mix well 1m
- 42 Centrifuge at 200 g for 3 min at 4°C 3m
- 43 Aspirate the supernatant (be careful not to aspirate the organoid pellet) 1m
- 44 Add 200 µL of ice-cold matrigel to the falcon tube and mix the organoids well 2m
- 45 Pipette 50 µL of the matrigel/organoid suspension into the center of each of four wells of a prewarmed 24-well plate to form domes in the center of each well 3m
- 46 Place the lid on the culture plate and quickly turn the plate upside down 1m
- 47 Incubate at 37°C for 10 minutes to set the matrigel 10m
- 48 Turn the plate back upside so that you can apply media to wells. Gently add 500 µL complete media to each well 2m
- 49 Add basal media to surrounding wells to prevent the wells from drying 2m
- 50 Exchange the media every 3-4 days or whenever the media turns yellow.



51 Split again in 6-8 days.

## Preparation of organoids for ELISA

38m

52 Organoids are grown for 4-6 days

53 Aspirate the supernatant (be careful not to aspirate the organoid pellet)

2m

54 Add 500  $\mu$ L of cold DPBS to a well and let it sit for about 15 seconds

1m

55 Gently pipette up and down ~10 times with P-1000 to remove matrigel (but not to break the organoids up) and move to a 15 mL falcon tube.

2m

56 Add 500  $\mu$ L of cold DPBS to the tube

2m

57 Put the tube on ice for 5 min

5m

58 Gently pipette up and down ~10 times with P-1000

2m

59 Add 10 mL cold DPBS and spin down at 200 g for 3 min at 4°C

3m

60 Aspirate the supernatant as much as possible

1m

61 Add 1 mL cold DPBS and gently pipette up and down ~10 times with P-1000

2m

62 Add 10 mL cold DPBS and spin down at 200 g for 3 min at 4°C

3m

63 Aspirate the supernatant as much as possible

1m



- 
- 64 Repeat the step #33-35 twice to completely wash the matrigel out 12m
- 65 Add 100  $\mu$ L DPBS + 0.1% BSA and move organoids to an eppendorf tube (it is important to add 0.1% BSA to prevent organoids from sticking to tubes) 1m
- 66 Remove the supernatant and add 100  $\mu$ L Ringer's (140 mM NaCl, 5 mM KCl, 2 mM  $\text{CaCl}_2$ , 2 mM  $\text{MgCl}_2$ , 10 mM D-glucose, and 10 mM HEPES-Na [pH 7.4])) buffer for ELISA assay (do not lyse the organoids as you are measuring released serotonin). 1m

---

## 5-HT ELISA

- 67 ELISA is performed according to the manufacturer's protocol ([https://www.ldn.de/wp-content/uploads/IFU-BA-E-5900R-V15.0\\_wz.pdf](https://www.ldn.de/wp-content/uploads/IFU-BA-E-5900R-V15.0_wz.pdf)).