



Oct 23, 2024

Mouse PBMC Isolation

DOI

dx.doi.org/10.17504/protocols.io.j8nlk9bo6v5r/v1

Malu G Tansey¹

¹University of Florida



Janna Jernigan Posey

Indiana University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk9bo6v5r/v1>

Protocol Citation: Malu G Tansey 2024. Mouse PBMC Isolation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlk9bo6v5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 23, 2024



Last Modified: October 24, 2024

Protocol Integer ID: 110749

Keywords: ASAPCRN, mouse pbmc isolation, mouse pbmc isolation this protocol, isolating peripheral mononuclear blood cell, peripheral mononuclear blood cell, pbmc, isolation

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP) initiative

Grant ID: ASAP-020527

Disclaimer

The **protocols.io** team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

Abstract

This protocol describes the process of isolating peripheral mononuclear blood cells (PBMCs) from the blood of a mouse via ACK lysis.

Materials

Centrifuge that can be cooled to 4-degree Celsius

Isoflurane

Forceps

Hemostat

Scissors

1ml syringe

Monoject Blood Collection Tube Glycerine Coated Lavender Stopper, EDTA (K3) 0.06 mL 7.5% Solution, 3 mL Draw, 10.25 mm x 64 mm (8881311248)

15ml tubes

Ammonium-Chloride-Potassium (ACK) lysis buffer

Gilbico™ HBSS -calcium/-magnesium/-phenol red (Fisher, 14175103)

1.5ml tubes

Troubleshooting

Safety warnings

⚠ When working with Isoflurane work in a chemical safety hood. Wear appropriate PPE and use caution around sharps.

Before start

Remove Ammonium-Chloride-Potassium lysis buffer from fridge and fast temp centrifuge to 4 degrees celsius.

Blood Collection

- 1 Anesthetize mouse with isoflurane. Keep animal under isoflurane (we put kim wipe that have some isoflurane on them in a 50ml conical and keep the mouse's head in the 50ml conical) and open the abdomen at the sternum.
- 2 Open up the sternum and make some later cuts to expose all of the diaphragm. Be careful not to cut other organs.
- 3 Cut the diaphragm, being careful not to cut the lungs or heart behind it. Cut the lateral sides of the ribs cage and fold up to expose the heart. Keep the ribs folded back by clamping them down with a hemostat.
- 4 Cut the right atrium of the heart and collect blood that pools into the cavity with the 1ml syringe and immediately add to vacutainer with EDTA; flick tube or invert to mix and prevent clotting. Each mouse usually gives ~300ul of blood (range 100ul-700ul/mouse)

PBMC isolation

- 5 Transfer blood into a 15ml conical on ice.
 - a) For flow cytometry transfer 250uL of blood
 - b) For RNA extraction transfer all blood collected
- 6 Add 5mL of ACK lysis buffer to blood samples in 15ml tubes and vortex briefly
- 7 Incubate 5 minutes on ice. Red blood cell lysis is complete when the solution goes from opaque to translucent.
- 8 Quench lysis with 5mLs of 1XHBSS (-/-) bringing total volume up to 10mLs.
- 9 Centrifuge 15ml tubes at 4 degree for 5 minutes at 350xg.
- 10 After the 15mL tubes are done spinning, assess the pellet color. If the pellet is white no further lysis is needed.
- 11 If the pellet is red then you will need to repeat the ACK lysis step once more. To do so aspirate off supernatant and resuspend cell pellet with 1ml of ACK lysis and incubate for 1 minute on ice and then quench with 5mLs of HBSS (-/-). Centrifuge sample after incubation at 4 degree for 5 minutes at 350xg.



- 12 Carefully place samples on ice to not disrupt the pellet and aspirate off supernatant, can leave up to 300ul of supernatant to make sure the pellet is not disrupted
- 13 Add 5ml of 1X HBSS (-/-) to each tube and gently resuspend cells for wash (do not vortex, flick or hand pipet to resuspend)
- 14 Centrifuge samples 4 degrees Celsius, 5 minutes, at 350xg
- 15 Remove supernatant and resuspend cells according to downstream protocol.