

Mar 21, 2024

Mouse stereotaxic surgery

DOI

dx.doi.org/10.17504/protocols.io.81wgbx373lpk/v1

madalynn.erb Erb¹

¹Van Andel Institute

Team Lee



Jane Balster

ASAP - Team Lee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgbx373lpk/v1>

Protocol Citation: madalynn.erb Erb 2024. Mouse stereotaxic surgery. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.81wgbx373lpk/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 09, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 95133

Keywords: ASAPCRN, mouse stereotaxic surgery this protocol, mouse stereotaxic surgery, stereotaxic surgery, mouse

Funders Acknowledgements:

Madalynn Erb

Grant ID: ASAP000592

Disclaimer

The **protocols.io** team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

Abstract

This protocol details the mouse stereotaxic surgery.

Attachments



[1003-2591.docx](#)

28KB



Materials

Materials:


- Anesthesia : Inhalant Isoflurane
 1. Induction: 3.0 - 3.5%
 2. Maintenance: 2.0 – 2.5%
- Analgesia :
 1. Marcaine (Local)
 2. Ethiqx XR / buprenorphine extended-release (Systemic)
- Sterile 0.9% saline
- Sterile ophthalmic ointment
- Electric hair clippers
- 70% ethanol
- 3% Hydrogen peroxide
- Motorized stereotaxic frame
- Heating pad
- Sterile surgical instruments
- Sterile gauze and swabs
- Surgical drill
- Hamilton syringe
- Steel needle
- Surgical tubing
- Pump
- Surgical glue

Troubleshooting





Procedure


3d 0h 6m 30s

- 1 Turn on the heating pad to and set to  36.9 °C .
- 2 Weigh mouse and record weight before starting surgery.
- 3 Place mouse in isoflurane chamber (3% isoflurane).
- 3.1 Record time of anesthesia.
- 4 After mouse is fully anesthetized, use hair clippers to shave head.
- 5 Place mouse back into isoflurane chamber (3% isoflurane) until fully anesthetized.
- 6 Transfer mouse onto stereotaxic frame and place the mouth / nose into nose cone with isoflurane.
- 6.1 Turn down isoflurane to 2%.
- 7 Apply sterile ophthalmic ointment to eyes using a sterile swab (this will prevent desiccation).
- 8 Use ear bars to secure mouse in stereotaxic frame ensuring head is level in all directions.
- 9 Using a sterile swab apply 70% EtOH to the scalp.
- 10 Using a sterile swab, apply povidone-iodine solution to the scalp.
- 10.1 Wait for povidone-iodine solution to dry before making surgical incision.



11 Perform a subcutaneous injection of Ethiq^a XR (buprenorphine extended-release) into the leg.

11.1 Use  0.05 mL per  20 g mouse.

12 Inject  30 μ L of Marcaine in 2-3 locations underneath the scalp near the incision site. This is a local analgesic.

12.1 Wait 30 seconds -  00:01:00 for the Marcaine to diffuse before performing incision.

1m

13 Using a sterile scalpel, make a surgical incision to expose the skull.

13.1 Minimize the size of the incision as much as possible.



Note

You will need to see bregma and have access to the injection site. For substantia nigra injections this will be located at the caudal region of the skull.

14 Position your injection needle at bregma and save location in the AP and XY axis.

15 Enter and save injection coordinates into the motorized stereotaxic frame.

15.1 Coordinates for right substantia nigra

- Anterior-posterior (AP):  -2.9 mm
- Medial-lateral (ML):  -1.3 mm
- Dorso-ventral (DV): -4.2

16 Raise the needle away from the skull slightly and move to AP and ML injection coordinates.

17 Slowly lower the needle to touch the skull. Raise the needle and drill a small hole where the needle touched the skull.



17.1 Slowly drill through the skull keeping the drill shallow enough to not damage brain tissue.

18 Once the hole is drilled lower the needle to the surface of the brain ensuring that the needle is not deflected by the skull.

Note

The needle should be completely straight.

19 Raise the needle ± 30 mm - ± 40 mm providing space to flush the needle and load with virus.

19.1 Flush the needle with sterile H_2O at Room temperature.

- Cover the mouse's head with sterile gauze to absorb the H_2O .

19.2 Draw up a $1 \mu\text{L}$ air bubble.

19.3 Load needle with $2.5 \mu\text{L}$ of virus.

- Slowly draw up virus and watch the air / liquid interface to determine volume.
- Needle / tubing should be marked with $1 \mu\text{L}$ intervals using a sharpie to help with this step.

20 Move the prepared needle down to the surface of the brain. Move the needle to the desired DV coordinate at a slow speed ($250 \mu\text{m} / \text{sec}$).

21 Wait $00:00:30$ and then start the virus injections at $0.2 \mu\text{L} / \text{min}$.

30s




22 After the injection is finished wait for $00:05:00$, leaving the needle in place and allowing the virus to diffuse away from injection site.

5m

23 Remove the needle at slow speed ($250 \mu\text{m} / \text{sec}$).

24 Move the needle up away from the skull and then release it and rotate it out of the way.



- 25 Close the incision site using sterile forceps and surgical glue.
- 26 Administer  0.5 mL sterile saline solution subcutaneously to avoid dehydration. 
- 27 Place the mouse back into the home cage and monitor recovery.
- 27.1 Record time at recovery.
- 28 Monitor post operation recover for  72:00:00 and record any observations of pain or distress onto surgery cards. 