

Jul 10, 2024

Mouse synapse imaging and analysis

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn6mzpl5d/v1>

Shiyi Wang¹

¹Duke University



Shiyi Wang

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn6mzpl5d/v1>

Protocol Citation: Shiyi Wang 2024. Mouse synapse imaging and analysis. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.14egn6mzpl5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 10, 2024



Last Modified: August 13, 2024

Protocol Integer ID: 103199

Keywords: ASAPCRN, analysis mouse synapse imaging, mouse, imaging

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP) initiative

Grant ID: ASAP-020607

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

The **protocols.io** team notes that research involving animals and humans must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Ethics Committee or Board.

Abstract

Mouse synapse imaging and analysis



1 ****Tissue Sectioning and Preparation****

- 1.1 - Prepare coronal sections of 30 μ m thickness containing the anterior cingulate cortex (ACC) and primary motor cortex (MOp) from WT and LRRK2 G2019Ski/ki mice.

2 ****Synaptic Staining****

- 2.1 - Perform synaptic staining using the following antibody combinations:

- 2.2 - Excitatory (Intracortical): VGluT1 (pre-synaptic) and PSD95 (post-synaptic).

- 2.3 - Inhibitory: VGAT (pre-synaptic) and Gephyrin (post-synaptic).

3 ****Antibody Incubation****

- 3.1 - Dilute primary antibodies (anti-VGluT1, anti-PSD95, anti-VGAT, anti-Gephyrin) appropriately in blocking solution.

- 3.2 - Incubate tissue sections overnight at 4°C with primary antibodies.

4 ****Secondary Antibody Incubation****

- 4.1 - Wash sections with 1x TBS containing 0.2% Triton X-100 (TBST).

- 4.2 - Block with 10% normal goat serum (NGS) diluted in TBST.

- 4.3 - Incubate sections with Alexa Fluor-conjugated secondary antibodies (Life Technologies) for 2-3 hours at room temperature.

5 ****Confocal Microscopy****

- 5.1 - Acquire high-magnification images using an Olympus FV 3000 inverted confocal microscope.
- 5.2 - Use a 60x objective with 1.64x optical zoom.
- 5.3 - Acquire z-stack images consisting of 15 optical sections spaced 0.34 μm apart.

6 ****Image Processing****

- 6.1 - Convert each z-stack into 5 maximum projection images (MPI), each corresponding to a 1 μm section of the z-plane, using FIJI.

7 ****Synaptic Puncta Analysis****

- 7.1 - Analyze synaptic puncta using FIJI plugin Puncta Analyzer¹⁵³.
- 7.2 - Identify synapses by the colocalization of pre and postsynaptic puncta (VGluT1/PSD95 for excitatory synapses, VGAT/Gephyrin for inhibitory synapses).

8 ****Data Collection****

- 8.1 - Analyze 15 MPis per mouse (5 MPis per tissue section, 3 tissue sections per mouse).
- 8.2 - Analyze between 4 and 5 age- and sex-matched mice per genotype and condition, as specified in the figure legends.

9 ****Data Handling****

- 9.1 - Ensure all animals appear healthy at the time of collection.



9.2 - Include all data collected without exclusions, as per experimental design.

10 Notes:

11 - Maintain consistency in tissue processing and staining protocols.

12 - Perform all image analysis using standardized settings and procedures to ensure data reproducibility.