

Jul 21, 2025

Mouse tissue lysate preparation

DOI

<https://dx.doi.org/10.17504/protocols.io.14egnyd3zv5d/v1>

Hely Rodriguez Cruz¹, Michael Hanna¹, Pietro De Camilli¹

¹Yale University



Michael G Hanna

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egnyd3zv5d/v1>

Protocol Citation: Hely Rodriguez Cruz, Michael Hanna, Pietro De Camilli 2025. Mouse tissue lysate preparation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egnyd3zv5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 05, 2025



Last Modified: July 21, 2025

Protocol Integer ID: 219603

Keywords: ASAPCRN, mouse tissue lysate preparation this protocol, preparation of lysate, lysate preparation, mouse tissue, lysate, detergent, preparation, tissue

Abstract

This protocol is for the preparation of lysate of mouse tissues using a detergent in non-denaturing conditions

Materials

- Ice
- Mouse tissue
- 10-15mL dounce homogenizer
- MilliQ H₂O
- 15mL conical tubes
- 1.5mL microcentrifuge tubes
- TRIS
- SDS
- NaCl
- Protease Inhibitor Cocktail
- Heat Block

Mouse tissue lysate preparation

- 1 Weigh tissue, keep on ice, and prepare 1mL of Homogenization Buffer for approximately every 100mg of tissue.

Homogenization Buffer:

- 25mM TRIS pH 7.5
- 150mM NaCl
- 1x Protease Inhibitor Cocktail

- 2 Using a tissue homogenizer, submerge your tissue in the appropriate amount of homogenization buffer (1mL/100mg) in the homogenizing vessel. Homogenize the tissue thoroughly. Be sure to keep the tissue grinder submerged while homogenizing to avoid frothing.
- 3 Transfer the homogenate to a 15mL conical tube and add 0.5ul of benzonase per mL of homogenate.
- 4 Using a concentrated stock solution of SDS, bring the homogenate up to 2% SDS. Use the equation below to determine the volume of SDS to add from your stock solution to your homogenate in order to bring it to 2%. The volume of stock solution to add is as follows:

$$V_{\text{stock}} = V_{\text{start}} / ([\text{stock}]/[\text{target}] - 1)$$

where V_{stock} is the volume of stock solution to add, V_{start} is the total amount of homogenate you start with, [stock] is the concentration of your stock solution, and [target] is your final desired or target concentration (in this case 2%)

- 5 Take a small sample for BCA analysis (plan to dilute ~5x for BCA).
- 6 Add concentrated loading sample buffer and bring to 1x. Aliquot the lysate into microcentrifuge tubes. Boil tubes @ 95°C for 3 minutes.
- 7 Spin samples briefly at max speed for 10 seconds to spin down condensation left on the caps. Load samples immediately in a gel or store @ -20°C.