



Aug 22, 2024

Mouse Tissue mtDNA Copy Number Protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn6jwzl5d/v1>

Livia Hecke Moraes¹, Jack Devine²

¹California Institute of Technology; ²Columbia University Irving Medical Center



Jack Devine

Columbia University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn6jwzl5d/v1>

Protocol Citation: Livia Hecke Moraes, Jack Devine 2024. Mouse Tissue mtDNA Copy Number Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.14egn6jwzl5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 10, 2024

Last Modified: August 22, 2024

Protocol Integer ID: 105075

Keywords: mouse tissue mtdna copy number protocol, mouse

Abstract

Mouse Tissue mtDNA Copy Number Protocol developed in the Picard lab at Departments of Psychiatry and Neurology, Robert N Butler Columbia Aging Center, Columbia University Irving Medical Center, New York, NY, USA

Attachments



Tissue mtDNA Copy

Nu...

150KB

Consumables, with catalog numbers:

- 1
 - BD Vacutainer Tubes (BD #363083)
 - Tris HCl (Sigma #T3253)
 - Tween 20, 10% (Sigma #P1379)
 - Nuclease free water (ThermoFisher Scientific #AM9939)
 - Proteinase K (20 mg/mL) (ThermoFisher Scientific #AM2548)
 - Stericup Quick Release-GV Sterile Vacuum Filtration System – 500 mL (Millipore Sigma S2GVU05RE)
 - BrandTech 96-well semi-skirted plate (BrandTech #781375/VWR #10141-434)
 - BrandTech 8-strip domed tube caps (BrandTech #781340/VWR #80087-132)
 - 2x TaqMan Universal MasterMix Fast (LifeTech #4444965)
 - MicroAmp™ Optical 384-Well Reaction Plate with Barcode (ThermoFisher #4309849) MicroAmp™
 - Optical Adhesive Film (ThermoFisher #4311971) Validated
 - mtDNA/nDNA Primer & Probe sequences (idtDNA.com)
 - Additional consumables not listed here: 1.5 mL Eppendorf tubes, pipette tips as needed.

2

Lysis buffer for mtDNA measurement

- 3 *Keep everything sterile. The below calculations provide sufficient lysis buffer for 1 L of complete lysis buffer, or about 5,000 reactions.*
- 4 **Tris-HCl (0.38M, pH 8.5)** (Sigma #T3253)
 - 18.168 g of Tris HCl (MW = 157.60 g/mol) in 200 mL of nuclease-free dH₂O
 - pH to 8.5 with 5M NaOH (about 12.5 mL, but measure pH continuously)
 - Note: Accurate pH is critical for proper extraction. Make sure your pH meter has been recently calibrated.
 - Make up to 300 mL with nuclease-free dH₂O
 - Sterile filter Tris-HCl buffer using a SteriCup (*Make 6 mL or 12 mL aliquots, depending on if you want to run 1 or 2 plates at once*)
- 5 **Tween 20, 10%** (Sigma #P1379)
 - 540 mL of ultrapure dH₂O
 - 60mL of Tween 20 (wash tip thoroughly or leave in solution)
 - Sterile filter Tween 10% using a SteriCup (*Make 12 mL or 24 mL aliquots, depending on if you want to run 1 or 2 plates at once*)



NB: Thawing these 10% Tween aliquots can be slow. I recommend using a 37°C bead bath.

- 6 **ddH₂O (Sterile, nuclease-free)** (Thermofisher Scientific #AM9939)
(Make 5 mL aliquots or use direct from bottle)
- 7 **Proteinase K (20 mg/mL)** (Thermofisher Scientific #AM2548)
- To be put in fresh from 20 mg/mL (1:100 dilution, final 200 µg/mL)

Store aliquots at -30°C

8

	Recipe for 1ml Lysis Buffer
	300 µL – Tris HCl
	600 µL – Tween 20
	90 µL – nuclease free H ₂ O
	10 µL – Proteinase K

9

	Examples of Lysis Buffer Calculations (including a bit of excess)			
	Component	48 Samples	96 Samples (1 plate)	2 Plates
	Tris	3 mL	6 mL	12 mL
	Tween	6 mL	12 mL	24 mL
	Proteinase K	100 µL	200 µL	400 µL


	Nuclease-free water	900 μ L	1800 μ L	3.6 mL
	Total Volume	10 mL	20 mL	40 mL

10

Lysis for Mouse Tissue

- Thaw all reagents and prepare lysis buffer as described above.
- Add 180 μ L of lysis buffer in each well of a 96-well semi-skirted plate. We use BrandTech #781375. Supplied in US by VWR #10141-434.
- Add 20 μ L of sample into each well.
- We recommend one extra well with 200 μ L lysis buffer only to use as a “no template control”.
- Seal wells with 8-strip domed caps. We use BrandTech #781340. Supplied in US by VWR #80087-132.
- To ensure proper mixing of the sample and buffer, vortex vigorously.
- To pull the sample/buffer mixture to the bottom of the wells, quickly centrifuge at 1,000g for 10 sec before lysis.
- Proceed to heat-activated lysis of the samples in a Thermocycler for 16 hours at 55°C, followed by heat inactivation for 10 minutes at 95°C.

11 o



This digested sample can be used directly as template DNA in qPCR for cf-DNA measurements.

- Before performing qPCR, vortex and centrifuge again the template DNA-containing tube(s)/plate(s) at 1,000g for 10 sec at room temperature. If not using within 24 - 48h, freeze the DNA samples at -80°C or colder.



6. qPCR data analysis and error checking

- 12 Compute the mean, standard deviation, and coefficient of variation (CV) for the cycle thresholds (Ct) across each triplicate for each sample.
- 13 If the CV across the 3 triplicates is >2%, check the triplicates to see if any of the wells has a Ct that is >1 unit different from the other two, which could indicate a qPCR failure in that well.
- 14 In the case where amplification has failed in a well, remove the outlier.
- 15 The average Ct for each sample should be used compute copies/mL.

Appendix A: qPCR Reagents and MasterMix Recipe

16

	Recipe for each well of qPCR	
	TaqMan Universal MasterMix Fast	10 µL
	ND1 Primers F+R	+ 0.6 µL
	ND1 Probe	+ 0.4 µL
	B2M Primers F+R	+ 0.6 µL
	B2M Probe	+ 0.4 µL
	Total Reagent Volume	= 12 µL
	+ Sample Volume	+ 8 µL
	Total Reaction Volume	= 20 µL

The reagents used may depend on the qPCR equipment available, the genomic locations of interest, or other experimental parameters.

Below, we provide a summary of possible reagents as well as other validated primers and probes.

- 17 **2x TaqMan Universal MasterMix Fast** (Life Tech #4444965)
MicroAmp™ Optical 384-Well Reaction Plate with Barcode (ThermoFisher #4309849)
MicroAmp™ Optical Adhesive Film (ThermoFisher #4311971)
Validated mtDNA/nDNA Primer & Probe sequences (idtDNA.com)

18

qPCR Target	Sequences (5'→3')
Mouse mt-COX1-F	ACCACCATCATTCTCCTT CTC
Mouse mt-COX1- R	CTCCTGCATGGGCTAGATT T
Mouse mt-COX1-Probe:	HEX/AAGCAGGAG/ZEN/CA GGAACAGGATGAA/3IABkF Q
Mouse B2M-F	GAGAATGGGAAGCCGAAC ATA
Mouse B2M-	CCGTTCTTCAGCATTGGA TTT
Mouse B2M-Probe	FAM/CGTAACACA/ZEN/GT TCCACCCGCCTC/3IABkFQ

- 19
- Reconstitute primers in appropriate volume of nuclease free water to achieve 100 µM stock concentration.
 - Combine 120 µL each of 100 µM forward and reverse primers with 960 µL nuclease free water to achieve 10 µM working concentration.
 - Store primers at -30°C or -80°C until use.



- Reconstitute probe in appropriate volume of nuclease free water to achieve 100 μM stock concentration.
- Dilute 100 μM stock probe 20x to achieve 5 μM working concentration.
- Store probes at -30°C or -80°C until use (avoid freeze-thaw) and protect probes from light.