




Oct 14, 2025

Version 1

mtscMultiome ATAC + GEX V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgqq2w3gk5/v1>

Marc Appis^{1,2,3}, Yu-Hsin Hsieh^{1,2}, Pauline Kautz^{1,2,4}, Lena Nitsch^{1,2,3}, Maren Salla^{1,2,3}, Klara Penkert^{1,2}, Caleb Lareau⁵, Patrick Maschmeyer^{1,2}, Leif S. Ludwig^{1,2}

¹Berlin Institute of Health at Charité - Universitätsmedizin Berlin (BIH);

²Max-Delbrück-Center for Molecular Medicine in the Helmholtz Association (MDC), Berlin Institute for Medical Systems Biology (BIMSB);

³Freie Universität Berlin, Department of Biology/Chemistry/Pharmacy;

⁴Technische Universität Berlin, Institute of Biotechnology;

⁵Computational and Systems Biology Program, Memorial Sloan Kettering Cancer Center

**Marc L Appis**

Max Delbrück Center, Berlin Institute of Health at Charité

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)OPEN  ACCESS**DOI:** <https://dx.doi.org/10.17504/protocols.io.36wgqq2w3gk5/v1>



Protocol Citation: Marc Appis, Yu-Hsin Hsieh, Pauline Kautz, Lena Nitsch, Maren Salla, Klara Penkert, Caleb Lareau, Patrick Maschmeyer, Leif S. Ludwig 2025. mtscMultiome ATAC + GEX. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.36wgqg2w3gk5/v1>

Manuscript citation:

Mimitou EP, Lareau CA, Chen KY, Zorzetto-Fernandes AL, Hao Y, Takeshima Y, Luo W, Huang TS, Yeung BZ, Papalexi E, Thakore PI, Kibayashi T, Wing JB, Hata M, Satija R, Nazor KL, Sakaguchi S, Ludwig LS, Sankaran VG, Regev A, Smibert P. Scalable, multimodal profiling of chromatin accessibility, gene expression and protein levels in single cells. Nat Biotechnol. 2021 Oct;39(10):1246-1258. doi: 10.1038/s41587-021-00927-2. Epub 2021 Jun 3. PMID: 34083792; PMCID: PMC8763625.

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

This protocol is an implementation of the protocol published in Nature Biotechnology (see literature reference) and is successfully used in our lab on a regular basis. However, adaptations for application in another lab or for different target cell types might be necessary.

Created: June 26, 2025

Last Modified: October 14, 2025

Protocol Integer ID: 221118

Keywords: Single cell multiomics, Mitochondrial genomics, scATAC-seq, scRNA-seq, mitochondrial genome into the 10x genomics multiome atac, capture of mitochondrial dna, mitochondrial genome, mitochondrial dna, gex mtscmultiome atac, mtdna, 10x genomics multiome atac, mtscmultiome atac, epi multiome atac, cytometric enrichment, atac library generation, gex libraries ready for submission, human peripheral blood mononuclear cell, gex library, cryopreserved cell, desired cell population

Funders Acknowledgements:

National Human Genome Research Institute

Grant ID: UM1 HG012076

Disclaimer

The authors of this protocol don't take any legal responsibility for usage of this protocol, neither they give any guaranty for successful execution.

Abstract

mtscMultiome ATAC + GEX integrates the capture of mitochondrial DNA (mtDNA) for variant calling from the mitochondrial genome into the 10x Genomics Multiome ATAC + GEX workflow (now **Epi Multiome ATAC + Gene Expression**). The protocol was developed with human peripheral blood mononuclear cells but can be adapted for other cell types, specimens, and species, which may require further optimization.

The protocol starts from cryopreserved cells and includes (1) thawing, washing, and staining of surface markers for flow cytometric enrichment of desired cell populations, (2) cell fixation, permeabilization, and tagmentation, (3) GEM generation and barcoding, (4) post-GEM incubation cleanup, (5) pre-amplification PCR, (6) ATAC library generation, (7) cDNA amplification, and final (8) GEX library construction to yield ATAC and GEX libraries ready for submission to NGS.

The original version of this protocol was published by Mimitou et al., Nat Biotechnol. 2021 (doi: **10.1038/s41587-021-00927-2**).

Attachments



CG000338_ChromiumN

eX...

4.8MB

Image Attribution

Modified based on Nitsch L, Lareau CA, Ludwig LS. Mitochondrial genetics through the lens of single-cell multi-omics. Nat Genet. 2024 Jul;56(7):1355-1365. doi: 10.1038/s41588-024-01794-8. Epub 2024 Jul 1. PMID: 38951641; PMCID: PMC11260401.

Guidelines

This protocol follows the 10x Genomics Chromium Next GEM Single Cell Multiome ATAC + GEX workflow. Please consult the corresponding user guide before. The present protocol is based on the *Chromium Next GEM Single Cell Multiome ATAC + Gene Expression User Guide* • **Rev F** (see supplementary information to this protocol).

Starting from the tagmentation, important deviations from the 10x workflow are highlighted in this Protocols.io version. However, it is particularly useful to check reagent handling instructions, troubleshooting, and Chip loading with the Chromium Multiome user guide.

Please note that initial sample preparation and handling depend very much on the desired starting material and thus, sample handling until tagmentation might be adapted. This version of mtscMultiome ATAC + GEX was adapted to analyze human primary peripheral blood mononuclear cells (**PBMCs**).

Materials

Equipment required for mtscMultiome ATAC + GEX

- Standard laboratory equipment (microliter pipets and tips, multichannel pipets (ideally Rainin), tubes (PCR tubes, 1.5 mL, 2 mL, 15 mL, 50 mL), table-top centrifuges, Vortexer/shaker, heating block, serological pipets & tips)
- Vortexer with 10X adapter
- FACS tubes with filter (e.g., Cell Strainer Snap Cap)
- Tissue Culture Hood (clean bench)
- Swinging bucket centrifuge (with inlets for 50 mL, 15 mL, FACS tubes)
- Water bath (up to 37°C)
- Heat block (for microreaction tubes, up to 65°C)
- Cell culture microscope
- Cell counting equipment: either manual (i.e., Neubauer counting chamber) or automated (e.g., Countess 3 FL Automated Cell Counter)
- FACS sorter (if L/D and negative selection of neutrophils shall be performed; panel needs to be adapted based on the used FACS sorter)
- Chromium iX/X (Controller)
- 10x Magnetic Separator
- Dry ice
- BioAnalyzer (Agilent) (*recommended*) or TapeStation (Agilent) + corresponding hsDNA kit
- Qubit hsDNA Kit (Invitrogen)

Reagents required for mtscMultiome ATAC + GEX

1) For Day 1 (Steps 1 - 3)

	Reagent	Manufacturer / Vendor	Product ID	Link	Can be replaced?	Amount per sample
	Penicillin-Streptomycin (10,000 U/ml; 10 mg/ml)	PAN Biotech	#P06-07100	https://shop.pan-biotech.de/en/Penicillin-Streptomycin-10.000-U-ml-Penicillin-10-mg-ml-Streptomycin/	Optional supplement for thawing media, can be replaced by P/S from different manufacturer	500 µL



	Reagent	Manufacturer / Vendor	Product ID	Link	Can be replaced?	Amount per sample
				P06-07100		
	FBS Standard	PAN Biotech	#P30-3306	https://shop.pan-biotech.de/en/FBS-Standard-South-American-origin-fetal-bovine-serum-0.2-m-sterile-filtered/P30-3306	Necessary part of thawing media, can be replaced by FBS from different manufacturer	5 mL
	RPMI 1640 Glutamax	Thermo Fisher	#61870010	https://www.thermofisher.com/order/catalog/product/61870010	Necessary main part of the thawing media, can be replaced by RPMI from different manufacturer; Can be replaced by other full media (HPLM, DMEM etc.)	50 mL
	BSA stock	Miltenyi Biotec	#130-091-376	https://www.miltenyibiotec.com/DE-en/products/macsbbsa-stock-solution.html	Required for generation of protein-rich buffers; Can be replaced by BSA from different manufacturer	440 µL
	PBS	Thermo Fisher	#14190094	https://www.thermofisher.com/order/catalog/product/14190094	Can be replaced by sterile, Nuclease-free PBS from different manufacturer	13 mL



	Reagent	Manufacturer / Vendor	Product ID	Link	Can be replaced?	Amount per sample
	Protector RNase Inhibitor (40 U/μl)	Merck Millipore	#3335399001	https://www.sigmaaldrich.com/DE/de/product/roche/rnainhibitor	Could be replaced by RNase inhibitor from different manufacturer (with same activity/concentration); however, replacement is not recommended	31 μL
	NP-40 (10%)	Merck Millipore	#11332473001	https://www.fishersci.de/shop/products/nonidet-p40-substitute-ultrapure-thermo-scientific/15875388?srsltid=AfmBOorag-yw-6maozyGoIM59qTAWc5CKth054ntYJf8m1BFAnau9Xb5V	Can be replaced by NP-40 from different manufacturer	1.1 μL
	Dithiothreitol (DTT) (100 mM)	Thermo Fischer	#707265ML	https://www.thermofisher.com/order/catalog/product/707265ML	Can be replaced by DTT from different manufacturer	12.1 μL
	UltraPure Nuclease-free water	Thermo Fisher	#10977-035	https://www.thermofisher.com/order/catalog/	Can be replaced by nuclease-free water from different manufacturer	460 μL



	Reagent	Manufacturer / Vendor	Product ID	Link	Can be replaced?	Amount per sample
				log/product/10977035		
	Tris-HCl	SigmaAldrich	#T2194-100ML	https://www.sigmaaldrich.com/DE/de/product/sigma/t2194	Can be replaced by Tris_HCl from different manufacturer	
	NaCl (5M)	SigmaAldrich	#59222C-500ML	https://www.sigmaaldrich.com/DE/de/product/sigma/59222c	Can be replaced by NaCl from different manufacturer	
	MgCl₂ (1M)	SigmaAldrich	#M1028-100ML	https://www.sigmaaldrich.com/DE/de/product/sigma/m1028	Can be replaced by MgCl ₂ from different manufacturer	
	Trypan Blue	ThermoFischer	#T10282	https://www.thermofisher.com/order/catalog/product/T10282	Can be replaced by Trypan Blue from different manufacturer or by similar L/D staining dye	20-30 µL
	Pierce™ 16 % Formaldehyde (w/v), methanol-free	Thermo Fisher	#28908	https://www.thermofisher.com/order/catalog/product/28908	Can be replaced by PFA from different manufacturer; should have molecular biology grade and be methanol-free	3 µL
	Glycine	SigmaAldrich	#50046-50G	https://www.sigmaaldrich.com	Can be replaced with glycine (molecular biology grade) from a	25 µL of a 2.5 M solution



	Reagent	Manufacturer / Vendor	Product ID	Link	Can be replaced?	Amount per sample
				om/DE/de/product/sigma/50046	different manufacturer	
	Glycerol (50% v/v)	Ricca Chemical Company	#3290-32	https://www.riccachemical.com/products/general-use-reagents/other-aqueous-solutions/3290-(-r3290000-)	Can be replaced by Glycerol 50% (v/v) from a different manufacturer	~5 mL

Reagents required for day 1 of the mtscMultiome protocol

The following table includes the components of the recommended panel for FACS-based sorting of live CD66b-PBMCs. It can be adapted based on the target cell type and lab-specific requirements.

	Antigen	Colour / channel	Recommended product	Link	Notes	Amount per sample
	CD14	FITC (= AF488)	BioLegend AlexaFluor 488 anti-human CD14 (clone 63D3) #367130	https://www.biolegend.com/en-us/products/alexa-fluor-488-anti-human-cd14-antibody-15066	Can be replaced by PE anti-human CD14 from another manufacturer (new titration might be necessary)	1 μ L
	CD66b	PE	Biolegend PE anti-human CD66b /clone G10F5) #305106	https://www.biolegend.com/en-us/products/pe-anti-human-cd66b-	Can be replaced by PE anti-humanCD66b from another manufacturer (new titration might be necessary)	1 μ L

	Antigen	Colour / channel	Recommended product	Link	Notes	Amount per sample
				antibody-6529		
	FcX block	n/a	Human TruStain FcX (BioLegend, #422301)	https://www.biolegend.com/en-us/products/human-trustain-fcx-fc-receptor-blocking-solution-6462	Recommended to prevent non-specific binding of antibodies; can be replaced by other FcX-receptor block	2.5 µL
	L/D	Sytox Blue (DAPI)	SYTOX Blue dead cell stain (Invitrogen, #S34857)	https://www.thermofisher.com/order/catalog/product/S34857	Can be replaced by other L/D dye in the same channel	<1 µL

Recommended FACS panel for sorting of viable non-CD66b PBMCs.

The following table gives an overview on the components of the 10X Chromium Next GEM Single Cell Multiome ATAC + GEX kit required for day 1 of the mtscMultiome ATAC + GEX protocol. We strongly suggest to follow the handling instructions as provided in the kit manual. Furthermore, these items should not be replaced.

	Item	10x Product number	Kit Box	Amount per 1 sample + 10%
	20X Nuclei Buffer	2000207	ATAC Kit A	2.75 µL
	ATAC Buffer B	2000193	ATAC Kit A	7 µL
	ATAC Enzyme B	2000265	ATAC Kit A	3 µL
	Barcoding Reagent Mix	2000267	GEM Kit A	54.45 µL
	Barcoding Enzyme Mix	2000266	GEM Kit A	8.25 µL
	Template Switch Oligo	3000228	GEM Kit A (after resuspension to be)	1.21 µL



Item	10x Product number	Kit Box	Amount per 1 sample + 10%
		stored at -80°C)	
Reducing Agent B	2000087	GEM Kit A	2.1 µL
Quenching Agent	2000269	GEM Kit A	5 µL
Single Cell Multiome Gel Beads A	2000261	Gel Bead Kit A	50 µL
Chromium Next GEM Chip J	2000264	Next GEM Chip J & Gaskets	1 piece (up to 8 samples)
Gasket	370017	Next GEM Chip J & Gaskets	1 piece (1 per Chip J)
Partitioning Oil	2000190	Chromium Partitioning Oil	2-3 mL

Reagents from the 10x Chromium Next GEM Single Cell Multiome kit required for the first part of the mtscMultiome ATAC + GEX protocol (fixation, lysis, tagmentation, GEM generation, RT-PCR & quenching)

2) For day 2 (steps 4 - 5)

Reagent	Manufacturer	Product ID	Link	Can be replaced	Amount per sample
Dynabeads MyOne SILANE	10x Genomics	2000048	https://cdn.10xgenomics.com/image/upload/v1660261285/support-documents/2000048_Dynabeads_MyOne_SILANE_sds.pdf	Can't be replaced	14.3 µL
Elution buffer (EB)	Qiagen	19086	https://www.qiagen.com/de/products/discovery-and-translational-research/lab	Might be replacable. However, replacement by a different elution buffer was not tested and is therefore	220 µL



	Reagent	Manufacturer	Product ID	Link	Can be replaced	Amount per sample
				- essentials/buffers-reagents/buffer-eb?catno=19086	not recommended.	
	10% Tween-20				Can be replaced by 10% Tween-20 from a different manufacturer	<1 µL
	Ethanol (HPLC grade)	Carl Roth Chemie	9065.1	https://www.carlroth.com/com/en/alcohols/ethanol/p/9065.1	Can be replaced by analytical/HPLC-grade ethanol from a different manufacturer	2 mL
	UltraPure Nuclease-free water	ThermoFisher	10977-035	https://www.thermofisher.com/order/catalog/product/10977035	Can be replaced by nuclease-free water from different manufacturer	<2 mL
	SPRIselect DNA Size Selection Reagent	Beckman Coulter Life Sciences	B23317	https://www.beckman.de/reagents/genomic/cleanup-and-size-selection/size-selection/b23317	Can't be replaced	250 µL

Reagents required for day 2 of the mtscMultiome protocol

The following table gives an overview on the components of the 10X Chromium Next GEM Single Cell Multiome ATAC + GEX kit required for day 2 of the mtscMultiome ATAC + GEX protocol. We strongly suggest to follow the handling instructions as provided in the kit manual. Furthermore, these items should not be replaced.

	Item	10x Product number	Kit box	Amount per 1 sample
	Cleanup Buffer	2000088	GEM Kit A	200.2 µL
	Reducing Agent B	2000087	GEM Kit A	<1 µL
	Recovery Agent	220016	Chromium Recovery Agent	125 µL

	Item	10x Product number	Kit box	Amount per 1 sample
	Pre-Amp Mix	2000274	Amp Kit A	50 µL
	Pre-Amp Primers	2000271	Amp Kit A	4 µL

Reagents from the 10x Chromium Next GEM Single Cell Multiome kit required for the second part of the mtscMultiome ATAC + GEX protocol (Post GEM incubation cleanup and pre-amplification PCR & cleanup)

3) For day 3 (step 6)

For ATAC library construction, freshly prepared 80% ethanol ( 500 µL per sample), SPRIselect reagent ( 155 µL per sample), nuclease-free water ( 15 µL per sample) and Elution buffer ( 21 µL per sample) are required. For reagent details see table at day 2. Additionally, Qubit & BioAnalyzer reagents are required for library QC.

	Item	10x Product Number	Kit box	Amount per 1 sample
	Amp Mix	2000047	Amp Kit A	50 µL
	SI-PCR Primer B	2000128	Amp Kit A	7.5 µL
	Single Index Plate N Set A	3000427	Single Index Kit N Set A	2.5 µL

Reagents from the 10x Chromium Next GEM Single Cell Multiome kit required for the third part of the mtscMultiome ATAC + GEX protocol (ATAC library construction)

4) For day 4 (step 7)

For cDNA amplification, the following materials are required (see table at day 2 for more details on the reagents):

- SPRIselect reagent:  60 µL per sample
- Freshly prepared 80% ethanol:  400 µL per sample
- Elution buffer:  40.5 µL per sample
- BioAnalyzer reagents



	Item	10x Product number	Kit box	Amount per 1 sample
	Amp mix	2000047	Amp Kit A	55 µL
	cDNA Primers	2000089	Amp Kit A	16.5 µL

Reagents from the 10x Chromium Next GEM Single Cell Multiome kit required for the fourth part of the mtscMultiome ATAC + GEX protocol (cDNA amplification)

5) Step 8

This step is identical to the 10x Chromium Next GEM Single Cell Multiome ATAC + GEX kit step 7 (Gene Expression Library Construction). Please refer to the user guide for required materials.



Protocol materials

☒ SYTOX[®]; Blue Dead Cell Stain, for flow cytometry **Thermo Fisher Catalog #S34857**

☒ Human TruStain FcX[™] **BioLegend Catalog #422301**

☒ NP-40 10% **Merck MilliporeSigma (Sigma-Aldrich) Catalog #11332473001**

☒ MACS BSA Stock Solution **Miltenyi Biotec Catalog # 130-091-376**

☒ DTT (100mM Solution) **Thermo Fisher Scientific Catalog #707265ML**

☒ Protector RNase Inhibitor (40 U/ul) 2000 units **Merck MilliporeSigma (Sigma-Aldrich) Catalog #3335399001**

☒ Nuclease-free Water

☒ 1X PBS

☒ Pierce 16% Formaldehyde (w/v), Methanol-free **Thermo Fisher Scientific Catalog #28908**

☒ Dynabeads MyOne Silane **10x Genomics Catalog #2000048**

☒ Sodium Chloride Solution 5 M **Merck MilliporeSigma (Sigma-Aldrich) Catalog #59222C**

☒ RPMI 1640 Medium, GlutaMAX[®]; Supplement **Thermo Fisher Catalog #61870010**

☒ FBS Standard, South America origin, fetal bovine serum, 0.2 µm sterile filtered **PAN Biotech Catalog #P30-3306**

☒ Tris-HCl (pH7.4) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T2194-100ML**

☒ Glycine **Merck MilliporeSigma (Sigma-Aldrich) Catalog #50046**

☒ Penicillin-Streptomycin, 10.000 U/ml Penicillin, 10 mg/ml Streptomycin **PAN Biotech Catalog #P06-07100**

☒ DPBS, no calcium, no magnesium **Thermo Fisher Catalog #14190094**

☒ MgCl₂ **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M1028-100mL**

☒ BSA Stock Solution **Miltenyi Biotec Catalog #130-091-376**

Safety warnings

⚠ Please be careful when handling (potentially) **infectious biological material**. Stick to local regulations when working with biohazards.

This protocol includes the use of the following **hazardous goods**. Check safety instructions and hazardous goods sheets before starting to work with them (see the "Materials" section for links to the manufacturer websites):

- Paraformaldehyde
- Ethanol
- Penicillin/Streptomycin
- Dithiothreitol (DTT)
- NP-40

Ethics statement

Please follow all the ethical regulations for the use of human donor and patient samples. It is the responsibility of the user to ensure that all legal obligations are followed.

Before start

It is highly recommended to follow the checklist outlined below before starting with the experiment:

1. Book tissue culture hood, FACS sorter, centrifuges, Chromium X/iX, and PCR cyclers in advance
2. Check the availability of all materials and equipment as outlined in the "Materials" section of this protocol
3. At least one day before starting with sample handling, the buffers & reagents as listed in step 2.1 of this protocol should be prepared. They can be stored for multiple weeks to months at indicated temperatures.
4. The buffers as listed in steps 2.2-2.5 must be prepared freshly on the same day when the experiment is conducted.
5. 80% ethanol must be prepared freshly for steps 4-8.



Step 0: Preparations of materials and equipment

55m

1 Prepare cell thawing

15m

1. Cool down a swinging-bucket centrifuge to 4°C .
2. Have dry ice ready.
3. Pre-warm thawing media to 37°C .

2 Prepare buffers

2.1 General buffers and media

20m

1. The buffers and media in this step can be stored for several weeks at 4°C .

1) Prepare cell thawing media: RPMI + 10% FBS + 1% P/S

- Requires:

FBS Standard, South America origin, fetal bovine serum, 0.2 μm sterile filtered **PAN Biotech Catalog #P30-3306**

RPMI 1640 Medium, GlutaMAX[®]; Supplement **Thermo Fisher Catalog #61870010**

- Optional:

Penicillin-Streptomycin, 10.000 U/ml Penicillin, 10 mg/ml Streptomycin **PAN Biotech Catalog #P06-07100**

- Add 50 mL FBS (and optionally 5 mL P/S) to 500 mL RPMI
- 50 mL per sample required
- Filter media through a $0.2\text{ }\mu\text{m}$ filter
- Store at 4°C

2) Prepare recovery buffer for FACS: PBS + 0.04% BSA (PBS/BSA)

- Requires: BSA Stock Solution **Miltenyi Biotec Catalog #130-091-376**

DPBS, no calcium, no magnesium **Thermo Fisher Catalog #14190094**

- Add 2 mL 10% volume BSA to 500 mL PBS
- $5\text{--}7\text{ mL}$ per sample required
- Store at 4°C

3) Prepare 2X ATAC-Resuspension buffer (2X ATAC-RSB)



- Requires: Nuclease-free Water

Tris-HCl (pH7.4) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T2194-100ML**

Sodium Chloride Solution 5 M **Merck MilliporeSigma (Sigma-Aldrich) Catalog #59222C**

MgCl₂ **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M1028-100mL**

- Add per 116.46 mL Nuclease-free water: 720 µL 1 Molarity (M) MgCl₂,
 2399 µL 1 Molarity (M) Tris-HCl and 479 µL 5 Molarity (M) NaCl
- 620 µL per sample required
- Store at 4 °C

4) Prepare 2.5 M Glycine stock

- 2.5 Molarity (M) glycine (75.07 g/mol)
- Requires 1X PBS and
 Glycine **Merck MilliporeSigma (Sigma-Aldrich) Catalog #50046**
- Dissolve 4.7 g glycine in 25 mL 1X PBS
- Filtrate sterile (0.2 µm)
- 25 µL per sample required
- Store at Room temperature

Optional (in case a new Template Switch Oligo (TSO) aliquot must be prepared):

5) Prepare Template Switch Oligo (TSO)

Centrifuge a new tube of TSO (in the GEM Kit A, stored at -20 °C) briefly, resuspend in 80 µL Low TE Buffer. Vortex 00:00:15 at maximum speed, centrifuge briefly, leave at room temperature for a minimum of 00:30:00 . After resuspension, store at -80 °C .

2.2 PBS/BSA/RI (PBR) buffer

(Prepare freshly, keep on ice until use)

Requires:

- 1X PBS

5m

-  MACS BSA Stock Solution **Miltenyi Biotec Catalog # 130-091-376**
-  Protector RNase Inhibitor (40 U/ul) 2000 units **Merck MilliporeSigma (Sigma-Aldrich) Catalog #3335399001**

	Reagent	Stock concentration	Final concentration	per 1 sample + 10% (X * 1.1)
	PBS	-	-	2460.7 µL
	BSA	10%	1%	275 µL
	RNase inhibitor	40 U/µl	0.2 U/µl	14.3 µL
			Total	2750 µL

Pipetting scheme for the PBR buffer (required for washing of cells after sorting)

2.3 Lysis buffer (LLL buffer)

(Prepare freshly, keep on ice until use)

5m

Requires:

- 2X ATAC-RSB (prepared in step 2.1)
-  NP-40 10% **Merck MilliporeSigma (Sigma-Aldrich) Catalog #11332473001**
-  MACS BSA Stock Solution **Miltenyi Biotec Catalog # 130-091-376**
-  DTT (100mM Solution) **Thermo Fisher Scientific Catalog #707265ML**
-  Protector RNase Inhibitor (40 U/ul) 2000 units **Merck MilliporeSigma (Sigma-Aldrich) Catalog #3335399001**
-  Nuclease-free Water

	Reagent	Stock conc.	Final conc.	Per 1 sample + 10% (X * 1.1)
	2x ATAC-RSB	2x	1x	55 µL
	NP-40	10%	0.1%	1.1 µL
	BSA	10%	1%	11 µL
	DTT	100 mM	1 mM	1.1 µL
	RNase inhibitor	40 U/µL	1 U/µL	2.75 µL



	Reagent	Stock conc.	Final conc.	Per 1 sample + 10% (X * 1.1)
	Nuclease-free water	-	-	39.05 µL
			Total	110 µL

Pipetting scheme for the lysis buffer (LLL), mild lysis to retrieve intact but permeabilized mitochondria

2.4 Wash buffer (WB)

5m

(Prepare freshly, keep on ice until use)

Requires:

- 2X ATAC-RSB (prepared in step 2.1)
-  MACS BSA Stock Solution **Miltenyi Biotec Catalog # 130-091-376**
-  DTT (100mM Solution) **Thermo Fisher Scientific Catalog #707265ML**
-  Protector RNase Inhibitor (40 U/ul) 2000 units **Merck MilliporeSigma (Sigma-Aldrich) Catalog #3335399001**
-  Nuclease-free Water

	Reagent	Stock conc.	Final conc.	Per 1 sample + 10% (X * 1.1)
	2x ATAC-RSB	2x	1x	550 µL
	BSA	10%	1%	110 µL
	DTT	100 mM	1 mM	11 µL
	RNase inhibitor	40 U/µL	1 U/µL	13.75 µL
	Nuclease-free water	-	-	415.25 µL
			Total	1100 µL

Pipetting scheme for the wash buffer, used to wash cells post-sorting.

2.5 1x Nuclei buffer (1x NB)

5m

(Prepare freshly, keep  On ice until use)

Requires:

- 20X Nuclei Buffer (from the 10x Genomics Epi Multiome ATAC + Gene Expression Kit, 10x-Product number: 2000207, box "ATAC Kit A")
-  Protector RNase Inhibitor (40 U/ul) 2000 units **Merck MilliporeSigma (Sigma-Aldrich) Catalog #3335399001**
-  DTT (100mM Solution) **Thermo Fisher Scientific Catalog #707265ML**
-  Nuclease-free Water

	Reagent	Stock conc.	Final conc.	Per 1 sample
	20x Nuclei buffer	20x	1x	2.5 µL
	RNase inhibitor	40 U/µL	1 U/µL	1.3 µL
	DTT	100 mM	1 mM	0.5 µL
	Nuclease-free water	-	-	45.7 µL
			Total	50 µL

Pipetting scheme for the 1X Nuclei buffer, used for recovery of fixed, permeabilized cells for subsequent tagmentation and barcoding

Step 1: Cell Preparation - Thawing, washing, staining, FACS

3h 5m 30s

3 Thawing of cryopreserved PBMCs from liquid nitrogen stocks.

3.1 For cryopreserved cells, retrieve the cryovial from liquid nitrogen and transport it on dry ice. 

3.2 Quickly thaw in a water bath at  37 °C until only a small ice crystal remains. Remove from the water bath, spray with ethanol and transfer to a clean tissue culture hood. 3m

3.3 Carefully transfer the thawed cells from the cryovial to a 50 mL conical tube. 1m



3.4 Rinse the cryovial with  1 mL pre-warmed thawing media (RPMI + 10% FBS) and add it dropwise to the 50 mL conical tube while gently shaking the tube.

2m

3.5 Sequentially dilute cells in the 50 mL conical tube by incremental 1:1 volume additions of pre-warmed media at a speed of 1 mL/3-5 sec to the tube and swirl. Repeat for a total of 5 times with ~1 min pauses between additions.

7m

3.6 Centrifuge  400 x g, 4°C, 00:05:00

5m 30s



3.7 Carefully remove the supernatant and resuspend the cell pellet in 10 mL PBS or FACS buffer.

2m

3.8 Centrifuge  400 x g, 4°C, 00:05:00

6m



3.9 Carefully remove the supernatant and resuspend the cell pellet in 1 mL PBS or FACS buffer. Transfer the cell suspension to a new 15 mL tube and top up to a 5 or  10 mL total volume.

2m

3.10 Determine the cell concentration and viability using a

3m

Equipment

Countess 3 FL Automated Cell Counter

NAME

Automated Cell Counter

TYPE

ThermoFisher scientific

BRAND

AMQAF2000

SKU

<https://www.thermoFisher.com/th/en/home/life-science/cell-analysis/cell-analysis-instruments/automated-cell-counters/models/countess-3-fl.html>

LINK

or a hemocytometer. Write down cell counts and cell suspension volume.

3.11 Centrifuge  400 x g, 4°C, 00:05:00

6m



- 3.12 Remove the supernatant without disrupting the cell pellet and resuspend in  1 mL PBS/BSA. Transfer to a 1.5 mL microcentrifuge tube. Rinse the 15 mL tube with  0.4 mL PBS/BSA and transfer the rinse to the 1.5 mL tube containing the cells.

2m

4 FACS staining

- 4.1 Centrifuge  400 x g, 4°C, 00:05:00 and remove the supernatant.

6m



- 4.2 Resuspend cells in  100 μ L PBS/BSA (max. 5 million cells per  100 μ L).

1m

- 4.3 Add  2.5 μ L Fc blocking reagent

6m

 Human TruStain FcX™ BioLegend Catalog #422301 (= 1:40, V:V dilution) and

incubate for  00:05:00 at  4 °C .

- 4.4 Add  1 μ L CD66b-PE ( 1 % (v/v)),  1 μ L CD14-FITC ( 1 % (v/v)).

1m

Note

Optional co-staining of CD66b and CD14 for FACS:

Fluorescence-activated cell sorting (FACS) ensures reducing the fraction of dead cells and excluding residual neutrophils from the sample (e.g., after Ficoll), as both can negatively impact sequencing data quality. To remove neutrophils, the exclusion of CD66b⁺ cells is sufficient; however, because monocytes may express basal levels of CD66b, co-staining with CD14 is recommended to accurately exclude CD66b⁺CD14⁻ neutrophils while retaining all CD14⁺ monocytes.

- 4.5 Incubate for  00:20:00  On ice .

20m

- 4.6 Wash 2x with  1 mL PBS/BSA (spin  400 x g, 4°C, 00:05:00).

15m





4.7 Resuspend in  400 μ L PBS/BSA with  0.1 % (v/v)

2m



SYTOX[®] Blue Dead Cell Stain, for flow cytometry **Thermo**
Fisher Catalog #S34857

5 Sorting

5.1 Pass the cell suspension through a 5 mL FACS Tube with Cell Strainer Snap Cap to generate a single cell suspension right before sorting.

5m

5.2 Sort live CD66b⁺ cells into a 15 mL conical tube with  5 mL  On ice PBS + 0.04% BSA (prepared in step 1.2 (FACS recovery buffer) of the preparations section "step 0").

1h 30m

On a BD FACS Aria III, the following settings are recommended:

- Collect cells into a cooled 15 mL tube, rinsed with protein-rich buffer (see above; e.g., FACS buffer or PBS + BSA)
- ND filter: ND 1.5 filter
- Nozzle: 100 μ m

Note that those settings might be adapted if another FACS sorter is used.

Note

Optional in case of sensitive cells:

In case of sensitive cell types, the cells can also be sorted into media to reduce stress. To minimize adhesion of cells to the tube walls, the media should contain at least 0.4% BSA, FBS or any other protein source.

Expected result

Required cell number post-sort:

It is recommended to retrieve at least 500k viable cells (of the population of interest) post-sort to have sufficient starting material for fixation and lysis to end up with 40'000 cells to load onto the Chromium Chip.

Step 2: Fixation, Lysis/Permeabilization, and Tagmentation

2h 38m 45s



6 Fixation (~0.5hr)

- 6.1 Centrifuge the conical tubes with sorted cells for 400 x g, 4°C, 00:05:00 . 6m
- 6.2 Remove the supernatant without disrupting the pellet. 1m
- 6.3 Resuspend cells in 450 µL Room temperature PBS/BSA/RI (PBR buffer) (prepared in step 2.2 of the preparations section ("Step 0")). 1m
- 6.4 *(Optional, in case of not having sorted the cells):* Pass the cell suspension through a 5 mL FACS Tube with a Cell Strainer Snap Cap to generate a single cell suspension. *
- 6.5 Transfer to a 1.5 mL tube and immediately proceed with cell fixation. 30s
- 6.6 Add 3 µL 6m
 Pierce 16% Formaldehyde (w/v), Methanol-free **Thermo Fisher Scientific Catalog #28908**
on top of the cells in 450 µL PBS/BSA/RI, quickly invert the tube several times to fix the cells in 0.1% formaldehyde. Incubate for 00:05:00 at room temperature, swirl and invert occasionally.
- 6.7 Quench with glycine solution to a final concentration of 0.125 Mass Percent . Add 25 µL 2.5 Mass Percent glycine solution into the 453 µL fixed cells. Quickly invert the tube several times. 1m
- 6.8 Wash 2x with 1 mL On ice PBS/BSA/RI, mix by inverting the tube 5 times and centrifuge at 400 x g, 4°C, 00:05:00 . 14m
- 6.9 Remove the supernatant carefully and proceed with cell lysis/permeabilization.

7 Cell lysis, permeabilization, and intracellular staining

7.1 Add  100 μL ice-cold lysis buffer to the cell pellet, pipet gently up and down 3 times. Recommended cellular concentrations are 300-500k cells/100 μL . 1m

7.2 Incubate on ice for  00:03:00 (primary cells). 3m

7.3 Wash 1x by adding  1 mL of chilled wash buffer and mix by inverting the tube 3 times. **(Do NOT mix by pipetting up & down!)**. Centrifuge at  500 x g, 4°C, 00:05:00  **(not 400g!)**. 7m

7.4 Discard the supernatant gently without disrupting the pellet. 1m

7.5 Dilute cells in chilled 1x Diluted Nuclei buffer with a volume depending on the number of cells to be recovered after sequencing (see table below). Taking into account the expected loss of material following the lysis and wash steps, it is recommended to initially resuspend in a small volume of 10-20 μL . 10m

Depending on how many cells you expect, an initial recommendation is to resuspend the cells in  12.5 μL 1x Nuclei buffer from which  2.5 μL are taken for cell counting.

Number of cells recovered	Cell stock concentration (cells/ μL)
10,000	3,230 - 8,060

The final number of cells recovered at the end of the experimental workflow/detected in the sequencing data depends on the number of cells loaded. The mtscMultiome ATAC + GEX protocol loads 5 μL suspension of tagmented single cells onto the Chromium chip. Loading maximum 8k cells/ μL is recommended to obtain an optimal number of cells post-sequencing.

Note

Tip in case of poorly singulated cells:

In case cells aggregated during processing and didn't yield a well separated single cell suspension, filtering through a cell strainer snap cap (35 μm pore size) is recommended.

7.6 Count an aliquot using Trypan Blue and a cell counter, e.g. 10m



Equipment

Countess 3 FL Automated Cell Counter

NAME

Automated Cell Counter

TYPE

ThermoFisher scientific

BRAND

AMQAF2000

SKU

<https://www.thermoFisher.com/th/en/home/life-science/cell-analysis/cell-analysis-instruments/automated-cell-counters/models/countess-3-fl.html>

LINK

If necessary, adjust the cell count by diluting the cells in 1X Nuclei Buffer. Typically, we target **1M 8000 cells/μL**.

7.7 Proceed immediately with the tagmentation reaction.

8 Tagmentation

8.1 Prepare the Transposition Mix on ice. For multiple samples, prepare the Transposition Mix per sample in an individual PCR tube of an 8-tube strip. Pipette mix 10 times, centrifuge briefly and maintain on ice.

2m

	Transposition Mix	B	1x (μL)
	Total		10.0
	ATAC Enzyme		3.0
	ATAC buffer B		7.0

Composition of the Tagmentation/Transposition Mix

8.2 Gently mix the cell stock by pipetting it slowly up and down. Add the **5 μL** of the cell stock to the tube containing the Transposition Mix. Gently pipette mix 6 times (pipette set

5m

to 10 µL). **DO NOT centrifuge.**

8.3 Incubate in a thermal cycler using the following protocol:

Lid Temperature  50 °C , Reaction Volume  15 µL , Run Time  01:00:00 .

1h



	Step	Temperature	Time
	Incubate	37°C	00:60:00
	Hold	4°C	Hold

PCR cyclor program for the tagmentation.

8.4 Immediately proceed to the next step: 30 mins after the initiation of tagmentation, thaw reagents needed for GEM Generation and Barcoding.

8.5 Cut one tube of beads (Single Cell Multiome Gel Beads, PN. 2000261, -80°C) from the strip per sample, equilibrate the beads to room temperature. **Do not keep them on ice!**

8.6 Thaw the Template Switch Oligo (TSO) (stored at  -80 °C), Reducing Agent B, Barcoding Reagent Mix at room temperature (upon thawing, place  On ice).

Step 3: GEM Generation and Barcoding

2h 16m

9 Preparations

9.1 Prepare the GEM Generation Master Mix  On ice . Pipette mix 10x and centrifuge briefly.

5m

	GEM Generation Master Mix	Per 1 sample (X * 1.1) in µL
	Barcoding Reagent Mix	54.45
	Template Switch Oligo (TSO)	1.21
	Reducing Agent B	2.09
	Barcoding Enzyme Mix	8.25
	Total	66

Pipetting scheme for the GEM Generation and Barcoding Master Mix

9.2 Prepare chip loading

15m

Place  On ice :

1. Samples after tagmentation
2. Master Mix
3. Pre-chilled 8-strip emulsion-safe PCR tubes (1 tube per sample).

 Room temperature :

1. Single Cell Multiome Gel Beads
2. Chromium Next GEM Chip J and 10x Gasket
3. Partitioning oil and 50% glycerol
4. Reagent reservoir with divider
5. Two lids of tip boxes (to cover the chip and reservoir while loading)
6. 200 μ L emulsion-safe Low-Retention Filter tip (e.g., Rainin Filter Tip)
7. 1 Emulsion-safe PCR tubes for the remaining beads.

10 Chip assembly, loading of the Chromium Controller, and transfer of GEMs

35m



- 10.1 ■ Follow the instructions as outlined in the 10x Genomics Epi Multiome ATAC + Gene Expression user guide, CG000338 Rev E: Step 2.1-2.4.

10.2 General remarks

1. Use a Rainin pipette
2. While pouring 50% glycerol and partitioning oil into the two compartments of the reagent reservoir, do not mix the two reagents!
3. Before taking the chip and gasket, spread your gloves with ethanol to prevent static electricity and only touch the side of chip/gasket!
4. Cover the chip and reagent reservoir with the lid of the tip box between steps.
5. Put the gel bead holder onto the 10x Genomics vortex adaptor before loading your sample.

- 10.3 Assemble the Chromium Next GEM Chip J according to the manufacturer's instructions.

- 10.4 Load the Chromium Next GEM Chip J according to the manufacturer's instructions.

- 10.5 Run the Chromium Controller according to the manufacturer's instructions.

10.6 Transfer GEMs from the chip to the PCR strip  On ice according to the manufacturer's instructions.

11 GEM Incubation

11.1 Incubate in a Bio-Rad C1000 Touch Thermal cycler with a 96-Deep Well Reaction Module with the following protocol:

1h 15m

Cycler settings:

- Lid temperature  50 °C
- Reaction volume  100 µL
- Total run time  01:15:00
- Program:

	Step	Temperature	Time
	1	37°C	00:45:00
	2	25°C	00:30:00
	3	4°C	Hold (do not leave over night!)

PCR cycler program for GEM incubation (RT-PCR)

Equipment

C1000 Touch™ Thermal Cycler with 96-Deep Well Reaction Module	NAME
Thermal Cycler	TYPE
BioRad	BRAND
1851197	SKU
https://www.bio-rad.com/en-us/sku/1851197-c1000-touch-thermal-cycler-with-96-ndash-deep-well-reaction-module?ID=1851197	LINK

11.2 Note: **After GEM incubation, proceed immediately to Quench!**

11.3 Retrieve the Quenching Agent (PN-2000269) from  -20 °C and equilibrate to  Room temperature while the PCR program is running.

11.4 Add  5 µL of the Quenching Agent to each sample to stop the reaction.

1m

11.5 Slowly pipette mix 10 times (pipette set to 90 µl). The solution will be viscous. Ensure that no liquid remains along the tube walls and pipette tips. If necessary, aspirate the entire volume and dispense slowly back into the tube.

5m

11.6 Store at  -80 °C for up to 4 weeks, or proceed to the next step.

II

Step 4: Post GEM Incubation Cleanup

1h 11m

12 Follow the instructions as outlined in the Chromium Next GEM Single Cell Multiome ATAC + Gene Expression user guide, CG000338 Rev F: Step 3.

13 Material and pre-prep

13.1 Thaw the Cleanup Buffer at  65 °C for  00:10:00 at maximum speed on a Thermomixer with shaking feature. Afterwards, cool to  Room temperature .

10m



13.2 Thaw the Reducing Agent B to  Room temperature (Upon thawing, place  On ice).

13.3 Prepare the Dynabeads Cleanup Mix

5m

- Vortex
- Maintain at  Room temperature .
- Requires  Dynabeads MyOne Silane **10x Genomics Catalog #2000048** *

Component	Per 1 sample + 10% (X * 1.1)
Cleanup Buffer	200.2 µL

Component	Per 1 sample + 10% (X * 1.1)
Dynabeads MyOne SILANE*	14.3 µL
Reducing Agent B	5.5 µL
Total	220 µL

Pipetting scheme for the Dynabeads Cleanup Mix.

* Store at  4 °C . Vortex thoroughly (≥30 sec) before adding the Dynabeads to the mix. Aspirate the full liquid volume with a pipette tip to verify that the beads have not settled in the bottom of the tube. If clumps are present, pipette mix to resuspend completely. **DO NOT centrifuge before use.**

13.4 Prepare the Elution Solution I

5m

- Add reagents in the order listed
- Vortex, spin down
- Maintain at  Room temperature

Component	Per 1 sample + 10% (X * 1.1)
Buffer EB	59.25 µL
10% Tween 20	0.605 µL
Reducing Agent B	0.605 µL
Total	60.5 µL

Pipetting scheme for the Elution Solution I

13.5 Prepare  80 % (v/v) ethanol ( 10 mL for 5 reactions) freshly.

2m

14 Dynabeads cleanup



14.1 Add  125 µL Recovery Agent to each sample at room temperature. **DO NOT pipette mix or vortex the biphasic mixture!**

1m

14.2 Gently invert the tube 10 times to mix. *Keep the lid actively closed with the fingers to prevent unintentional leaking of the tube during inversion!*

1m



Centrifuge briefly.

- 14.3 Slowly remove and discard  125 μL Recovery Agent/Partitioning Oil (pink) from the bottom of each tube without aspirating any aqueous sample. 2m
- 14.4 Vortex Dynabeads Cleanup Mix and add  200 μL to each sample. Pipette mix 5x with the pipette set to 200 μL . **Do NOT close the cap of the PCR tube**, otherwise the sample will spill out!  2m
- 14.5 Cover the PCR tube with the lid of a tip box and incubate for  00:10:00 at  Room temperature  10m
- 14.6 At the end of 10 min incubation, place on the 10x Magnetic Separator, high position (magnet•**High**) until the solution clears. 2m
- 14.7 Discard the supernatant. 1m
- 14.8 While on the magnet, add  300 μL freshly prepared  80 % (v/v) ethanol to each sample without disrupting the pellet. 30s
- 14.9 Wait  00:00:30 and then remove the ethanol. 1m
- 14.10 While the sample is on the magnet, add  200 μL  80 % (v/v) ethanol to each sample without disrupting the pellet. 30s
- 14.11 Wait  00:00:30 and then remove the ethanol. 30s
- 14.12 Remove tubes from the magnet and centrifuge briefly. Place again on the magnet•**Low**. 1m
- 14.13 Remove any remaining ethanol. 30s
- 14.14 Remove the tube strip from the magnet. Immediately add  50 μL Elution Solution I to each sample. 1m
- 14.15 Pipette mix (pipette set to 30 μL) without introducing bubbles.



14.16 Incubate  00:01:00 at room temperature.

1m

1m

14.17 Centrifuge briefly. Place on the magnet•**Low** until the solution clears.

1m

14.18 Transfer  50 μL sample to a new tube strip and proceed to the SPRI cleanup as per protocol.

1m

15 **SPRIselect cleanup**

15.1 Vortex the SPRIselect reagent until fully resuspended. Add  90 μL (1.8X) SPRIselect reagent to each sample. Pipette thoroughly.

1m

15.2 Incubate  00:05:00 at  Room temperature .

5m

15.3 Centrifuge briefly. Place on the magnet•**High** until the solution clears.

2m

15.4 Discard the supernatant.

1m

15.5 While keeping the tube(s) on the magnet•**High**, add  200 μL  80 % (v/v) ethanol to each sample without disrupting the pellet. Wait  00:00:30 and then remove the ethanol.

1m

15.6 Repeat the wash for a total of 2 washes.

3m

15.7 Centrifuge briefly and place on the magnet•**Low**.

1m

15.8 Remove any remaining ethanol.

1m

- 15.9 Remove the tube strip from the magnet and immediately add  46.5 µL Buffer EB. 1m
- 15.10 Pipette mix (pipette set to 45 µL) without introducing bubbles. 1m
- 15.11 Incubate for  00:02:00 at  Room temperature . 2m
- 15.12 Centrifuge briefly and place on the magnet. **Low** until the solution clears. 1m
- 15.13 Transfer  46 µL sample to a new tube strip. 1m

Step 5: Pre-Amplification PCR & Clean-up

53m 30s

16 Material and pre-prep

- 16.1 Equilibrate the Pre-Amp Primers to  Room temperature . Upon thawing, vortex and centrifuge briefly. Maintain  On ice .
- 16.2 Thaw the Pre-Amp Mix on ice. Pipette mix and centrifuge briefly.

17 ATAC/GEX Pre-amplification PCR

(as outlined in the Chromium Next GEM Single Cell Multiome ATAC + Gene Expression user guide, CG000338 Rev F: Step 4)

- 17.1 Prepare the Pre-Amplification Mix  On ice . Pipette mix 10X and centrifuge briefly. 3m

	Component	10X product number (PN)	Per sample (µL)
	Pre-Amp Mix	2000274	50
	Pre-Amp Primers	2000271	4

Pipetting scheme for the **Pre-Amplification Mix**.

17.2 Add sample ( 46 µL from step 15.13) to the Pre-Amplification Mix. Pipette mix and centrifuge briefly.

2m

17.3 Incubate in a Bio-Rad C1000 Touch Thermal cycler with a 96-Deep Well Reaction Module with the following protocol.

30m



- Lid temperature:  105 °C
- Reaction volume:  100 µL
- Total run time:  00:30:00

	Step	Temperature	Run time
	1	72°C	00:05:00
	2	98°C	00:03:00
	3	98°C	00:00:20
	4	63°C	00:00:30
	5	72°C	00:01:00
	6	Go to step 3, repeat 6x (7 cycles total)	
	7	4°C	Hold

PCR program for pre-amplification PCR

17.4 **[STOP POINT]** Store at  4 °C for up to 18 h or proceed to the next step.



18 Post Pre-Amplification SPRI cleanup

! Important: Differences to the 10x Genomics protocol !

(Different volumes will be used for library elution and splitting.)



18.1 Vortex to resuspend SPRIselect reagent. Add  160 µL SPRIselect reagent to each sample. Pipette mix 15 times. (pipette set to 200 ul).

2m

18.2 Incubate  00:05:00 at room temperature.

5m

18.3 Centrifuge briefly. Place on the magnet • **High** until the solution clears.

2m



- 18.4 Discard the supernatant. 30s
- 18.5 While sample on the magnet, add  300 μL freshly prepared [M] 80 % (v/v) ethanol without disrupting the pellet. 30s
- 18.6 Wait  00:00:30 and then remove the ethanol. 30s
- 18.7 While on the magnet, add  200 μL [M] 80 % (v/v) ethanol to the sample without disrupting the pellet. 1m
- 18.8 Wait 30 sec and then remove the ethanol. 1m
- 18.9 Centrifuge briefly. Place on the magnet • **Low**. 30s
- 18.10 Remove any remaining ethanol. 30s
- 18.11 **! Difference to the 10x Genomics protocol !**
Remove the tube strip from the magnet. Immediately add  100.5 μL Buffer EB. 30s 
- 18.12 Pipette mix (pipette set to 90 μl) without introducing bubbles. 1m
- 18.13 Incubate  00:02:00 at  Room temperature . 2m
- 18.14 Centrifuge briefly. Place on the magnet • **High** until the solution clears. 1m
- 18.15 **! Difference to the 10x Genomics protocol !**
Transfer  100 μL sample to a new tube strip and proceed to library construction as per protocol. The 100 μL pre-amplified, SPRI-cleaned sample includes barcoded transposed DNA fragments and barcoded cDNA. 30s 
- 18.16 **! Difference to the 10X protocol !** 

The sample is divided and used as input for three separate steps:

1.  25 μL sample (25%) is used for ATAC Library Construction (step 6)
2.  35 μL sample (35%) is used for cDNA Amplification (step 7)
3. Store the remaining Pre-amplification product at  $-20\text{ }^{\circ}\text{C}$ for long term.

18.17 **[STOP POINT]** Store at  $4\text{ }^{\circ}\text{C}$ for up to 72 h or  $-20\text{ }^{\circ}\text{C}$ for long-term storage.



Step 6: ATAC Library Construction

1h 16m

19 As outlined in the Chromium Next GEM Single Cell Multiome ATAC + Gene Expression user guide, CG000338 Rev E: Step 5.

! Important: Differences to the 10x Genomics protocol !

20 Material and pre-prep

20.1 Equilibrate the Single Index Plate N set A to room temperature. Upon thawing, centrifuge at  1000 x g, 4°C , 00:01:00 . Maintain  On ice .

20.2 Thaw the SI-PCR primer B  On ice . Vortex and centrifuge briefly.

20.3 Thaw the Amp Mix  On ice . Pipette mix and centrifuge briefly.

21 ATAC-Sample Index (SI) PCR

21.1 ! Difference to the 10x Genomics protocol !

Prepare the Sample Index PCR Mix  On ice .

4m



Note: due to the reduced resuspension volume in the previous step, water has to be added to the mix to fill the input for the ATAC Index PCR mix ad $40\text{ }\mu\text{L}$ total volume ($25\text{ }\mu\text{L}$ library + $15\text{ }\mu\text{L}$ water).

	Component	10X product number (PN)	1 sample (μL)
	Amp Mix	2000047	50

Component	10X product number (PN)	1 sample (μL)
SI-PCR Primer B	2000128	7.5
Nuclease-free water	-	15
Total		72.5

Pipetting scheme for the Sample Index PCR Mix

21.2 ! Difference to the 10x Genomics protocol !

Add  25 μL pre-amplified sample (from 5.2.16.a) to 72.5 μL the Sample Index PCR Mix. Pipette mix and centrifuge briefly.

21.3 Add  2.5 μL of an individual Single Index N Set A to each well. Record the barcode/well that was used for each sample. Pipette mix and centrifuge briefly.

21.4 Incubate in a Bio-Rad cycler with 96-Deep Well Reaction Module with the following protocol.
(For mtscATAC/mtscMultiome ATAC library construction, we typically add 1-2 PCR cycles compared to what is recommended by 10x Genomics for their standard single-cell ATAC-seq protocol.)

- Lid temperature:  105 °C
- Reaction volume:  100 μL
- Run time: ~

Step	Temperature	Time
1	98°C	00:00:45
2	98°C	00:00:20
3	67°C	00:00:30
4	72°C	00:00:20
5	Go to step 2, see table below for number of cycles	
6	72°C	00:01:00
7	4°C	Hold

PCR protocol for the ATAC sample index PCR

2m

5m



40m



	Number of loaded cells	Recommended total PCR cycles
	<2,000	10
	2,001 - 6,000	9
	6,001 - 10,000	8

Cycle number optimization table. The number of PCR cycles needs to be adapted based on the number of cells per μL loaded on the Chip J as determined in step 7.6 of this protocol (step 1.0 in the 10x Genomics protocol).

21.5 **[STOP POINT]** Store at  4 °C for up to 72 h or proceed to the next step.



22 Post Sample Index Double Sided Size Selection – SPRIselect

22.1 Vortex to resuspend the SPRIselect reagent. Add  60 μL (0.6X) SPRIselect reagent to each sample. Pipette mix 15 times. (pipette set to 150 μL).

1m 30s

22.2 Incubate  00:05:00 at  Room temperature .

5m

22.3 Place on the magnet•**High** until the solution clears.

1m

22.4 Transfer  150 μL supernatant to a new strip tube. **DO NOT discard the supernatant.**

1m



22.5 Vortex to resuspend the SPRIselect reagent. Add  95 μL SPRIselect reagent (1.55X) to each sample. Pipette mix 15 times. (Pipette set to 180 μL).

1m

22.6 Incubate  00:05:00 at  Room temperature .

5m

22.7 Place on the magnet•**High** until the solution clears.

1m



- 22.8 Discard the supernatant. 30s
- 22.9 While the sample is on the magnet, add  300 μL freshly prepared  80 % (v/v) ethanol without disrupting the pellet. 30s
- 22.10 Wait 30 sec and then remove the ethanol. 1m
- 22.11 While the sample is on the magnet, add  200 μL  80 % (v/v) ethanol without disrupting the pellet. 1m
- 22.12 Wait 30 sec and then remove the ethanol. 1m
- 22.13 Centrifuge briefly and place on the magnet • **Low**. 30s
- 22.14 Remove remaining ethanol. 30s
- 22.15 Remove from the magnet and **immediately** add  20.5 μL Buffer EB. Pipette mix. 1m
- 22.16 Incubate  00:02:00 at  Room temperature . 2m
- 22.17 Centrifuge briefly. Place on the magnet • **Low** until the solution clears. 1m
- 22.18 Transfer  20 μL sample to a new tube strip. 30s
- 22.19 **[STOP POINT]** Store at  4 °C for up to 72 h or  -20 °C for long-term storage. II
- 23 **ATAC Post Library QC & Quantification**




Follow the instruction as described in the Chromium Next GEM Single Cell Multiome ATAC + Gene Expression user guide, CG000338 Rev F: Step 5.3.

Expected result

For the scATAC libraries from a mtscMultiome ATAC + GEX run, we typically expect concentrations between [M] 40-100 ng/ μ L and molarities of [M] 100-400 nanomolar (nM)

This is how the BioA profile of a typical scATAC library from a mtscMultiome run should look like:



Example of an expected DNA fragment size profile of an ATAC library from a mtscMultiome ATAC + GEX run recorded with a BioAnalyzer (Agilent)

Step 7: cDNA Amplification

3h 15m 30s

- 24 Follow the steps as outlined in Chromium Next GEM Single Cell Multiome ATAC + Gene Expression user guide, CG000338 Rev E: Step 6.

25 cDNA Amplification

50m

For the cDNA amplification PCR, use the following PCR program and PCR cycles.



	Step	Temperature	Time
	1	98°C	00:03:00
	2	98°C	00:00:15
	3	63°C	00:00:20
	4	72°C	00:01:00
	5	Go to step 2, see table below for the number of total cycles	
	6	72°C	00:01:00
	7	4°C	Hold

PCR program for cDNA amplification.

	Nr. of loaded cells per μL	Total PCR cycles
	>2,0001	9
	2,001-6,000	7
	6,001-10,000	6

Cycle number optimization table. The number of PCR cycles needs to be adapted based on the number of cells per μL loaded on the Chip J as determined in step 7.6 of this protocol (step 1.0 of the 10x Genomics protocol). This table is identical to the table "Recommended starting point for cycle number optimization" in step 6.1 of the 10x Genomics protocol.

[STOP POINT]

After the PCR, amplified cDNAs can be stored at  4 °C for up to 72 h or at  -20 °C for up to 7 days or proceed to the next step.

26 cDNA Cleanup – SPRIselect

[STOP POINT] After the clean-up, the libraries can be stored at  4 °C for up to 72 h or at  -20 °C for up to a month.

27 cDNA QC & Quantification

2h

Follow the instructions as described in Chromium Next GEM Single Cell Multiome ATAC + Gene Expression user guide, CG000338 Rev F: Step 6.3. **DNA concentration quantification via Qubit is not required**, 1 μL undiluted amplified cDNA can directly be



used for the determination of peak sizes & distribution via BioAnalyzer. However, quantification via Qubit can **is recommended** to obtain a quick estimation of yield and to ensure that no dilution is required before loading the sample onto the BioA chip.

Step 8: GEX Library Construction

6h

- 28 Follow the instruction as described in Chromium Next GEM Single Cell Multiome ATAC + Gene Expression user guide, CG000338 Rev F: Step 7.

6h

Acknowledgements

We thank Paul Mendez and Moritz Kochanek for critically reviewing the draft document of this protocol.