


Apr 22, 2020

Version 1

Multi-Seq: my notes from the lab V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.bejyjcpw>



Luciano Martelotto¹

¹University of Adelaide



Luciano G Martelotto

Waters Corporation (former BD Biosciences)

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bejyjcpw>

Protocol Citation: Luciano Martelotto 2020. Multi-Seq: my notes from the lab. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bejyjcpw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol in our group and it is working. Optimization is needed.

Created: April 03, 2020

Last Modified: April 22, 2020

Protocol Integer ID: 35160

Keywords: notes from the lab, note, lab

Abstract

Please note these are just my notes and by no means are ground truth of how Multi-seq needs to be performed. I WELCOME to comments, suggestions and amendments.



Guidelines

Notes are for nuclei/CMO; however, for cells/LMO the protocol is the basically the same but using 0.04% BSA because BSA sequesters LMOs (See Multi-Seq paper, McGinnis *et al.*, Nat Med, 2019).

You may scale up or down LMO/CMO concentrations while maintaining the same ratio LMO/CMO per cell/nuclei used in paper.

-For LMOs/CMOs on live cells that's 8 pmol/100K or 40 pmol/500K, in a reaction volume of 200 uL, which gives you 40 nM (100K cells) and 200 nM (500K cell), respectively.

-For LMOs/CMOs on nuclei they seem to have used 500 uM (500K nuclei), that is 100 pmol/500K, in (likely) a reaction volume of 200 uL, so it'll be 20 pmol/100K nuclei (this is 2.5 times higher than for lives cells). We have used CMOs at 200 nM in a 200 uL reaction and works fine. Both version will be given in the step-by-step protocol, but the one inferred from paper (using 500 nM) will be *[in between brakets and italicized]*.

-I would recommend to determine what's the optimal amount of CMO or LMO/number of cells or nuclei in your hands. For GM12878 cell line I used 200 nM for nuclei and managed to get very good libraries. You may also want to test 200 nM versus 500 nM and see which is best for your sample. From Supplementary Figure 1, panels A and C, 200 nM seems to be enough.

-We have also replaced polyA for Capture sequence 2 (Feature), libraries look very good as well. This was done to avoid competition between CMO/LMOs and endogenous RNAs for polyA oligos.

Sequence below:

5'-CCTTGGCACCCGAGAATTCCAXXXXXXXXXXGCTCACCTATTAGCGGCTAAGG-3'

Materials

Materials are those in Multi-Seq paper, McGinnis *et al.*, Nat Med, 2019. You will also need either DNA binding dies like DAPI, 7-AAD or DRAQ-7 if you are using FANS. For sucrose gradient use this NUC201-1KT (Sigma).

Safety warnings

 Every below step was done 4°C.



Before start

-I strongly recommend to read Multi-Seq paper, McGinnis *et al.*, Nat Med, 2019.

-Multi-seq is a great tech and not difficult, and yet it needs a bit optimisation to get it to work consistently. Chris and David are amazingly helpful blokes and always keen to help if you have any questions. Don't hesitate in getting in touch...I did, a lot!

-Be mindful these are just my notes and by no means are ground truth of how Multi-seq needs to be performed. I WELCOME to comments, suggestions and amendments.

- 1 First thing to do is to get **very clean nuclei preps**. You may start with less clean nuclei if you like but it is likely that CMO will bind to the debris, and that is not good (just a thought). So, if you work with cell lines, iPSCs, PBMCs or some of the 'liquid' cancers then nuclei will be pretty clean, but depending on the solid tissues, as you probably have experienced already, you may end up getting dirty nuclei preps.
- 2 To get clean nuclei you can use either a sucrose gradient or FANS based on DNA content staining DAPI/7-AAD/DRAQ-7^(*). I prefer FANS (but this is just me and because I have the sorter in the lab), and for most snRNA-Seq experiments I use DAPI or 7-AAD or DRAQ-7 (I used DAPI more often). After FANS nuclei usually looks super clean and start from there. Sucrose gradient works also quite well, specially to remove myelin from brain tissue, but haven't used it regularly enough. In any case, starting with clean nuclei (when possible) is good practice to try to have the best and reproducible labelling as possible.
(*). These are also dyes are known to not distort chromatin (or very little, as opposed to PI which it does) so we keep this door open for when we do snRNA and snATAC from same sample.
Note: FANS based on DNA content is not just an expensive way to clean up nuclei, it can help you identify potential sub-population with different ploidy, in particular for cancer samples.
- 3 It is important that you visually inspect the nuclei preps under the optical microscope (brightfield) before and after sorting. Check for signs of nuclei damage or partially lysed cells. Nuclei sizes and shape vary a lot, so I don't take this too much into account. I do look at the edges of the nuclei - usually intact nuclei have clear edges. Damaged nuclei are easy to identify with 40x-60x magnification. Sometimes disorganised chromatin popping our out of the nucleus is visible, which is bad sign. Nuclei blebbing is also bad sign.
- 4 Mix Anchor and Barcoded oligos in 1:1 molar ratio in PBS (no $\text{Ca}^{2+}/\text{Mg}^{2+}$) at [IM] 0.8 micromolar (μM) ^(*), that is 10 times solution (of the duplex) and incubate at  Room temperature for  00:05:00 . We have been preparing this mix fresh. (*)
The working concentration is optional.

[Mix Anchor and Barcoded oligos in 1:1 molar ratio in PBS (no $\text{Ca}^{2+}/\text{Mg}^{2+}$) at

[IM] 2 micromolar (μM) ^(), that is 10 times solution (of the duplex) and incubate at*

 Room temperature for  00:05:00 . We have been preparing this mix fresh. ()*

The working concentration is optional.]

What (I think) it's key is to maintain the ratio  40 pmol [or  100 pmol] of CMOs (Anchor-Barcoded Oligo and Co-Anchor) per 500,000 cells ([IM] 200 nanomolar (nM)

[or 500 nanomolar (nM)] final in 200 μ L. So for 5x less nuclei, namely 100,000 nuclei, we use 8 pmol [or 20 pmol] CMO, which is 5x less.

- 5 The number of nuclei you use of your choice, for example ~100,000 nuclei^(*) in Washing Buffer [PBS (no Ca²⁺/Mg²⁺) + 1 % BSA (MACS BSA Stock Solution from Miltenyi, Cat. 130-091-376)]. Note: remember that for cells you use 0.04 % BSA or no BSA in your Washing Buffer.
If FANS is not an option, then directly resuspend nuclei in 180 μ L Washing Buffer and continue with Step 7 below.
We usually chose FANS for our nuclei preps clean up step prior to snRNA-Seq (see *Frankenstein* protocol for snRNAseq using fresh and frozen samples on 10x Platform).
(*)Our sorter offset is ~40%, so we sort ~140,000 events to get ~100,000 nuclei. The sorting is done in 50 μ L of Washing Buffer in 96-well plate (culture plate, round bottom, but LoBind 1.5 mL eppis work too). I like culture plates because I always inspect nuclei under the microscope after sorting, so these are ideal.
- 6 After sorting, transfer the nuclei to a 1.5 mL LoBind eppis tube and complete with chilled Wash Buffer to ~ 180 μ L 1.5 mL (Note: when limited nuclei number, I usually use the additional volume to 'wash' the well where the nuclei come from and then transfer this wash the respective nuclei in the eppi).
- 7 Add 8 pmol Anchor-Barcode oligo (8 pmol = 10 μ L of 0.8 micromolar (μ M)) and gently but thoroughly mix by pipetting using wide-bore tip. So far we've maintained the 40 pmol CMO/500,000 nuclei (8 pmol CMO/100,000, that's 5x less). For 500,000 nuclei add 20 μ L of 2 micromolar (μ M) CMO for a final concentration of 200 nanomolar (nM) in 200 μ L final volume).
Here I am accounting for 5x less cells, so I add 10 μ L of 0.8 micromolar (μ M) CMO in a final volume of ~ 190 μ L (Note: at this stage the Anchor-Barcoded oligo concentration is 42.105 nanomolar (nM) , but it'll become 40 nanomolar (nM) when I add the Co-Anchor below (bringing the final volume from 190 μ L to 200 μ L). In any case, for 100,000 nuclei Anchor-Barcoded Oligo and Co-Anchor will at the same molarity in the end, that is 40 nanomolar (nM) or 8 pmol each in



200 μL final volume, which is 5x less of what's used for 500,000 nuclei, namely

200 nanomolar (nM)).

[Add 20 pmol Anchor-Barcode oligo (20 pmol = 10 μL of 2 micromolar (μM)) and gently but thoroughly mix by pipetting using wide-bore tip.

So far we've maintained the 100 pmol CMO/500,000 nuclei (20 pmol

CMO/100,000, that's 5x less). For 500,000 nuclei add 20 μL of

2 micromolar (μM) CMO for a final concentration of 500 nanomolar (nM) in

200 μL final volume).

Here I am accounting for 5x less cells, so I add 10 μL of 2 micromolar (μM)

CMO in a final volume of ~ 190 μL (Note: at this stage the Anchor-Barcoded oligo

concentration is 105.26 nanomolar (nM) , but it'll become 100 nanomolar (nM)

when I add the Co-Anchor below (bringing the final volume from 190 μL to

200 μL). In any case, for 100,000 nuclei Anchor-Barcoded Oligo and Co-Anchor will

at the same molarity in the end, that is 100 nanomolar (nM) or 20 pmol each in

200 μL final volume, which is 5x less of what's used for 500,000 nuclei, namely

500 nanomolar (nM)).]

8 Incubate for 00:05:00 at 4 $^{\circ}\text{C}$.

9 Add 8 pmol [or 20 pmol] Co-Anchor and gently but thoroughly mix.

10 Incubate for 00:05:00 at 4 $^{\circ}\text{C}$.

11 Add 1 mL of chilled Wash Buffer and mix by inverting 3x.

12 500 x g, 4 $^{\circ}\text{C}$, 00:05:00

13 Remove supernatant without disturbing the pellet. For this, I always leave behind ~

20 μL - 30 μL of supernatant. This is key to avoid losing the nuclei.

- 14 Repeat wash for 2 more times using  1.2 mL of chilled Wash Buffer. **Do not resuspend** nuclei after every wash – this help reducing nuclei loss.
- 15 In the final wash resuspend in  80 μ L of chilled Wash Buffer supplemented with  0.2 U/ μ l -  0.5 U/ μ l of RNase Inhibitor (Protector RNase Inhibitor, Cat. RNAINH-RO, Merck/Roche). One may use RNase inhibitor all the way, but it becomes expensive. I use RNase Inhibitor during cell lysis, and in this final resuspension and RNA/cDNA quality looks good.
- 16 If possible, count using Countess II Automated Cell Counter using Trypan Blue (tend to precipitate so it's a pain. I spin top speed every time I use) or Ethidium Homodimer. Alternatively, proceed with pooling.
- 17 Combine nuclei from multiple samples to multiplex and filter with a 40- μ m Flowmi® Cell Strainer.
- 18 Count using Countess II Automated Cell Counter.
- 19 Adjust nuclei concentration if needed.
- 20 Proceed with 10 $\times$ 3' v3 protocol.
- 21 After GEM formation and RT proceed as per 3' v3.x workflow until cDNA amplification. Elute cDNA in  35 μ L .
- 22 Transfer  34 μ L of cDNA onto cDNA Amplification mix supplemented with  1 μ L of  2.5 micromolar (μ M) of Multi-Seq Additive primer.
- 23 Perform cDNA amplification and 0.6X SPRI clean-up according to 3' v3.x workflow.
- 24 Endogenous transcript cDNA is bound to the SPRI beads. Continue protocol for endogenous transcripts without change.

- 25 Add 1.4X SPRI to supernatant to obtain a final SPRI volume of 2X SPRI. For example, the cDNA PCR reaction volume is  100 μL , therefore add another  140 μL of SPRI to the supernatant which already contains 60 μL of SPRI. The final volume will be  300 μL containing 2X SPRI. Transfer entire volume into a low-bind 1.5 mL tube and incubate for  00:05:00 at room temperature.
- 26 Place tube on magnet and wait ~  00:02:00 until solution is clear and then carefully remove and discard the supernatant.
- 27 Add 400 μL 80% Ethanol to the tube without disturbing the pellet and stand for  00:00:30 (only one 80% Ethanol wash). Carefully remove and discard the ethanol wash.
- 28 Centrifuge tube briefly and return it to magnet. Remove and discard any remaining ethanol.
- 29 Resuspend in beads in  50 μL water.
- 30 Perform another round of 2X SPRI purification by adding  100 μL SPRI reagent directly onto resuspended beads. Mix by pipetting and incubate  00:10:00 at room temperature.
- 31  00:02:00 Place tube on magnet and wait ~  00:02:00 until solution is clear. Carefully remove and discard the supernatant.
- 32 Add  200 μL 80% Ethanol to the tube without disturbing the pellet and stand for  00:00:30 (1st Ethanol wash). Carefully remove and discard the ethanol wash.
- 33 Add  200 μL 80% Ethanol to the tube without disturbing the pellet and stand for  00:00:30 (2nd Ethanol wash). Carefully remove and discard the ethanol wash.
- 34 Centrifuge tube briefly and return it to magnet. Remove and discard any remaining ethanol and allow the beads to air dry for  00:02:00 . **Do not over dry beads.**



35 Resuspend beads in  40 μL water or EB. Pipette mix vigorously and incubate at room temperature for  00:02:00 .

36 Place tube on magnet and transfer clear supernatant into one PCR tube.

Note

Quantification of barcode DNA concentration using Qubit (typical range is 0.5 - 5 ng/ μL) is recommended but it's optional.

-If quantifying the use  3.5 ng of material for indexing PCR.

-Alternatively use  5 μL of material for indexing.

-The use of the whole material is also possible but some primer optimization may be needed as you PCR bubbles due to over-cycling is likely (i.e. primers get exhausted). When this happen, run a P5/P7 reaction for one-2 cycles.

37 Set up the PCR reaction as follows:

 50 μL Kapa HiFi HotStart ReadyMix (2X)

 5 μL 10 μM SI Primer (e.g. 3' v3.1, PN 2000095)

 5 μL 10 μM RPIx primer (choose unique RPI for each sample from 10X lane)

 3.5 ng barcode DNA or  5 μL or entire volume (~  40 μL)

Nuclease-free water to  50 μL final volume.

Cycling conditions:

1.  95 $^{\circ}\text{C}$ x  00:05:00

2.  98 $^{\circ}\text{C}$ x  00:00:20

3.  60 $^{\circ}\text{C}$ x  00:00:30

4.  72 $^{\circ}\text{C}$ x  00:00:30

Go to 2 for ~6-10 cycles

5.  72 $^{\circ}\text{C}$ x  00:05:00

6. Hold at  4 $^{\circ}\text{C}$.



- 38 Add  160 μL (1.6X) SPRI to each PCR product, pipette mix thoroughly. Incubate at room temperature for  00:05:00 .
- 39 Place tube on magnet (High), wait for solution to clear. Remove and discard supernatant.
- 40 Wash beads twice on magnet (High) with  200 μL of 80% ethanol and allow to stand for  00:00:30 between washes.
- 41 After second wash, briefly centrifuge beads and place on magnetic rack (Low). Remove remaining ethanol with P20 micropipette. Air-dry beads on magnet for  00:02:00 .
Do not over dry beads.
- 42 Remove from magnet, resuspend beads in  25 μL buffer EB and pipette mix thoroughly to resuspend. Incubate at room temperature for  00:02:00 .
- 43 Return to magnet (Low), wait for solution to clear, and transfer supernatant to PCR strip tube.
- 44 Quantify barcode library concentration (1:5 - 1:10) using Bioanalyzer High Sensitivity DNA analysis. See representative electropherogram below.




Figure 1. This is part of Supplementary Figure 9 of paper:

(A) Bioanalyzer traces following cDNA amplification and MULTI-seq barcode enrichment using 3.2X SPRI with 1.8X 100% isopropanol exhibits two distinct peaks. Bioanalyzer traces are representative of all datasets presented in this study (n = 4). The first peak (p1) is an average of 65-70bp in length and likely corresponds to barcodes amplified via the MULTI-seq additive primer. The second peak (p2) is an average of 100bp in length and likely corresponds to barcodes that successfully underwent MMLV-RTase template switching and were subsequently amplified by the standard 10X Genomics Single Cell V2 primer.

(B) Bioanalyzer analysis following library preparation PCR exhibits one distinct peak (p3) with an average length of 173bp, matching expectations. Bioanalyzer traces are representative of all datasets presented in this study (n = 4).




Figure 2. Bioanalyzer traces. These libraries were prepared using the entire volume from step 36. Left panel inset shows PCR bubble.

- 45 Barcodes can be sequenced independently or as fraction of endogenous cDNA library. Aim for ~3000-5000 read pairs per cell.
- 46 For analysis go to <https://github.com/chris-mcginnis-ucsf/MULTI-seq>.