



Oct 01, 2025

Multiome sequencing murine ovaries

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg317z1l25/v1>

David Arboleda-Prado¹, Benjamin Cosgrove¹

¹Cornell University



David Arboleda-Prado

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg317z1l25/v1>

Protocol Citation: David Arboleda-Prado, Benjamin Cosgrove 2025. Multiome sequencing murine ovaries. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kqdg317z1l25/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 30, 2025



Last Modified: October 01, 2025

Protocol Integer ID: 228592

Keywords: mouse ovary, sequencing murine, protocol for single nuclei, single nuclei

Abstract

Protocol for single nuclei multiome sequencing - mouse ovaries.

Ovary collection and nuclei isolation

- 1 Female mice were staged for estrous cycles, and both ovaries were collected at diestrus.
- 2 Ovaries were immediately flash frozen in liquid nitrogen and stored at -80°C until use.
- 3 For single-nuclei isolation, frozen ovaries were processed using the Chromium Nuclei Isolation Kit (10x Genomics, Pleasanton, CA) following the manufacturer's protocol for frozen tissue.
- 4 Nuclei were washed, counted using an automatic hemocytometer, and resuspended in Diluted Nuclei Buffer (10x Genomics) to the concentration recommended for Multiome library preparation.

Single-nuclei Multiome library preparation and sequencing

- 5 Nuclei were loaded on the Chromium Next GEM Chip J using the Chromium Single Cell Multiome ATAC + Gene Expression kit (10x Genomics, Pleasanton, CA), following the manufacturer's instructions.
- 6 The targeted recovery was 10000 nuclei per sample.
- 7 Libraries were sequenced on a NovaSeq X platform (Illumina, San Diego, CA). Gene expression (GEX) and ATAC libraries were run across NovaSeqX-25B lanes.