

Sep 19, 2022

Multiplex fluorescent immunostaining of thin, fixed mouse brain tissue sections to characterize human iPSC-derived cell xenografts

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l275mxg1y/v1>

Benjamin Trist¹, Louise Cottle¹

¹The University of Sydney



Benjamin Trist

The University of Sydney

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l275mxg1y/v1>

Protocol Citation: Benjamin Trist, Louise Cottle 2022. Multiplex fluorescent immunostaining of thin, fixed mouse brain tissue sections to characterize human iPSC-derived cell xenografts. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.4r3l275mxg1y/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 08, 2022

Last Modified: June 01, 2024

Protocol Integer ID: 69719

Keywords: Human iPSC, Immunohistochemistry, Fluorescent, Human-to-mouse xenograft, ASAPCRN, fixed mouse brain tissue section series, mouse brain tissue section series, multiplex fluorescent immunohistochemistry protocol, fixed mouse brain tissue section, mouse brain tissue section, maturation of human ipsc, brain of athymic mice, brain tissue, human ipsc, derived cell, derived cell xenograft, immunohistochemistry, cell xenografts this protocol, living brain, ipsc, cell, mice, athymic mice

Funders Acknowledgements:

Michael J Fox Foundation

Grant ID: ASAP-000497

Abstract

This protocol describes our multiplex fluorescent immunohistochemistry protocol used to identify human iPSC-derived cells within thin, fixed mouse brain tissue section series'. We apply this workflow for post-mortem assessment of the survival, growth and maturation of human iPSC-derived cells which have been transplanted into the living brain of athymic mice.

Attachments



[519-1074.docx](#)

102KB

Materials

Equipment:

- Horizontal rocker
- Vortex
- Microcentrifuge
- Glass petri dish
- Oven

Consumables:

- 20mL scintillation vials
- Paint brushes
- Gelatin-Chrom Alum-coating microscope slides
 1. See related protocol - *Coating superfrost microscope slides with gelatin-chromium potassium sulfate*
- Microscope slide coverslips (no. 1.5 thickness, 22×50mm)
- glass pipettes
- rubber teats
- transfer pipettes

Key reagents:

- Optimal Cutting Temperature (OCT) compound
- Bovine Serum Albumin (BSA)
- Casein
- Sodium citrate
- Tween-20 and Triton X-100
- Ethanol
- ProLong Diamond Antifade mountant

Primary antibodies:

	A	B	C	D
	Target	Species	Dilution	Company, Catalog #
	alpha-Synuclein	Mouse	1:1000	Santa Cruz, sc-12767
	CTIP	Rat	1:1000	Abcam, ab18465
	FOXG1	Rabbit	1:500	Abcam, ab18259

	A	B	C	D
	HNA	Mouse	1:1000	Novus, #NOVNB-313912
	NCAM/CD56 (ERIC-1)	Mouse	1:100	Santa Cruz, sc-106
	NeuN	Chicken	1:1000	Merck-Millipore, ABN91
	Pax6	Rabbit	1:1000	ThermoFisher Scientific, #42-660
	SOX2	Mouse	1:500	R&D Systems, AF2018
	Tbr1	Rabbit	1:1000	Abcam, ab31940
	Tyrosine Hydroxylase	Rabbit	1:1500	Pel-freeze, P40101-0

✕ Anti- α -synuclein Antibody (211) **Santa Cruz Biotechnology Catalog #sc-12767**

✕ Anti-Ctip2 antibody [25B6] (ab18465) **Abcam Catalog #ab18465**

✕ Anti-FOXG1 antibody (ab18259) **Abcam Catalog #ab18259**

✕ Anti-NCAM/CD56 Antibody (ERIC 1) **Santa Cruz Biotechnology Catalog #sc-106**

✕ Anti-NeuN Antibody **Merck MilliporeSigma (Sigma-Aldrich) Catalog #ABN91**

✕ PAX6 Polyclonal Antibody **Thermo Fisher Scientific Catalog #42-6600**

✕ Human/Mouse/Rat SOX2 Antibody **R&D Systems Catalog #AF2018**

✕ Anti-TBR1 antibody **Abcam Catalog #ab31940**

✕ Rabbit Tyrosine Hydroxylase **Pel-Freez Catalog #P40101-0**

Secondary antibodies:

	A	B	C	D
	Target	Species	Dilution	Company, Catalog #
	anti-mouse AF647	Donkey	1:500	Life Technologies, A31571
	anti-rabbit AF647	Goat	1:500	Life Technologies, A11008
	anti-chicken CF594	Goat	1:500	Merck, SAB4600094
	anti-goat CF488A	Donkey	1:500	Merck, SAB4600032

✕ Donkey anti-Mouse IgG (H L) Highly Cross-Adsorbed Secondary Antibody Alexa Fluor™ 647 **Thermo Fisher Scientific Catalog #A-31571**

Goat-anti-rabbit-Alexafluor 488 **Thermo Fisher Scientific Catalog #A11008**

Anti-Chicken IgG (H L) highly cross-adsorbed CF594 antibody produced in donkey **Merck MilliporeSigma (Sigma-Aldrich) Catalog #SAB4600094**

Anti-Goat IgG (H L) highly cross-adsorbed CF™ 488A antibody produced in donkey **Merck MilliporeSigma (Sigma-Aldrich) Catalog #SAB4600032**

Solutions:

- 1x PBS, pH 7.4

- Antigen retrieval (AR) buffer

2.94 g (10 millimolar (mM)) sodium citrate, 500 μL (0.05%) Tween-20, up to 1 L with dH_2O , pH 6

- 1x PBST

- 500 μL (0.05%) Tween-20 in 1 L 1x PBS

- Blocking solution

- 1 g (1% w/v) casein, 250 μL (0.25% v/v) Triton X-100, 1.5 g (1.5% w/v) glycine, 5 g (5% w/v) BSA up to 100 mL with 1x PBS

Material input (animal, cell, tissue, fraction details)

Thin, fixed athymic mouse brain tissue sections prepared from whole mouse brains grafted with human iPSC-derived neural progenitor cells.



Day 1 (~3-4 hrs)

- 1 Pre-heat oven and Antigen Retrieval (AR) buffer to 70 °C .
- 2 Label scintillation vials to match labels on section storage plates (mouse and section series IDs, name, date etc.).
- 3 Transfer sections into scintillation vials using a transfer pipette or fine paintbrush.

- 4 Remove anti-freeze solution and perform 3× 7 min washes in 1x PBS at Room temperature with gentle agitation.

Note

- Slow shaking on an orbital rocker recommended for washes/incubations to ensure even contact with solutions.
- Use a glass pipette and rubber teat to remove solution during wash changes.
- Anti-freeze solution must be rinsed off prior to immunostaining.

- 4.1 Remove anti-freeze solution and perform 3x 00:07:00 washes in 1x PBS at Room temperature with gentle agitation (1/3).
- 4.2 Remove anti-freeze solution and perform 3x 00:07:00 washes in 1x PBS at Room temperature with gentle agitation (2/3).
- 4.3 Remove anti-freeze solution and perform 3x 00:07:00 washes in 1x PBS at Room temperature with gentle agitation (3/3).

Antigen retrieval (AR) - Day 1

- 5 Incubate sections in pre-heated AR buffer for 00:30:00 at 70 °C .
- 6 After antigen retrieval, allow sections to cool for 00:30:00 before proceeding with staining.



7 Perform 3x 7 min washes in 1x PBST with gentle agitation.



7.1 Perform 3x  00:07:00 washes in 1x PBST with gentle agitation (1/3).

7m

7.2 Perform 3x  00:07:00 washes in 1x PBST with gentle agitation (2/3).

7m

7.3 Perform 3x  00:07:00 washes in 1x PBST with gentle agitation (3/3).

7m

Blocking step - Day 1

8 Incubate sections in blocking solution for  01:00:00 at  Room temperature with

1h



Primary antibody step - Day 1

9 Incubate sections with desired primary antibody combinations (to characterise either cortical or ventral midbrain derived xenografts) diluted in blocking buffer  Overnight at  4 °C with gentle agitation.

1h



Day 2 (~4hrs)

10 Perform 3 × 7 min washes in 1x PBST with gentle agitation.



10.1 Perform 3 x  00:07:00 washes in 1x PBST with gentle agitation (1/3).

7m

10.2 Perform 3 x  00:07:00 washes in 1x PBST with gentle agitation (2/3).

7m

10.3 Perform 3 x  00:07:00 washes in 1x PBST with gentle agitation (3/3).

7m



Secondary antibody step

1d 4h 7m

- 11 Incubate sections in secondary antibodies, that match the species of chosen antibodies, diluted in blocking buffer for 02:00:00 at Room temperature . Incubate in a dark room or cover the vials with foil to protect the fluorophores from light. 2h
- 12 Perform 3× 7 min washes in 1x PBST with gentle agitation keeping vials protected from the light.
- 12.1 Perform 3x 00:07:00 washes in 1x PBST with gentle agitation keeping vials protected from the light (1/3). 7m
- 12.2 Perform 3x 00:07:00 washes in 1x PBST with gentle agitation keeping vials protected from the light (2/3). 7m
- 12.3 Perform 3x 00:07:00 washes in 1x PBST with gentle agitation keeping vials protected from the light (3/3). 7m
- 13 **Mount tissue sections** in a dimly lit room (to protect the fluorophores) onto super-frost slides pre-coated with gelatin-chrome alum and allow to dry in the dark at Room temperature for 01:00:00 - 02:00:00 . 3h
- 14 Remove ProLong Diamond Antifade mountant from storage at 4 °C and allow the reagent to equilibrate to Room temperature before use.
- 15 **Coverslip slides with ProLong Diamond Antifade mountant** and allow to set at Room temperature for at least 24:00:00 in the dark before proceeding with microscopy. 1d
- 16 Image sections using fluorescent microscopy for subsequent xenograft characterization.
- 17 Slides can be stored at 4 °C or -20 °C .