



Feb 02, 2024

## N\_terminal protein labeling

DOI

[dx.doi.org/10.17504/protocols.io.rm7vzxpbrgx1/v1](https://dx.doi.org/10.17504/protocols.io.rm7vzxpbrgx1/v1)

Patricia Yuste Checa<sup>1</sup>, F Ulrich Hartl<sup>1</sup>

<sup>1</sup>Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Patricia Yuste-Checa

MPI Biochemistry

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzxpbrgx1/v1>

**Protocol Citation:** Patricia Yuste Checa, F Ulrich Hartl 2024. N\_terminal protein labeling. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.rm7vzxpbrgx1/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** January 24, 2024



**Last Modified:** May 31, 2024

**Protocol Integer ID:** 94602

**Keywords:** ASAPCRN, n-terminal protein, using clusterin protein, clusterin protein as example, clusterin, protein, labeling, terminus, label, protocol detail

**Funders Acknowledgements:**

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

## Abstract

This protocol details how to efficiently label a protein at the N-terminus using Clusterin protein as example.

## Attachments





[973-2530.docx](#)

19KB

# Materials

## Buffers

- Labeling buffer: 1M 0.1 Molarity (M) sodium bicarbonate buffer pH 8.3
- 1x PBS pH 7.2

-  Alexa 488 NHS ester **Thermo Fisher Scientific Catalog #A20000**
-  pHrodo Red Succinimidylester **Thermo Fisher Scientific Catalog #P36600**

## Equipment

Nap-5 column	NAME
Cytiva NAP™	BRAND
45-000-151	SKU
<a href="https://www.fishersci.com/shop/products/illustra-nap-columns-nap-5-2/45000151">https://www.fishersci.com/shop/products/illustra-nap-columns-nap-5-2/45000151</a> <sup>LINK</sup>	

## Troubleshooting





## N-terminal protein labeling

1h 30m

- 1 Exchange the protein buffer to Labeling buffer using a Nap-5 column (Thermo Fisher Scientific, 45-000-151). Equilibrate the column with 10 column volumes (CV: 1 mL). Load the protein onto the column and elute the protein with the corresponding amount of Labeling buffer following the column manufacturer's instructions. Collect each eluted drop into an Eppendorf low binding tube. Measure the protein concentration in each fraction by nanodrop and pool the protein-containing fractions.

### Note

The protein will be diluted and partially lost in the process. Therefore, it is recommended to start with relatively high protein concentration. Estimated labeled Clusterin yield: around

 300  $\mu\text{L}$  at  30 micromolar ( $\mu\text{M}$ ) labeled Clusterin if starting with  100  $\mu\text{L}$  purified Clusterin at  200 micromolar ( $\mu\text{M}$ ) .

- 2 Dissolve the dye (Alexa488 NHS ester(refer materials section); pHrodo Red Succinimidylester in DMSO. With a pipette tip gently touch the dye powder which will stick to the tip. Immerse the tip in some DMSO previously dispensed in a tube. Repeat the procedure several times until the solution reaches the desired color.

### Note

If labeling a protein for uptake assays analyzed by flow cytometry, it is recommended to use A488 because the 488 nm signal outside the cell can be easily quenched by adding Trypan blue right before measurement. pHrodo red dye is a pH sensitive dye which fluoresces brightly only in acidic environments and therefore can be used to specifically monitor phagocytosis and endocytosis.

- 3 Quantify the diluted dye concentration by nanodrop. Dilute the dye with water for measurement in order to reach an absorbance of  $\lambda < 1$  for an accurate measurement.

### Note

Physical characteristics of the dyes to be set in the nanodrop:

Alexa488: Absorbance maximum ( $\lambda_{\text{max}}$ ): 495 nm; Extinction coefficient ( $\epsilon$ ): 71,000  $\text{cm}^{-1}\text{M}^{-1}$ ; Correction factor at 280 nm ( $\text{CF}_{280}$ ): 0.11; Correction factor at 260 nm ( $\text{CF}_{260}$ ): 0.3.

pHrodo Red: Absorbance maximum ( $\lambda_{\text{max}}$ ): 560nm; Extinction coefficient ( $\epsilon$ ): 65,000  $\text{cm}^{-1}\text{M}^{-1}$ ; Correction factor at 280 nm ( $\text{CF}_{280}$ ): 0.12; Correction factor at 260 nm ( $\text{CF}_{260}$ ): 0.36.




4 Add the corresponding amount of diluted dye to the eluted protein to reach a final protein:dye ratio of 1:4-1:10.



5 Incubate  01:30:00 at  Room temperature in the dark.

1h 30m



6 Remove free dye by using a Nap-5 column, pre-equilibrated with 1x PBS  7.2 buffer or desired final buffer. Collect each eluted drop in an Eppendorf low binding tube. Measure protein concentration in each fraction by nanodrop and pool the fractions containing protein.

#### Note

Avoid pooling the very last fractions because they may contain free dye.

7 Measure labeling efficiency (dye molarity/protein molarity) by nanodrop. Quantify independently dye concentration and protein concentration.

#### Note

If labeling Clusterin, it should be noted that it is a heterodimer containing 2 N-termini and therefore the labeling efficiency may be higher than 1.