

Aug 12, 2019

NADH Oxidase Activity

DOI

dx.doi.org/10.17504/protocols.io.3jcgkiw



Eva Feldman¹

¹University of Michigan - Ann Arbor

Diabetic Complications Consortium
Tech. support email: rmcindoe@augusta.edu



Lili Liang

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.3jcgkiw

External link: <https://www.diacomp.org/shared/document.aspx?id=36&docType=Protocol>

Protocol Citation: Eva Feldman 2019. NADH Oxidase Activity. protocols.io <https://dx.doi.org/10.17504/protocols.io.3jcgkiw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 31, 2019

Last Modified: August 12, 2019

Protocol Integer ID: 23876

Keywords: NADH Oxidase activity, Cardiovascular, Nephropathy, Neuropathy, Pediatric Endocrinology, Retinopathy, Uropthy, Wound-Healing

Abstract

Summary:

Describes assay to quantitate NADH Oxidase activity from tissues.

Diabetic Complications:



Cardiovascular



Nephropathy



Neuropathy



Pediatric Endocrinology



Retinopathy



Uropathy



Wound-Healing



Materials

MATERIALS

⊗ 50mM Tris **Gibco - Thermo Fisher Scientific Catalog #15504-012**

⊗ 50mM MES buffer **Merck MilliporeSigma (Sigma-Aldrich)**

⊗ 150μM NADH **Merck MilliporeSigma (Sigma-Aldrich) Catalog #N 6879**

Reagent Preparation:

Tris-MES buffer (pH 7.0):

NADH

Tris-MES buffer (pH 7.0): Prepare 50mM Tris buffer solution and pH to 7.0 with 50mM MES. 50mM Tris - 302.85mg Tris in 50ml deionized H₂O. 50mM MES - 319.89mg MES in 50ml deionized H₂O.

NADH: Prepare 600μM solution. 2.55mg in 6mL de-ion H₂O, enough for whole plate.

Note:

Fisher Scientific, RRID:SCR_008452

Sigma-Aldrich, RRID:SCR_008988

Sample Preparation:

- 1 Turn on Multiskan, set temp to 37°C and set up plate layout.
- 2 Sonicate tissue on ice in 20mM PB pH 7.4 with PMSF inhibitor or thaw prepared samples on ice.
- 3 Remove 25µL for protein analysis.
- 4 Prepare NADH, enough for whole plate.
- 5 Dilute samples 1:5 with de-ionized H₂O.
- 6 Using a clear plate: Add 50µL sample to wells and 50µL diluted sample to wells in duplicate.
- 7 Add 50µL buffer to 3 wells for blanks for positive control.
- 8 Add 100µL Tris-Mes to each sample and blanks.
- 9 For negative control add 200µL Tris-Mes to 3 wells.
- 10 Place plate in Multiskan and add 50µL 600µM NADH to the sample and positive blanks.
Do not add NADH to the 3 negative control wells.
- 11 Press start and read at 340nm for 10 minutes @ 1 minute intervals.



- 12 Save raw data as an Excel file into the NADHx data folder. Use the naming convention NAXXXX.xls, where XXXX is the date in mmdd format.