

Mar 03, 2025

Nanoparticle Tracking Analysis (NTA)

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gp9761vzp/v1>

Emily Tabaie¹

¹University of California, Riverside



Emily Tabaie

University of California, Irvine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gp9761vzp/v1>

Protocol Citation: Emily Tabaie 2025. Nanoparticle Tracking Analysis (NTA). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.dm6gp9761vzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 03, 2025



Last Modified: March 03, 2025

Protocol Integer ID: 123713

Keywords: nanoparticle tracking analysis, nanoparticle, concentration of particle, nta, isolation of extracellular vesicle, extracellular vesicle, particle, tracking analysis, concentration

Abstract

This protocol is designed to identify the average size and concentration of particles using a nanoparticle tracking analysis (NTA). These steps are to be done immediately after the isolation of extracellular vesicles.

Nanoparticle Tracking Analysis (NTA)

- 1 Turn on the machine and tighten in the 4 beige screws into the glass chamber
- 2 Place the glass chamber into the machine and move up the red lever to lock into place
- 3 Grab a new 1 mL syringe and wash the machine with 1 mL of distilled water 2-3 times until no particles are visible
- 4 Attach the syringe to the injection port and fix it in the injection port
- 5 Go to capture (left top corner) → click start camera
- 6 Go to interface → hardware (on the right) → syringe pump (on the top) → click infuse
- 7 Scroll up and down to see the brightness of the particles passing on the screen
- 8 Once the wash is done grab that same 1 mL syringe and take up all 1 mL of your sample
- 9 Insert the syringe and state infusion
- 10 Wait until the volume on the syringe reaches 300-400 ul and start the analysis
- 11 Go SOP (left) → recent measurement (top) → double click on generic triplicate (will show the script panel on the left) → click on file name → click on D drive → Ziting → Emily → create a folder with the date
- 12 Hit run → creator name (Emily) → diluent (100) → sample description (whatever the sample is) → okay → settings okay



- 13 Hit settings okay → export → click on the folder → drive → folder → results → pdf → export to USB drive
- 14 Take the syringe and inject into the waste and throw the syringe into the red biohazard bag
- 15 Get a new syringe and insert 1 mL of distilled water into the syringe pump to wash the pump and repeat once
- 16 Go to hardware → syringe pump → infuse
- 17 When all the samples are complete wash the chamber with a total of 3 mL distilled water
- 18 Clean the glass chamber with ethanol, turn off the machine, and close the program