

Nov 04, 2025

Nanopore Sequencing of Environmental Microbial Metagenomes and Metatranscriptomes Using the Native Barcoding Kit 24 V14

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gpzz9plzp/v1>

Thuc Trong Nguyen¹, Runar Stokke¹

¹Department of Biological Sciences and Centre for Deep Sea Research, University of Bergen



Thuc Trong Nguyen

Department of Biological Sciences, University of Bergen, Cen...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpzz9plzp/v1>

Protocol Citation: Thuc Trong Nguyen, Runar Stokke 2025. Nanopore Sequencing of Environmental Microbial Metagenomes and Metatranscriptomes Using the Native Barcoding Kit 24 V14. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.dm6gpzz9plzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 31, 2024

Last Modified: November 04, 2025

Protocol Integer ID: 101013

Keywords: Nanopore sequencing, Direct cDNA sequencing, Native barcoding, Metagenomic, Metatranscriptomic, Biofilm, nanopore sequencing of environmental microbial metagenome, nanopore sequencing, environmental microbial metagenome, streamlined workflow for nanopore, environmental biofilm metagenome, oxford nanopore technology, nanopore, studying microbial community, microbial communities in environment, read sequencing, rna extraction, barcode ligation, sequencing output, native barcoding kit 24 v14, metatranscriptome

Funders Acknowledgements:

Research council of Norway

Grant ID: 315427

Disclaimer

FOR RESEARCH PURPOSES ONLY; USE AT YOUR OWN RISK

This protocol is provided for informational purposes only and does not constitute legal, medical, clinical, or safety advice. Content added to **protocols.io** is not peer-reviewed and may not have undergone formal approval. The information in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action taken or omitted based on this protocol is solely at the user's own risk.

Abstract

This protocol presents a streamlined workflow for **Nanopore sequencing of environmental biofilm metagenomes and metatranscriptomes** using the **Native Barcoding Kit 24 V14**. It details the DNA/RNA extraction, rRNA depletion, polyadenylation, cDNA synthesis, barcoding, and library preparation steps, optimized for long-read sequencing. The workflow incorporates magnetic bead-based purification and barcode ligation to enable multiplexed sequencing on an Oxford Nanopore Technologies (ONT) platform. Designed to maximize DNA/cDNA yield, minimize purification losses, and ensure high-quality sequencing output, this method is well-suited for studying microbial communities in environments, including deep-sea hydrothermal vents and other extreme habitats.

Attachments



[r2002_zymbiomics_d...](#)

351KB



[ligation-sequencing-...](#)

1.3MB



[ligation-sequencing-...](#)

1.6MB



[highprep_rna_elite_p...](#)

54KB



[ribopool-protocol.pd...](#)

2.2MB



[Ampure XP beads.pdf](#)

565KB

Guidelines

Using the **Native Barcoding Kit 24 V14**

Key Considerations

- **Sample Quality:** Ensure proper collection and storage of environmental samples to prevent degradation.
- **RNA Handling:** Use RNase-free reagents and clean the workspace with RNase AWAY to avoid contamination.
- **Minimizing Sample Loss:** Be mindful of purification steps, as RNA losses can be significant.
- **Bead-Based Purification:** Magnetic bead clean-ups must be performed carefully to maximize yield.
- **Barcode Selection:** Choose unique barcodes (NB01–24) for multiplexed sequencing.
- **Library Preparation:** Follow ONT guidelines for adapter ligation and sequencing setup.

Materials

- ZymoBIOMICS DNA/RNA Miniprep Kit R2002 (Zymo Research)
- DNase I E1010 (Zymo Research)
- AMPure XP Beads (Beckman Coulter)
- HighPrep RNA Elite RC-90005 (MagBio)
- Pan-Prokaryote riboPOOL Kit (siTOOLs Biotech)
- Qubit RNA HS Assay Kit (Thermo Fisher)
- Qubit RNA IQ Assay Kit (Thermo Fisher)
- dNTP Mix (10 mM each) R0193 (Thermo Fisher)
- RNase Cocktail Enzyme Mix AM2286 (Thermo Fisher)
- RNaseOUT Recombinant Ribonuclease Inhibitor 10777019 (Thermo Fisher)
- ATP Solution (10 mM) PV3227 (Thermo Fisher)
- Maxima Reverse Transcriptase (200 U/μL) EP0741 (Thermo Fisher)
- UltraPure BSA (50 mg/mL) (Thermo Fisher)
- VN Primer (2 μM) (/5phos/ACTTGCCTGTCGCTCTATCTTCTTTTTTTTTTTTTTTTTTTTTTVN)
- Strand-switching Primer (10 μM) (TTTCTGTTGGTGCTGATATTGCTmGmGmG)
- PR2 Primer (10 μM) (/5Phos/TTTCTGTTGGTGCTGATATTGC)
- E. coli Poly(A) Polymerase M0276S (New England Biolabs)
- Blunt/TA Ligase Master Mix M0367S (New England Biolabs)
- NEBNext Quick Ligation Reaction Buffer B6058S (New England Biolabs)
- NEBNext Companion Module for Oxford Nanopore Technologies Ligation Sequencing E7180S (New England Biolabs)
- Native Barcoding Kit 24 V14 (SQK-NBD114.24) (Oxford Nanopore Technologies)
- ONT Flow Cell R10.4.1 (Oxford Nanopore Technologies)



Safety warnings

- ❗ The loss of RNA from extraction to cDNA synthesis is very high (~98.11%). Start with **≥1000 ng of total RNA per sample** to retain enough cDNA for sequencing.

Safety

- Use PPE when handling chemicals and biological samples.
- Follow the latest version of protocols on **The Nanopore Community** for the best library prep and sequencing results.

Before start

Clean the workbench and all the equipment with **RNase AWAY®** before starting any work.



DNA/RNA extraction using ZymoBIOMICS DNA/RNA Miniprep Kit

33m 30s

- 1 In a **ZR BashingBead Lysis Tube (0.1 & 0.5 mm)**, add up to 250 mg of biofilm sample and 750 µL DNA/RNA Shield . Ensuring a total volume of approximately 1000 µL per tube. 1000 µL 5m
- 2 Close the tube caps tightly. Homogenize the sample using a Homogenizer Bead Mill or FastPrep-24 at a speed of 6 m/s for two rounds, each lasting 1 minute. Allow the sample tubes to cool On ice for 5 minutes between rounds. 10m
- 3 Centrifuge at 16000 x g, 00:00:30 30s
- 4 Carefully transfer the entire supernatant into a 1.5 mL Eppendorf DNA LoBind tube, ensuring that the pellet remains undisturbed. 2m
- 5 Add an equal volume of **DNA/RNA Lysis Buffer** to the supernatant at a 1:1 ratio (e.g., 400 µL of Lysis Buffer for 400 µL of sample). Mix thoroughly by vortexing, then briefly spin down. 2m
- 6 Transfer the sample into a **Spin-Away Filter** (yellow) placed in a **Collection Tube**. Centrifuge at 16000 x g, 00:00:30 . Carefully transfer the **flow-through** to a clean 1.5 mL Eppendorf DNA LoBind tube for **RNA** collection, and retain the filter for DNA purification. 30s
- 7 Repeat the previous step until the entire mixture from **Step 5** has been filtered. 10m

STEP CASE

Metatranscriptomic sequencing 153 steps

RNA extraction

55m 30s

- 8 Add an equal volume of absolute ethanol to the flow-through and mix thoroughly by pipetting. 3m



- 9 Transfer the mix into a **Zymo-Spin IIICG Column** (green) placed in a **Collection Tube**. Centrifuge at 16000 rpm, 00:00:30 . Discard the flow-through. If necessary, reload the column and repeat the centrifugation.

30s

9.1 DNase I Treatment in-column

30s

Add 400 μ L DNA/RNA Wash Buffer to the column. Centrifuge at 16000 x g, 00:00:30 . Discard the flow-through

- 9.2 Prepare the following mix in a clean Eppendorf tube, adjusting the quantities based on the number of columns you use.

2m

Reaction mix	Volume
DNase I reconstituted (1U/ μ L)	5 μ L
DNA digestion buffer	75 μ L
Total volume per reaction	80 μL

Mix by gentle inversion, then briefly spin down. **Do not vortex.**

- 9.3 Add 80 μ L DNase I reaction Mix (from step 9.2) directly into the column matrix. Incubate at Room temperature for 00:15:00

15m



- 10 Add 400 μ L DNA/RNA Prep Buffer to the column. Centrifuge at 16000 rpm, 00:00:30 . Discard the flow-through

30s

- 11 Add 700 μ L DNA/RNA Wash Buffer to the column. Centrifuge at 16000 rpm, 00:00:30 . Discard the flow-through

30s

- 12 Add 400 μ L DNA/RNA Wash Buffer to the column. Centrifuge at 16000 x g, 00:02:00 . Carefully transfer the column into a new nuclease-free Eppendorf tube.

2m

- 13 Add 100 μ L ZymoBIOMIC DNase/Rnase-Free Water directly to the column matrix. Incubate at Room temperature for 00:05:00 . Centrifuge at

5m 30s





 16000 x g, 00:00:30 to elute RNA from the column.

Note

Alternatively, for highly concentrated RNA, use ≥ 50 mL elution

- 14 Place a **Zymo-Spin III-HRC Filter** in a Collection Tube. Add  600 μ L ZymoBIOMICS HRC Prep Solution . Centrifuge at  8000 x g, 00:03:00 .

3m

- 15 Transfer the eluted RNA from **Step 13** into a prepared **Zymo-Spin III-HRC Filter (Step 14)** placed in a nuclease-free Eppendorf DNA LoBind tube. Centrifuge at  16000 x g, 00:03:00 .

3m



Note

The filtered RNA can be used immediately or stored frozen at  -80°C .

Quantify  1 μ L of the eluted sample using Qubit fluorometer (Qubit RNA HS).

Total RNA Purification using HighPrep RNA Elite Beads

49m 30s

- 16 Thoroughly shake the **HighPrep RNA Elite** reagent before use.

1m

Note

Allow the **HighPrep RNA Elite** reagent to reach  Room temperature for at least  00:30:00 before use.

- 17 Transfer the total RNA sample into a clean 1.5 mL Eppendorf tube.

1m



- 18 Add the appropriate volume of **HighPrep RNA Elite** reagent based on the volume of your total RNA sample.

2m

	Volume of RNA (μL)	Volume of HighPrep RNA Elite 1.8X (μL)
	10	18
	50	90
	100	180

- 19 Mix thoroughly by pipetting up and down 6-8 times.

1m

- 20 Incubate on a Hula mixer for 00:05:00 at Room temperature .

5m



- 21 Place the mixture on a magnetic separation device for 00:15:00 or until the solution becomes clear.

15m

- 22 While keeping the Eppendorf tube **on the magnet**, carefully pipette off and discard the supernatant.

1m

- 23 With the tube **on the magnet**, add 200 μL 70% ethanol to the Eppendorf tube. Incubate at Room temperature for 00:00:30

30s

- 24 While keeping the Eppendorf tube **on the magnet**, carefully pipette off and discard the supernatant.

1m

Note

If the pellet is disturbed, allow the beads to settle again before removing the ethanol.

- 25 Repeat steps 23-24 for **three** 70% ethanol **washes**.

5m



- 26 Dry the beads by incubating the Eppendorf tube **on the magnetic separation device** at Room temperature for up to 00:10:00 .

10m

**Note**

Ensure that all traces of ethanol are completely removed. However, avoid overdrying the beads, as a cracked dry pellet may reduce the yield.

- 27 Remove the tube from the magnetic separation device. Add 40 μ L elution buffer (e.g., nuclease-free water) to the reaction tube and mix by pipetting up and down.

1m

Note

Prewarming the nuclease-free water to 55 $^{\circ}$ C can help increase the yield.

- 28 Incubate the tube on a Hula mixer for 00:02:00 at Room temperature .

2m

- 29 Briefly spin down the tube. Place it back on the magnetic separation device and wait for 00:03:00 or until the beads have completely cleared from the solution.

3m

- 30 Transfer the eluate (cleared supernatant) to a clean 1.5 mL Eppendorf tube for storage or subsequent applications.

1m

**Note**

Quantify 1 μ L of the eluted sample using Qubit fluorometer (Qubit RNA HS).

rRNA depletion using riboPOOL Kit

59m 30s



31 Add the following reagent to a 1.5 mL Eppendorf tube

2m

	Reagent	Volume
	Hybridization Buffer (HB)	5 μ L
	Resuspended RP	1 μ L
	RNA sample (containing 100ng - 5 ug of total RNA)	14 μ L
	RNaseOUT™	0.5 μ L
	Total	20.5 μL

32 Vortex thoroughly and briefly spin down.

30s

33 Incubate at 68 °C for 00:10:00 to denature RNA.

10m

34 Allow the sample to cool gradually from 68 °C to 37 °C for optimal hybridisation.

15m

Note

To achieve this, turn off the heating block and allow the temperature to naturally decrease to 37°C. If using temperature-controlled ramping, cool down at a rate of 3°C per minute.

35 Resuspend the **streptavidin-coated magnetic beads (SMB)** by carefully vortexing the tube at medium speed.

30s

36 Transfer 90 μ L of bead suspension per sample into a fresh tube.

30s

Note

For batch washing of beads for multiple samples, aliquot the bead suspension into a single tube: use up to 540 μ L for 6 samples or 1080 μ L for 12 samples.



- 37 Place the tube containing **SMB beads** on a magnetic rack and wait for  00:01:00 . 
- 38 Remove and discard the supernatant carefully. 
- 39 Add  80 μL of **Depletion Buffer (DB)** per sample (e.g.,  480 μL for 6 samples and  960 μL for 12 samples). Agitate the tube well to fully resuspend the beads. 
- 40 Repeat **Steps 37-39**. 
- 41 Briefly centrifuge the tube containing  20.5 μL of hybridized riboPOOL and total RNA (from **Step 34**). 
- 42 Pipette  80 μL of the prepared beads (from **Step 40**) into the tube containing hybridised riboPOOL-RNA solution. Agitate the tube thoroughly to ensure proper resuspension. 
- 43 Incubate the tube at  37 °C for  00:15:00 , followed by  50 °C for  00:05:00 in a thermal cycler. 

- 44 Briefly spin down to collect any droplets. 
- 45 Place on the magnet for  00:02:00 . Then, carefully transfer the supernatant to a new 1.5 mL Eppendorf tube. 
- 46 Place the tube on the magnet for  00:01:00 to remove any trace amounts of beads (optional). Carefully transfer the supernatant to a new tube. 
 

Note

At this stage, the elution should primarily contain **mRNA**. The RNA concentration is expected to decrease by approximately **90%** from the initial amount due to the depletion of rRNA from the total RNA sample.

mRNA Purification using HighPrep RNA Elite

50m

- 47  [go to step #16](#) Follow the same protocol to **purify total RNA** precisely (from **Step 16** to **Step 30**). The volume of elution buffer (**Step 27**) can be adjusted to  30 μL to achieve a higher RNA concentration.

Note

Quantify  1 μL of the eluted RNA using Qubit fluorometer (Qubit RNA HS)
Qualify  1 μL of the eluted RNA using either Bioanalyzer RNA 6000 Pico Kit (Agilent) or Qubit fluorometer (Qubit RNA IQ)

Polyadenylating RNA with *E. coli* poly(A) polymerase

39m 30s

- 48 In a 1.5 mL Eppendorf tube, set up a 3' polyadenylation reaction as follows

2m



Reagent	Volume
Non-polyadenylated RNA (from step 47)	X μL (e.g., 10 μL)
10X <i>E. coli</i> Poly(A) Polymerase Reaction Buffer	2 μL
ATP (10mM)	2 μL
<i>E. coli</i> Poly(A) Polymerase	1 μL
Nuclease-free water	14.5 μL - X μL (e.g., 14.5 - 10 = 4.5 μL)
RNaseOUT™	0.5 μL
Total	20 μL



Gently agitate by pipetting up and down. **Do not vortex.**

- 49 Incubate reaction at 37 °C for 00:10:00 . 10m
- 50 Stop the reaction by adding 5 µL of 50 millimolar (mM) EDTA (to achieve a final concentration of 10 millimolar (mM) EDTA). The final reaction volume will be 25 µL . 30s
- 51 Add 45 µL of resuspended **HighPrep RNA Elite Beads** to the reaction tube. 30s
- 52 Incubate on a Hula mixer for 00:05:00 at Room temperature 5m
- 53 Spin down the sample, then pellet it using a magnetic separation device. **Keep the tube on the magnet**, carefully pipette off and discard the supernatant. 5m
- 54 Keep the tube **on the magnet** and wash the beads with 200 µL of freshly prepared 70% ethanol. Gently rotate the tube 180° twice in the rack to ensure thorough washing of the pelleted beads. Carefully pipette off and discard the supernatant. 2m
- 55 Repeat **Step 54** to perform a total of two washes. 2m
- 56 Briefly spin down the tube and place it back **on the magnet**. Carefully pipette off any residual ethanol. Allow the beads to dry for 00:00:30 , ensuring the pellet does not over-dry or crack. 30s
- 57 Remove the tube from the magnetic rack and resuspend the pellet in 12 µL of nuclease-free water. Incubate On ice for 00:05:00 . 5m
- 58 Place the tube **on the magnet** until the beads pellet and the eluate becomes clear and colourless. 5m
- 59 Carefully remove and retain 12 µL of the eluate containing the 3'-polyadenylated RNA in a clean 1.5 mL Eppendorf tube. 2m



Note

Quantify  1 μL of the eluted sample using Qubit fluorometer (Qubit RNA HS)

Reverse transcription and strand-switching

2h 4m

60 Prepare the following reaction in the 0.2 mL PCR tube:

2m

Reagent	Volume
RNA (100 ng Poly(A)+ RNA)	7.5 μL
VN Primer diluted to 2 uM	2.5 μL
10 mM dNTPs	1 μL
Total	11 μL

61 Mix gently by flicking the tube, then briefly spin down.

30s

62 Incubate at  65 $^{\circ}\text{C}$ for  00:05:00 and then snap cool on a pre-chilled freezer block for  00:01:00

6m

63 In a separate tube, mix the following:

2m

Reagent	Volume
5x RT Buffer	4 μL
RNase OUT	1 μL
Nuclease-free water	1 μL
Strand-switching Primer diluted to 10 uL	2 μL
Total	8 μL

63.1 Mix gently by flicking the tube, and spin down.

30s

64 Add  8 μL of the **strand-switching reagents** (prepared in **Step 63**) to the  11 μL of snap-cooled mRNA (from **Step 62**). Mix gently by flicking the tube, then briefly spin down. 1m

65 Incubate at  42 °C for  00:02:00 in the thermal cycler. 2m

66 Add  1 μL of **Maxima H Minus Reverse Transcriptase**, bringing the total volume to  20 μL .

67 Mix gently by flicking the tube, and spin down.

68 Incubate using the following protocol using a thermal cycler: 1h 50m

	Cycle step	Temperature	Time	No. of cycles
	Reverse transcription and strand-switching	42°C	90 mins	1
	Heat inactivation	85°C	5 mins	1
	Hold	4°C	∞	



RNA degradation and second-strand synthesis

1h 36m 30s

69 Add  1 μL of **RNase Cocktail Enzyme Mix** to the reverse transcription reaction (end product from **Step 68**). 1m

70 Incubate the reaction for  00:10:00 at  37 °C in the thermal cycler. 10m



71 Resuspend the **AMPure XP beads (AXP)** by vortexing thoroughly. 1m



- 72 Transfer the sample to a clean 1.5 mL Eppendorf DNA LoBind tube. 1m
- 73 Add  17 μL of resuspended **AMPure XP beads (AXP)** to the reaction and mix gently by flicking the tube. 1m
- 74 Incubate on a Hula mixer for  00:05:00 at  Room temperature . 5m
- 75 Spin down the sample, then pellet it using a magnetic separation device until the supernatant is clear and colourless. Keep the tube **on the magnetic** device and carefully pipette off the supernatant. 5m
- 76 Keep the tubes on the magnetic device and wash the beads with  200 μL of freshly prepared 80% ethanol, ensuring the pellet remains undisturbed. Carefully remove and discard the ethanol using a pipette. 2m
- Note

If the pellet is disturbed, allow the beads to re-pellet before removing the ethanol.
- 77 Repeat **Step 76** for a total of two washes. 2m
- 78 Spin down the tube and place it back on the magnetic device. Carefully pipette off any residual ethanol. Allow the beads to dry for approximately  00:00:30 , ensuring the pellet does not over-dry or crack. 30s
- 79 Remove the tube from the magnetic rack and resuspend the pellet in  20 μL nuclease-free water. 1m
- 80 Incubate on a Hula mixer for  00:10:00 at  Room temperature 10m
- 81 Briefly spin down the tube, then pellet the beads on the magnet until the eluate is clear and colourless, for at least  00:01:00 1m



82 Carefully remove and retain  20 μ L of eluate into a clean 1.5 mL Eppendorf DNA LoBind tube.

1m

83 Prepare the following reaction in a 0.2 mL thin-walled PCR tube:

3m

Reagent	Volume
2x LongAmp Taq Master Mix	25 μ l
PR2 Primer diluted to 10 μ M	2 μ l
Reverse-transcribed sample from above (Step 82)	20 μ l
Nuclease-free water	3 μ l
Total	50 μl

Mix gently by flicking the tube

84 Incubate using the following protocol in a thermal cycler:

20m

Cycle step	Temperature	Time	No. of cycles
Denaturation	94 °C	1 min	1
Annealing	50 °C	1 min	1
Extension	65 °C	15 mins	1
Hold	4 °C	∞	

85 Resuspend **the AMPure XP beads (AXP)** by vortexing thoroughly.

1m

86 Transfer the sample (end product of **Step 84**) to a clean 1.5 mL Eppendorf DNA LoBind tube.

1m

87 Add  40 μ L of resuspended **AMPure XP beads (AXP)** to the reaction and mix gently by flicking the tube.

1m

88 Incubate on a Hula mixer for  00:05:00 at  Room temperature .

5m





89 Spin down the sample, then pellet it using a magnetic separation device until the supernatant is clear and colourless. Keep the tube **on the magnetic device** and carefully pipette off the supernatant. 5m

90 Keep the tubes **on the magnetic device** and wash the beads with  200 μL of freshly prepared 80% ethanol, ensuring the pellet remains undisturbed. Carefully remove and discard the ethanol using a pipette. 2m

Note

If the pellet is disturbed, allow the beads to re-pellet before removing the ethanol.

91 Repeat **Step 90** for a total of two washes. 2m

92 Spin down the tube and place it back **on the magnetic device**. Carefully pipette off any residual ethanol. Allow the beads to dry for approximately  00:00:30, ensuring the pellet does not over-dry or crack. 1m

93 Remove the tube from the magnetic rack and resuspend the pellet in  21 μL nuclease-free water. 1m

94 Incubate on a Hula mixer for  00:10:00 at  Room temperature. 10m

95 Briefly spin down the tube, then pellet the beads on the magnet until the eluate is clear and colourless, for at least  00:01:00. 1m

96 Carefully remove and retain  21 μL of eluate in a clean 1.5 mL Eppendorf DNA LoBind tube. 2m

Note

Quantify  1 μL of the eluted sample using Qubit fluorometer (Qubit DNA HS)

cDNA repair and end-prep

34m 30s



97 Combine the following reagents in a 0.2 m PCR tube:

3m

Reagent	Volume
cDNA sample (Step 96)	20 µl
Nuclease-free water	30 µl
Ultra II End-prep reaction buffer	7 µl
Ultra II End-prep enzyme mix	3 µl
Total	60 µl

Note

Ensure the reaction is thoroughly mixed by gently pipetting up and down, then briefly spin down.

98 Using a thermal cycler, incubate at 20 °C for 00:05:00 , followed by 65 °C for 00:05:00 .

10m



99 Resuspend the **AMPure XP Beads (AXP)** by vortexing thoroughly.

1m

100 Transfer the sample (**Step 98**) to a clean 1.5 mL Eppendorf DNA LoBind tube.

1m

101 Add 60 µL of resuspended **AMPure XP Beads (AXP)** to the end-prep reaction (**Step 100**) and mix gently by flicking the tube.

2m

102 Incubate on a Hula mixer for 00:05:00 at Room temperature .

5m

103 Spin down the sample, then pellet it using a magnetic separation device until the supernatant is clear and colourless. Keep the tube **on the magnetic** device and carefully pipette off the supernatant.

5m

104 Keep the tubes **on the magnetic device** and wash the beads with 200 µL of freshly prepared 80% ethanol, ensuring the pellet remains undisturbed. Carefully remove and discard the ethanol using a pipette.

2m

Note

If the pellet is disturbed, allow the beads to re-pellet before removing the ethanol.

105 Repeat **Step 104** for a total of two washes

2m

106 Spin down the tube and place it back **on the magnetic device**. Carefully pipette off any residual ethanol. Allow the beads to dry for approximately  00:00:30, ensuring the pellet does not over-dry or crack.

30s

107 Remove the tube from the magnetic rack and resuspend the pellet in  20 μ L nuclease-free water. Incubate for  00:02:00 at room temperature.

2m



108 Pellet the beads on a magnetic rack until the eluate is clear and colourless, for at least  00:01:00

1m

109 Remove and retain  20 μ L of eluate into a clean 1.5 mL Eppendorf DNA LoBind tube.

II

Note

Quantify  1 μ L of the eluted sample using a Qubit fluorometer (Qubit DNA HS)

At this point, cDNA can be stored at  -20 $^{\circ}$ C for up to a month or at  -80 $^{\circ}$ C for longer storage.

Native barcode ligation

1h 6m 30s

110 For each sample to be run together on the same flow cell, select a unique barcode from the **Native Barcodes (NB01-24)**.

1m

111 In clean 0.2 mL PCR tubes, add the reagents in the following order per tube:

2m

Reagent	Volume
End-prepped DNA (from Step 109)	7.5 μ l



Reagent	Volume
Native Barcode (NB01-24)	2.5 µl
Blunt/TA Ligase Master Mix	10 µl
Total	20 µl

Note

Thoroughly mix the reaction by gently pipetting up and down, then briefly spin down.

112 Incubate for  00:20:00 at  Room temperature .

20m



113 Add the appropriate volume of EDTA to each tube, selecting the correct version based on the bottle cap colour. Mix thoroughly by pipetting, then briefly spin down.

2m

EDTA cap colour	Volume per tube
For clear cap EDTA	2 µl
For blue cap EDTA	4 µl

EDTA is added at this step to stop the reaction.

Note

The concentration of EDTA may vary depending on the version of your Nanopore reaction kit, as detailed below.

- The EDTA vials with a **clear cap** contain 0.5M EDTA (old version concentration).
- The EDTA vials with a **blue cap** contain 0.25M EDTA (new version concentration).

The new EDTA concentration is indicated in all batch numbers formatted as NBD1424.**20**.0001 and NBD1424.**20**.0002 (or higher). The old format was used in batch numbers following the format NBD1424.**10**.0001.

114 **Pool** all the barcoded samples in a 1.5 mL Eppendorf DNA LoBind tube.

2m



-	Volume per sample	For 6 samples	For 12 samples	For 24 samples
Total volume of preps using clear cap EDTA	22 µl	132 µl	264 µl	528 µl
Total volume of preps using blue cap EDTA	24 µl	144 µl	288 µl	576 µl

Note**IMPORTANT:**

Ensure that the liquid volumes at the base of all tubes or wells in the plate are consistent before and after pooling to confirm that all liquid has been transferred correctly.

115 Resuspend the **AMPure XP Beads (AXP)** by vortexing.

1m

116 Add **AMPure XP Beads (AXP)** to the pooled reaction at a 0.4X ratio and mix thoroughly by pipetting.

1m

-	Volume per sample	For 6 samples	For 12 samples	For 24 samples
Volume of AXP for preps using clear cap EDTA	9 µl	53 µl	106 µl	211 µl
Volume of AXP for preps using blue cap EDTA	10 µl	58 µl	115 µl	230 µl

117 Incubate on a Hula mixer for  00:10:00 at  Room temperature

10m





118 Spin down the sample, then pellet it using a magnetic separation device until the supernatant is clear and colourless. Keep the tube **on the magnetic device** and carefully pipette off the supernatant. 5m

119 Keep the tubes **on the magnetic device** and wash the beads with  700 μL of freshly prepared 80% ethanol, ensuring the pellet remains undisturbed. Carefully remove and discard the ethanol using a pipette. 2m

Note

If the pellet is disturbed, allow the beads to re-pellet before removing the ethanol.

120 Repeat **Step 119** for a total of two washes 2m

121 Spin down the tube and place it back **on the magnetic device**. Carefully pipette off any residual ethanol. Allow the beads to dry for approximately  00:00:30, ensuring the pellet does not over-dry or crack. 30s

122 Remove the tube from the magnetic rack and resuspend the pellet in  35 μL nuclease-free water by gently flicking. 1m

123 Incubate for  00:10:00 at  37 °C. Every  00:02:00, gently flick the sample for  00:00:10 to promote DNA elution. 10m



124 Pellet the beads **on a magnetic rack** until the eluate is clear and colourless. 5m

125 Carefully remove and retain  35 μL of eluate into a clean 1.5 mL Eppendorf DNA LoBind tube. 2m

Note

Quantify  1 μL of the eluted sample using a Qubit fluorometer (Qubit DNA HS)



Adapter ligation and clean-up

1h 8m

- 126 Thaw all reagents at Room temperature . Mix gently by pipetting, then briefly spin down. Place all reagents On ice .

3m

In a 1.5 ml Eppendorf LoBind tube, mix in the following order:

Reagent	Volume
Pooled barcoded sample	30 μ L
Native Adapter (NA)	5 μ L
NEBNext Quick Ligation Reaction Buffer (5X)	10 μ L
Quick T4 DNA Ligase	5 μ L
Total	50 μL

Note

IMPORTANT

The **Native Adapter (NA)** used in this kit and protocol is not interchangeable with other sequencing adapters.

Do **NOT** vortex the **Quick T4 DNA Ligase**.

- 127 Thoroughly mix the reaction by gently pipetting up and down, then briefly spin down.

1m

- 128 Incubate the reaction for 00:20:00 at Room temperature

20m



- 129 Resuspend the **AMPure XP Beads (AXP)** by vortexing.

1m

- 130 Add 20 μ L of resuspended **AMPure XP Beads (AXP)** to the reaction and mix thoroughly by pipetting.

1m

- 131 Incubate on a Hula mixer for 00:10:00 at Room temperature

10m



132 Spin down the sample, then pellet it on the magnetic rack. Keep the tube **on the magnet** and carefully pipette off the supernatant.

5m

133 Wash the beads by adding  125 μL **Short Fragment Buffer (SFB)**. Flick the tube to resuspend the beads, then briefly spin down. Return the tube to the magnetic rack and allow the beads to pellet. Carefully remove and discard the supernatant using a pipette.

5m



Note

IMPORTANT:

Use **Short Fragment Buffer (SFB)** rather than 80% ethanol to wash the beads, as ethanol can be detrimental to the sequencing reaction.

134 Repeat **Step 133** for a total of two washes.

5m

135 Spin down the tube and place it back **on the magnetic rack**. Carefully pipette off any residual supernatant.

2m

136 Remove the tube from the magnetic rack and resuspend the pellet in  15 μL of **Elution Buffer (EB)**.

2m

137 Spin down the tube and incubate for  00:10:00 at  37 °C . Every  00:02:00 , gently flick the sample for  00:00:10 to promote DNA elution.

10m



138 Pellet the beads **on a magnet** until the eluate is clear and colourless, for at least  00:01:00

1m

139 Carefully remove and retain  15 μL of the eluate containing the DNA library in a clean 1.5 mL Eppendorf DNA LoBind tube.

2m

Note

Quantify  1 μL of the eluted sample using a Qubit fluorometer (Qubit DNA HS)
The prepared library is now ready for loading onto the flow cell. Store the library on ice until it is ready to be loaded.

Priming and loading on the SpotON flow cell

27m 30s

140 In a 1.5 mL Eppendorf LoBind tube, mix in the following order:

2m


Reagent	Volume per flow cell
Flow cell flush (FCF)	1170 μL
Bovine serum albumin (BSA) at 50 mg/mL	5 μL
Flow cell tether (FCT)	30 μL
Total	1205 μL

Note

Mix by pipetting at room temperature

141 Place the flow cell on the sequencing device and follow the manufacturer's instructions for setup and operation.

1m



MinION Flow Cell

- 141.1 Complete a flow cell check to assess the number of pores available before loading the library. 5m
- 142 Slide the flow cell priming port cover clockwise to open the priming port. 30s
- 143 After opening the priming port, check for a small air bubble under the cover. Drawback a small volume to remove any bubbles by following these sub-steps 30s
- Note**

Visually check for a continuous buffer from the priming port across the sensor array.
- 143.1 Set a P1000 pipette to 200 μ L 30s
- 143.2 Insert the tip into the priming port 30s
- 143.3 Turn the wheel until the dial shows 220-230 μ L , to draw back 20-30 μ L or until you can see a small volume of buffer entering the pipette tip. 30s



Note

IMPORTANT

Take care when drawing back the buffer from the flow cell. **Do not** remove more than **20-30 µl**, and make sure that the array of pores is covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.

144 Load  800 µL of the priming mix into the flow cell via the priming port, avoiding introducing air bubbles. Wait for  00:05:00 . During this time, prepare the library for loading by following the steps below.

5m



145 Thoroughly mix the contents of the **Library Beads (LIB)** by pipetting.

1m

146 In a new 1.5 m Eppendorf DNA LoBind tube, prepare the library for loading as follows:

2m

Reagent	Volume per flow cell
Sequencing Buffer (SB)	37.5 µL
Library Beads (LIB) mixed immediately before use	25.5 µL
DNA library	12 µL
Total	75 µL

147 **Complete the flow cell priming:**

30s

147.1 Gently lift the SpotON sample port cover to make the **SpotON** sample port accessible.

30s

147.2 Load  200 µL of the priming mix into the flow cell **priming port** (not the SpotON sample port), avoiding introducing air bubbles.

30s

148 Mix the prepared library gently by pipetting up and down just before loading.

30s

149 Add  75 µL of the prepared library to the flow cell via the **SpotON sample port** in a **dropwise fashion**. Ensure each drop flows into the port before adding the next.

2m



- 150 Gently replace **the SpotON sample port cover**, ensuring the bung enters the SpotON port. Close the priming port. Place the light shield onto the flow cell. Close the device lid and set up a sequencing run. Start sequencing.

5m

Protocol references

ZymoBIOMICS DNA/RNA Miniprep Kit: Microbiome DNA and RNA from any sample. Catalog no. R2002

Pan-Prokaryote riboPOOL Kit Manual v1.1 (2019)

HighPrep RNA Elite Manual v1.01. Catalog No. RC-90050

Ligation Sequencing gDNA Native Barcoding V14 SQK-NBD114.24 (2022) ONT

Ligation Sequencing v14: Direct cDNA Sequencing DCS_9187_v114 (2023) ONT

Agencourt AMPure XP Information For Use Guide PN B37419AB (2016)

Acknowledgements

The research's funding comes from the Research Council of Norway as part of the DeepSeaSequence project (**Project Number:** 315427).