

Jul 20, 2022

Version 1

Nanopore (SQK-LSK109) without barcode V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvob837l4o/v1>

Wen-Ting Zeng¹

¹Kaohsiung Medical College



Wen-Ting Zeng

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvob837l4o/v1>

Protocol Citation: Wen-Ting Zeng 2022. Nanopore (SQK-LSK109) without barcode. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5qpvob837l4o/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: May 21, 2022



Last Modified: September 13, 2023

Protocol Integer ID: 62940

Keywords: barcode nanopore, nanopore, lsk109

Abstract

Nanopore (SQK-LSK109) without barcode



DNA repair and end-prep

1h

- 1 In a 200µl PCR tube, mix the following:

A	B	C
Reagent	Volume	Color
49~100fmol sample DNA + DI Water	50µl	
End Repair&A-Tailing Buffer (KAPA)	7µl	Purple
End Repair&A-Tailing Enzyme Mix (KAPA)	3µl	Purple
Total	60µl	

- 1.1 Ensure the components are thoroughly mixed by pipetting, and spin down.

- 1.2 Using a thermal cycler, incubate at 20 °C for 00:30:00 and 65 °C for 00:30:00 . Hold at 4 °C

1h

1.2

Adapter ligation and clean-up

30m

- 2 Spin down the Adapter Mix (AMX) and Quick T4 Ligase, and place on ice.

Note

Although the recommended 3rd party ligase is supplied with its own buffer, the ligation efficiency of Adapter Mix (AMX) is higher when using Ligation Buffer supplied within the Ligation Sequencing Kit.

- 2.1 Thaw Ligation Buffer at room temperature, spin down and mix by vortex.
*Place on ice immediately after thawing and mixing.

2.2 In a 200µl PCR tube, mix in the following order:

A	B	C
Reagent	Volume	Color
End repaired & A-tailing product	60µl	
Adapter Mix (AMX)	5µl	Green
DI Water	5µl	
Ligation Buffer (KAPA)	30µl	Yellow
NEBNext Quick T4 DNA Ligase (KAPA)	10µl	Yellow
Total	110µl	

2.3 Ensure the components are thoroughly mixed by pipetting, and spin down.

2.4 Using a thermal cycler, incubate at  20 °C  00:30:00 .

30m

Note

During incubation, take out LFB, SQB, FLT, FLB from the fridge 30 minutes earlier to thaw on ice.

3

A	B	C
Reagent	Volume	Color
Ligation product	110µl	
Agencourt AMPure XP beads	44µl	Stored in 4°C
Long Fragment Buffer (LFB) or Short Fragment Buffer (SFB)	250µl (twice)	Orange & Grey
Elution Buffer from the Oxford Nanopore kit (EB)	14µl	Black

- 3.1 Transfer the DNA sample (110ul) to a clean 1.5 ml tube.
- 3.2 Resuspend the AMPure XP beads by vortexing.
- 3.3 Add 44 μ l (0.4X) of resuspended AMPure XP beads to the adapter-ligated reaction and mix by pipetting.
- 3.4 Incubate on a Hula mixer (rotator mixer) for 5 minutes at room temperature.
- 3.5 Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.

Note

Pellet the beads on a magnet until the eluate is clear and colorless, for at least 1 minute.

- 3.6 Wash the beads by adding either 250 μ l Long Fragment Buffer (LFB) or 250 μ l Short Fragment Buffer (SFB).

Note

Depending on the wash buffer (LFB or SFB) used, the clean-up step after adapter ligation is designed to either enrich for DNA fragments of >3 kb, or purify all fragments equally.

- 3.7 Flick the beads to resuspend, spin down, then return the tube to the magnetic rack and allow the beads to pellet.

Note

Pellet the beads on a magnet until the eluate is clear and colorless, for at least 1 minute.

- 3.8 Remove the supernatant using a pipette and discard.



3.9 Repeat the previous step 3.6.

3.10 Spin down and place the tube back on the magnet. Pipette off any residual supernatant.

Note

Pellet the beads on a magnet until the eluate is clear and colourless, for at least 1 minute.

3.11 Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.

3.12 Remove the tube from the magnetic rack and resuspend the pellet in 14 µl Elution Buffer (EB).

3.13 Spin down and incubate for 10 minutes at room temperature.

Note

For high molecular weight DNA, incubating at 37°C can improve the recovery of long fragments

3.14 Remove and retain 14 µl of eluate containing the DNA library into a clean 1.5 ml Eppendorf DNA tube.

3.15 Quantify 1 µl of eluted sample using a Qubit fluorometer.

Note

The prepared library is used for loading into the flow cell. Store the library on ice until ready to load.

Priming and loading the SpotON flow cell

4 Thaw the flow cell, Sequencing Buffer (SQB), Loading Beads (LB), Flush Tether (FLT) and one tube of Flush Buffer (FB) at room temperature.

4.1 Open the GridION and slide the flow cell under the clip.

Note

Press down firmly on the flow cell to ensure correct thermal and electrical contact.




4.2 Check flow cell (Check the pore)

Note

This takes about 10 minutes

4.3 Mix the Sequencing Buffer (SQB), Flush Tether (FLT) and Flush Buffer (FB) tubes by pipetting and spin down at room temperature.

4.4 Remove the flow cell from the machine and slide the **priming port cover** clockwise to open the priming port.



4.5 After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles (a few μl):

1. Set a P1000 pipette to 200 μl
2. Insert the tip into the priming port
3. Turn the wheel until the dial shows 220-230 μl , or until you can see a small volume of buffer entering the pipette tip

Note

Visually check that there is continuous buffer from the priming port across the sensor array.

4.6 To prepare the flow cell priming mix, add 30 μl of thawed and mixed Flush Tether (FLT) directly to the tube of thawed and mixed Flush Buffer (FB), and mix by vortexing at room temperature.

	A	B	C
	Reagent	Volume	Color
	Flush Tether (FLT)	30 μl	Purple
	Flush Buffer (FB)	New one	Blue



	A	B	C
	Total	1.2ml	

- 4.7 Load 800 μ l of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for 5 minutes. During this time, prepare the library for loading by following the steps below.

Note

Closed the priming port

- 4.8 Thoroughly mix the contents of the Loading Beads (LB) by flick.

- 4.9 In a new tube, prepare the library for loading as follows:


	A	B	C
	Reagent	Volume	Color
	Sequencing Buffer (SQB)	37.5 μ l	Red
	Loading Beads (LB), mixed immediately before use	25.5 μ l	Pink
	DNA library	12 μ l	
	Total	75 μ l	

- 4.10 Complete the flow cell priming:
1. Gently lift the SpotON sample port cover to make the SpotON sample port accessible.
 2. Load 200 μ l of the priming mix into the flow cell via the priming port (not the SpotON sample port), avoiding the introduction of air bubbles.

Note

Load the library as soon as possible after this step.

- 4.11 Mix the prepared library gently by pipetting up and down just prior to loading.
- 4.12 Add 75 μ l of sample to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next.
- 4.13 Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the GridION.



Flow Cell Wash


30m**5**

Place the tube of Wash Mix (WMX) on ice. Do not vortex the tube.

- 5.1 Thaw one tube of Wash Diluent (DIL) at room temperature.
- 5.2 Mix the contents of Wash Diluent (DIL) thoroughly by vortexing, spin down briefly and place on ice.
- 5.3 In a clean 1.5 ml Eppendorf DNA tube, prepare the following Flow Cell Wash Mix:

	A	B	C
Component	Volume	Color	
Wash Mix (WMX)	2µl	Brown	
Wash Diluent (DIL)	398µl	Brown	
Total	400µl		


- 5.4 Mix well by pipetting, and place on ice. Do not vortex the tube.
- 5.5 Stop or pause the sequencing experiment in MinKNOW, and leave the flow cell in the device.
- 5.6 Ensure that the priming port cover and SpotON sample port cover are in the positions indicated in the figure below.
- 5.7 Using a P1000, remove all fluid from the waste channel through Waste port 1. As both the priming port and SpotON sample port are closed, no fluid should leave the sensor array area.



- 5.8 Rotate the flow cell priming port cover clockwise so that the priming port is visible.
- 5.9 Check for air between the priming port and the sensor array. If necessary, using a P1000 draw back a small volume to remove any air (a few μ ls):
1. Set a P1000 pipette to 200 μ l
 2. Insert the tip into the priming port
 3. Turn the wheel until the dial shows 220-230 μ l, or until you can see a small volume of buffer/liquid entering the pipette tip.
 4. Visually check that there is continuous buffer from the priming port across the sensor array.
- 5.10 Load 400 μ l of the prepared Flow Cell Wash Mix into the flow cell via the priming port, avoiding the introduction of air.
- 5.11 Close the priming port and wait for 30 minutes.
- 5.12 Ensure that the priming port cover and SpotON sample port cover are in the positions indicated in the figure below.
- 5.13 Using a P1000, remove all fluid from the waste channel through Waste port 1. As both the priming port and SpotON sample port are closed, no fluid should leave the sensor array area.

Note

It is vital that the flow cell priming port and SpotON sample port are closed to prevent air from being drawn across the sensor array area, which would lead to a significant loss of sequencing channels.



Store the MinION/GridION flow cell for later use


30m

- 6 Thaw one tube of Storage Buffer (S) at room temperature.
- 6.1 Mix contents thoroughly by pipetting and spin down briefly.
- 6.2 Rotate the flow cell priming port cover clockwise so that the priming port is visible.
- 6.3 Check for air between the priming port and the sensor array. If necessary, using a P1000 draw back a small volume to remove any air (a few μ ls):
 1. Set a P1000 pipette to 200 μ l
 2. Insert the tip into the priming port
 3. Turn the wheel until the dial shows 220-230 μ l, or until you can see a small volume of buffer/liquid entering the pipette tip.
 4. Visually check that there is continuous buffer from the priming port across the sensor array.
- 6.4 Slowly add 500 μ l of Storage Buffer (S) through the priming port of the flow cell.
- 6.5 Close the priming port.

- 6.6 Using a P1000, remove all fluid from the waste channel through Waste port 1. As both the priming port and SpotON sample port are closed, no fluid should leave the sensor array area.

Note

It is vital that the flow cell priming port and SpotON sample port are closed to prevent air from being drawn across the sensor array area, which would lead to a significant loss of sequencing channels.



- 6.7 The flow cell can now be stored at 4-8°C.
- 6.8 When you wish to reuse the flow cell, remove the flow cell from storage, and allow it to warm to room temperature for ~5 minutes. You will need to perform a Flow Cell Check before loading the next library.



Note

We recommend storing libraries in Eppendorf DNA LoBind tubes at 4°C for short term storage or repeated use, for example, re-loading flow cells between washes.

For single use and long term storage of more than 3 months, we recommend storing libraries at -80°C in Eppendorf DNA LoBind tubes.

For further information, please refer to the DNA library stability Know-How document.