

Apr 12, 2023

NASC-seq2 Protocol

 [Nature Cell Biology](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr43nogmk/v1>



Gert-Jan Hendriks¹, Daniel Ramsköld¹, Anton Larsson¹, Juliane Mayr¹, Christoph Ziegenhain¹, Leonard Hartmanis¹, rickard.sandberg¹, Michael Hagemann-Jensen¹

¹Karolinska Institute Stockholm



Gert-Jan Hendriks

Karolinska Institute Stockholm

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr43nogmk/v1>

External link: <https://doi.org/10.1038/s41556-024-01486-9>

Protocol Citation: Gert-Jan Hendriks, Daniel Ramsköld, Anton Larsson, Juliane Mayr, Christoph Ziegenhain, Leonard Hartmanis, rickard.sandberg, Michael Hagemann-Jensen 2023. NASC-seq2 Protocol. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.6qpvr43nogmk/v1>

**Manuscript citation:**

Daniel Ramsköld, Gert-Jan Hendriks, Anton J. M. Larsson, Juliane V. Mayr, Christoph Ziegenhain, Michael Hagemann-Jensen, Leonard Hartmanis, Rickard Sandberg (2024) Single-cell new RNA sequencing reveals principles of transcription at the resolution of individual bursts. *Nature Cell Biology* doi: [10.1038/s41556-024-01486-9](https://doi.org/10.1038/s41556-024-01486-9)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 21, 2022

Last Modified: May 18, 2026

Protocol Integer ID: 73034

Keywords: seq2 protocol insights into transcriptional bursting kinetic, transcriptional bursting kinetic, analyses of static rna count, static rna count, time nascent rna imaging, rna in cell, rna, cell profiling, single cell, bursting dynamic, cell, seq2 protocol insight, seq2, nasc

Funders Acknowledgements:

GJH - Human Frontier Science Program

Grant ID: LT000155/2017-L

CZ - European Molecular Biology Organization

Grant ID: ALTF 673-2017

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Insights into transcriptional bursting kinetics and regulation have emerged from real-time nascent RNA imaging and analyses of static RNA counts over cells. Here, we developed sensitive single-cell profiling of newly transcribed (or new) RNA in cells (NASC-seq2) that can easily be applied to tens of thousands of single cells to help shed new light on bursting dynamics and coordination.

Image Attribution

Illustration by Beata Edyta Mierzwa (www.BeataScienceArt.com, social media: @beatascienceart)

Materials

EQUIPMENT

- Single-channel and multichannel pipettes
- Contact liquid handler for sample transfer and Vapor-Lock dispensing (Such as Agilent Bravo)
- Nanodispenser (such as Dispensix I.DOT, I.DOT Mini or Formulatrix Mantis)
- Plate centrifuge
- Tube centrifuge
- Fluorescence plate reader (compatible with Quantifluor dsDNA dye or similar)
- PE200 Sequencing (such as MGI G400 or Illumina NovaSeq)
- Magnetic Rack (384-well plate)
- Magnetic Rack (1.5ml Eppendorf tubes)
- 384-well compatible thermal cyclers

MATERIALS

General Protocol

⊗ Armadillo PCR Plate, 384-well, clear, clear wells **Thermo Fisher Catalog #AB2384**

⊗ Vapor-Lock **Qiagen Catalog #981611**

⊗ Adhesive sealing sheets **Thermo Scientific Catalog #AB0558**

⊗ Filter Tips

⊗ dimethylsulfoxide (DMSO) **Merck MilliporeSigma (Sigma-Aldrich)**

⊗ Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787-50ML**

⊗ UltraPure DNase/RNase-Free Distilled Water **Thermo Fisher Scientific Catalog #10977023**

⊗ DNA LoBind Tubes **Eppendorf Catalog ##022431021**

⊗ 4-thiouridine (4sU) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T4509**

⊗ Tris-HCL (pH 8.4 1M) **Merck MilliporeSigma (Sigma-Aldrich)**

⊗ Iodoacetamide (single-use vial 56mg) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A3221**

⊗ Tris-HCl, pH 8.0 (UltraPure) **Thermo Fisher Scientific Catalog #15568025**

⊗ Sodium Chloride (5M) **Invitrogen - Thermo Fisher Catalog #AM9760G**

⊗ GTP (Tris buffered solution 100mM) **Thermo Scientific Catalog #R1461**

⊗ Magnesium Chloride (1M Solution) **Invitrogen - Thermo Fisher Catalog #AM9530G**

⊗ Poly Ethylene Glycol (PEG) 8000 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #89510-250G-F**

⊗ Dithiothreitol (DTT) **Thermo Fisher Scientific Catalog #707265ML**

- ⊗ Recombinant RNase Inhibitor **Takara Bio Inc. Catalog #2313A**
- ⊗ Maxima H Minus Reverse Transcriptase (200 U/μL) **Thermo Fisher Catalog #EP0751**
- ⊗ KAPA HiFi Hotstart PCR kit **Roche Catalog #KK2502**
- ⊗ dNTP Set 100 mM Solutions **Thermo Fisher Scientific Catalog #R0182**
- ⊗ QuantiFluor(R) dsDNA System **Promega Catalog #E2670**
- ⊗ Agilent High Sensitivity DNA Kit **Agilent Technologies Catalog #5067-4626**
- ⊗ Tris-HCl pH 7.5
- ⊗ Nextera XT sample prep kit, 96 samples **Illumina, Inc. Catalog #FC-131- 1096**
- ⊗ NN-Dimethylformamide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D4551**
- ⊗ SDS, 10% Solution **Life Technologies Catalog #AM9822**
- ⊗ Phusion High-Fidelity DNA Polymerase (2 U/μL) **Thermo Fisher Catalog #F530L**

When using AmpureXP beads

- ⊗ Agencourt AMPure XP **Beckman Coulter Catalog #A63880**

When using DIY beads

- ⊗ Sera-Mag Speed Beads **GE Healthcare Catalog #65152105050250**
- ⊗ Sodium Azide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S2002-100G**
- ⊗ IGEPAL-CA630 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #I3021 SIGMA-ALDRICH**
- ⊗ EDTA (0.5 M), pH 8.0 **Life Technologies Catalog #AM9260G**

PRIMERS

All primers are the same as used in Smart-seq3 and were HPLC purified, with the TSO below being RNase-free HPLC purified.

- Smartseq3_OligodT30VN
/5Biosg/ACGAGCATCAGCAGCATACGATTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTT
- Smartseq3_N8_TSO
/5Biosg/AGAGACAGATTGCGCAATGNNNNNNNNrGrGrG
- Fwd_PCR_primer
TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGATTGCGCAA*T*G
- Rev_PCR_primer
ACGAGCATCAGCAGCATAC*G*A



Citation

Hagemann-Jensen M, Ziegenhain C, Chen P, Ramsköld D, Hendriks GJ, Larsson AJM, Faridani OR, Sandberg R (2020). Single-cell RNA counting at allele and isoform resolution using Smart-seq3.

<https://doi.org/10.1038/s41587-020-0497-0>

LINK

OPTIONAL

- Molecular Spike synthetic RNAs that were *in vitro* transcribed in the presence of 4sUTP (NU-1156S, Jena Bioscience).

Citation

Christoph Ziegenhain, Gert-Jan Hendriks, Michael Hagemann-Jensen, Rickard Sandberg (2022)
. Molecular spikes: a gold standard for single-cell RNA counting. Nature Methods.

<https://doi.org/10.1038/s41592-022-01446-x>

LINK

Prepare lysis plates

- 1 Using an automated pipetting platform or multichannel pipette, dispense 3uL of Vapor-Lock (Qiagen) to each well of a 384-well plate.

Note:

Do not use a non-contact dispenser (such as Dispensix IDOT or Formulatrix Mantis) for this step as the Vapor-Lock may severely damage the dispenser. We do this using an Agilent Bravo platform.

- 2 Prepare the following lysis buffer on ice.

Reagent	Final concentration	1 reaction	384-well plate
Recombinant RNase Inhibitor (40 U/uL)	2.5 u/uL	0.0188 uL	8.66 uL
4sU containing spUMI pool (OPTIONAL, 0.01 ng/uL)	0.04 pg/uL	0.0012 uL	0.553 uL
Triton-X100 (2 %)	0.1%	0.015 uL	6.91 uL
Nuclease-free water	-	0.265 uL	122.1 uL
Total		0.3 uL	138.2 uL

 On ice

- 3 Use a nanodispenser to distribute 0.3 uL of freshly prepared lysis buffer to each well of vapor-lock containing 384-well PCR plate.

 0.3 µL Lysis buffer

 On ice

- 4 Spin down at >3,000 G for 10 seconds.

10s

3000 x g, 00:00:10

5 Optional: Store lysis plate.

If not immediately continuing with the next step, prepared lysis plates can be sealed with aluminium seals and stored at -80°C.

Cell culture and FACS sorting

6 Grow cells in the presence of 50 µM 4sU. It can be helpful to collect untreated cells as a control.

Note

- 4sU is light-sensitive, and direct light (i.e. light in the cell culture hood, etc) should be reduced to a minimum.
- Do not refreeze leftovers from 4sU aliquot.
- Labeling times can vary for different celltypes and biological applications, but in general, short labeling times (i.e. less than 30 minutes) may result in poor separation of new and old molecules. Depending on your downstream analysis, this may affect your results.

7 Stop the labeling by transferring your cells to 15-ml falcon tubes on ice and wash them with cold PBS.

On ice

8 Distribute single cells into each well of the lysis plate by FACS sorting.

9 Spin down and seal the plates with aluminium seals and store them at -80°C.

Alkylation

10 Prepare the alkylation mix at room temperature.

Reagent	Final concentration	1 reaction	384-well plate
Tris-HCL (pH 8.4, 1 M)	50 mM	0.03 uL	13.8 uL

Reagent	Final concentration	1 reaction	384-well plate
DMSO (100 %)	45 %	0.24 uL	110.6 uL
IAA in DMSO (200 mM)	10 mM	0.03 uL	13.8 uL
Total		0.3 uL	138.2 uL

Please note that the calculations in this mix are calculated to the final volume in the alkylation reaction (i.e. 0.6uL). The DMSO final concentration listed above also accounts for the DMSO that is added with the IAA.

🔥 Room temperature

Note

To avoid problems with iodoacetamide stability and variability of alkylating potential, we use single-use vials (Sigma A3221-10VL) of iodoacetamide that are dissolved in DMSO at 200mM right before they are used to prepare the alkylation mix above and the remainder is discarded.

Safety information

Iodoacetamide should be handled in a fume hood and in accordance with local environmental health and safety regulations.

- 11 Using a nanodispenser, distribute 0.3 uL of freshly prepared alkylation mix to all wells of the 384-well plate containing cells and lysis-buffer.
- 12 Spin the 384-well plate down at >3,000 x G for 10 seconds.
- 13 Incubate the plate at 50C for 15 minutes. While this is running, prepare the quenching mix (see the next section).

🔥 50 °C 15 minutes

Quenching and Denaturation

**14** Prepare the quencing mix on ice.

Reagent	Final concentration	1 reaction	384-well plate
DTT (1 M in H ₂ O)	35 mM	0.035 uL	16.13 uL
dNTPs (10 mM each)	0.5 mM each	0.2 uL	92.16 uL
Oligo-dT (100 uM)	0.6 uM	0.024 uL	11.06 uL
Recombinant RNase Inhibitor (40 U/uL)	0.4 u/uL	0.04 uL	18.43 uL
Nuclease-free water	-	0.101 uL	46.54 uL
Total		0.4 uL	184.32

Note that the dNTPs, oligo-dT and RRI concentrations above are calculated to the final volume in the Reverse Transcription mix (4uL) and not the Quenching mix (1uL).

15 Using a nanodispenser, distribute 0.4 uL of freshly prepared quenching and denaturation mix to all wells of the 384-well plate.**16** Spin the 384-well plate down at >3,000 x G for 10 seconds**17** Incubate the plate for 5 minutes at room temperature, followed by 10 minutes at 72C and a final hold at 4C. While this is running, prepare the Reverse Transcription mix (see the next section)

Reverse Transcription

18 Prepare the Reverse Transcription mix as below

Reagent	Final concentration	1 reaction	384-well plate
Tris-HCL (1 M, pH 8.0)	25 mM	0,1 uL	46.08 uL
NaCl (1 M)	35 mM	0,14 uL	64.51 uL
GTP (100 mM)	1 mM	0,04 uL	18.43 uL
MgCl ₂ (1 M)	2.5 mM	0.01 uL	4.61 uL
PEG (40 %)	5 %	0.5 uL	230.4 uL



Reagent	Final concentration	1 reaction	384-well plate
DTT (1 M)	2 mM + carry-over from quenching and denaturation mix	0.008 uL	3.69 uL
Recombinant RNase Inhibitor (40 U/uL)	0.4 U/uL + carry-over from lysis as well as quenching and denaturation mix	0.04 uL	18.43 uL
Template Switching Oligo (1 mM)	2 uM	0.008 uL	3.69 uL
Maxima H-minus RT enzyme (200 U/uL)	2 U/uL	0.04 uL	18.43 uL
H2O	-	2.11 uL	974.13 uL
Total		3 uL	1382.4 uL

- 19 Dispense 3 uL of the freshly prepared Reverse Transcription mix to all wells of the 384-well plate.
- 20 Spin the 384-well plate down at >3,000 x G for 10 seconds
- 21 Place the plate in the thermal cycler and start the Reverse Transcription program

Preamplification PCR

- 22 For the remainder of the protocol we are using standard Smart-seq3 reaction conditions.

Citation

Hagemann-Jensen M, Ziegenhain C, Chen P, Ramsköld D, Hendriks GJ, Larsson AJM, Faridani OR, Sandberg R (2020). Single-cell RNA counting at allele and isoform resolution using Smart-seq3. Nature biotechnology.

<https://doi.org/10.1038/s41587-020-0497-0>

LINK

A detailed description of the Smart-seq3 protocol is also available on protocols.io.

Protocol



NAME

Smart-seq3 Protocol

CREATED BY

Michael Hagemann-Jensen

Preview

For convenience, the preamplification PCR protocol is also included below.

- 23 Prepare the preamplification PCR mix as below

	Reagent	Final concentration	1 reaction	384-well plate
	Kapa HiFi HotStart buffer (5 X)	1 X	2.0 uL	820 uL
	dNTPs (25 mM/each)	0.3 mM/each	0.12 uL	49.2 uL
	MgCl ₂ (100 mM)	0.5 mM	0.05 uL	20.5 uL
	Fwd Primer (100 uM)	0.5 uM	0.05 uL	20.5 uL
	Rev Primer (100 uM)	0.1 uM	0.01 uL	4.1 uL
	Polymerase (1 U/uL)	0.02 U/uL	0.2 uL	82 uL
	Nuclease Free Water	-	3.57 uL	1463.7 uL
	Total		6 uL	2460 uL

Preamplification PCR mix from Hagemann-Jensen et al.

- 24 Dispense 6 uL of the freshly prepared preamplification PCR mix to all wells of the 384-well plate.

25 Spin the 384-well plate down at $>3,000 \times G$ for 10 seconds

26 Place the plate in the thermal cycler and start the following PCR program

	A	B	C	D
	Step	Temperature	Time	Cycles
	Initial denaturation	98 °C	3 min	1x
	Denaturation	98 °C	20 sec	20-25x
	Annealing	65°C	30 sec	
	Elongation	72 °C	4 min	
	Final Elongation	72 °C	5 min	1x
	Hold	4 °C	Hold	

Purification, QC and Tagmentation

27 For details on the purification of the amplified cDNA, the quality control and the tagmentation to produce the final sequencing libraries, please refer to the Smart-seq3 protocols.io.

Protocol



NAME

Smart-seq3 Protocol

CREATED BY

Michael Hagemann-Jensen

Preview



Sequencing and data analysis

- 28 We have chosen to perform relatively long short-read sequencing with most data produced for NASC-seq2 being sequenced PE200. We highly recommend this as it will increase your ability to call individual molecules as either new (labeled) or old (unlabeled).

For details on data processing and analysis, please see the lab github:

<https://github.com/sandberg-lab/NASC-seq2>

Citations

Christoph Ziegenhain, Gert-Jan Hendriks, Michael Hagemann-Jensen, Rickard Sandberg. Molecular spikes: a gold standard for single-cell RNA counting

<https://doi.org/10.1038/s41592-022-01446-x>