

May 13, 2023

Version 2

NCEM Drop - Conventional (TM-014) V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.36wgqjx95vk5/v2>

sandra.crameri¹

¹CSIRO ACDP



sandra.crameri

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgqjx95vk5/v2>

Protocol Citation: sandra.crameri 2023. NCEM Drop - Conventional (TM-014). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.36wgqjx95vk5/v2> Version created by [sandra.crameri](#)

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 13, 2023



Last Modified: May 13, 2023

Protocol Integer ID: 81821

Keywords: standard routine conventional negative contrast em, ncem drop, ncem, conventional

Abstract

Standard routine conventional negative contrast EM.

Protocol materials



nano-W **Catalog** ##2018-5ML



HEADER

1 SAN:

SPEC No:

OPERATOR:

Conventional

11m

2 Adsorb  10 μ L sample to grid  00:10:00 , inspect to ensure sample does not dry out.

10m

3 Drain excess sample from grid using filter paper, leave wet.

4 Stain  nano-W **Catalog ##2018-5ML**  00:01:00

1m

5 Drain & dry using filter paper