

Dec 14, 2022

nCoV-2019 Illumina Miniseq sequencing protocol (2,000bp amplicon)



Forked from [nCoV-2019 sequencing protocol \(RAPID barcoding, 1200bp amplicon\)](#).

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvrpjbgbmk/v1>

Bruno Gomez-Gil¹, Juli²

¹CIAD AC Mazatlan; ²CIAD

Juli : Julissa Enciso Ibarra;



Bruno Gomez-Gil

CIAD

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvrpjbgbmk/v1>

Protocol Citation: Bruno Gomez-Gil, Juli 2022. nCoV-2019 Illumina Miniseq sequencing protocol (2,000bp amplicon). protocols.io <https://dx.doi.org/10.17504/protocols.io.6qpvrpjbgbmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: February 28, 2021

Last Modified: December 14, 2022

Protocol Integer ID: 47758

Keywords: sequencing protocol, illumina miniseq platform, nextera xt, amplicon, tagmentation with nextera xt, good sequence, protocol, much of this protocol

Funders Acknowledgements:

CONACYT

Grant ID: 321122

Disclaimer

It work for us, any modification is up to you.

Abstract

This is a fork of the protocol <https://dx.doi.org/10.17504/protocols.io.bh7hj9j6> but modified for tiled 2000bp amplicons, tagmentation with Nextera XT, indexing, and sequencing with the Illumina Miniseq platform.

It has already produced very good sequences.

Much of this protocol is base on this paper: <https://doi.org/10.1093/biomethods/bpaa014>

Guidelines

Tested with high viral copy numbers (<30 Ct).

Materials

cDNA

- GoScript™ Reverse Transcriptase Kit by Promega.

Multiplex DNA

-

Library preparation

- Nextera XT DNA Library Preparation Kit

-

Protocol materials

 GoScript™ Reverse Transcriptase Kit **Promega Catalog #A5001**

 KAPA Taq HotStart **Merck MilliporeSigma (Sigma-Aldrich) Catalog #KK1510**

 Nextera XT DNA Library Preparation Kit **Illumina, Inc. Catalog #FC-131-1024**

Safety warnings

 Please follow standard health and safety guidelines when working with COVID-19 patient samples.



Sample preparation

10m

- 1 Dilute the sample depending on the Ct values, this will reduce the likelihood of PCR-inhibition.

Ct range	Sample	Water
----------	--------	-------

12-15	2 µL	198 µL
16-18	2 µL	18 µL
>18	no dilution	

cDNA preparation

2h

- 2 We use the GoScript™ Reverse Transcriptase Kit **Promega Catalog #A5001**
Mix (pipetting) the following components in Eppendorf tube.

5m

Component	Volume
-----------	--------

Nuclease-free water	4 µL
GoScript™ Reaction Buffer, Random Primer	2 µL
GoScript™ Enzyme Mix	0.4 µL
Final volume	10 µL

Note

It is good practice to carry a **negative control** (e.g. water) through the entire process from cDNA preparation to sequencing.

Note

The mastermix should be made up in the **mastermix cabinet** and aliquoted into PCR strip tubes. Tubes should be wiped down when entering and leaving the mastermix cabinet.

- 3 Prepare the **mastermix** on ice, mix by pipetting.

**Component****Volume**

GoScript™ Reverse Transcription Mix

10 µL

RNA

3 µL

Final volume

13 µL

4 Incubate the reaction as follows:

1h 21m

Anneal primers

25 °C

for

00:05:00

Extension

42 °C

for

01:00:00

Inactivation

70 °C

for

00:15:00

Snap cool in a prechilled metal rack or on ice 00:01:00

Note

A quick cooling step using a PCR cooling block or ice helps to inhibit secondary structure formation and can decrease variation in overall coverage.

Note

A mastermix should be made up in the **mastermix cabinet** and added to the denatured RNA in the **extraction and sample addition cabinet**. Tubes should be wiped down when entering and leaving the mastermix cabinet.

Primer pool preparation

5 **PRIMERS** for this protocol are described in protocol [nCoV-2019 sequencing protocol \(RAPID barcoding, 1200bp amplicon\)](https://dx.doi.org/10.17504/protocols.io.6gpvrdpjbghmk/v1) and are the 2,000 bp option.

We selected for 2,000 bp amplicons so they could be more easily tagged in the library preparation.

Two pools are prepared, Pool1 with 18 primers and Pool2 with 16 primers.

POOL1. Average T_m 61.01 °C



A	B	C	D	E	F	G	H	I	J	K	L
Name	Sequence	Direction	Start	End	Length	Product Size	Tm	%G C	Hairpin Tm	Pair Dimers Tm	Self Dimers Tm
SARSC oV_200 0_01_LEFT	ACCAACC AACTTTC GATCTCT TGT	forward	31	54	24	20 49	60. 7	41. 7	No ne	No ne	No ne
SARSC oV_200 0_01_RIGHT	ACACCAC CTGTAATG TAGGCCA	reverse	20 58	20 79	22	20 49	61. 4	50	39	No ne	5
SARSC oV_200 0_03_LEFT	TCGCACA AATGTCTA CTTAGCT GT	forward	37 72	37 95	24	181 4	60. 6	41. 7	No ne	0.7	No ne
SARSC oV_200 0_03_RIGHT	GTGTGCC CATGTACA TAACAGC T	reverse	55 63	55 85	23	181 4	61. 2	47. 8	42. 7	0.7	8.3
SARSC oV_200 0_05_LEFT	CAATCAT GCAATTG TTTTTCA GCTATTTT G	forward	72 99	73 28	30	182 5	60. 4	30	33	No ne	7.4
SARSC oV_200 0_05_RIGHT	CGTGTGT CAGGGCG TAACTTT	reverse	91 02	91 23	22	182 5	61. 6	50	No ne	No ne	No ne
SARSC oV_200 0_07_LEFT	GGACGTA CCATATTG GGTAGTG C	forward	10 88 6	10 90 8	23	183 4	60. 8	52. 2	43. 8	No ne	No ne
SARSC oV_200 0_07_RIGHT	TCTGTCTG TAGTGCA ACAGGAC T	reverse	126 98	127 19	22	183 4	61. 3	50	41. 4	No ne	14
SARSC oV_200 0_09_LEFT	ACCACTT CAGAGAG CTAGGTG T	forward	14 47 7	14 49 8	22	192 5	60. 5	50	No ne	No ne	No ne
SARSC oV_200 0_09_RIGHT	ACAACCT GGAGCAT TGCAAAC A	reverse	16 38 0	16 40 1	22	192 5	61. 5	45. 5	No ne	No ne	11.9

	A	B	C	D	E	F	G	H	I	J	K	L
	SARSC oV_200 0_11_LE FT	TGGCATA CCTAAGG ACATGAC CT	for wa rd	181 68	181 90	23	181 4	60. 9	47. 8	38. 3	No ne	11.2
	SARSC oV_200 0_11_RI GHT	CAGTGAG TGGTGCA CAAATCG T	rev ers e	19 96 0	19 98 1	22	181 4	61. 6	50	38. 7	No ne	20. 4
	SARSC oV_200 0_13_LE FT	TCCTCAG TTTTACAT TCAACTC AGGA	for wa rd	216 95	217 20	26	193 7	60. 2	38. 5	45. 2	No ne	0.9
	SARSC oV_200 0_13_RI GHT	TGACTAG CTACACTA CGTGCCC	rev ers e	23 61 0	23 63 1	22	193 7	61. 5	54. 5	37	No ne	No ne
	SARSC oV_200 0_15_LE FT	AGGAGTC AAATTACA TTACACAT AAACGAA	for wa rd	25 36 0	25 38 9	30	180 5	60. 1	30	No ne	No ne	No ne
	SARSC oV_200 0_15_RI GHT	ACTGCTA CTGGAAT GGTCTGT GT	rev ers e	271 42	271 64	23	180 5	61. 6	47. 8	No ne	No ne	No ne
	SARSC oV_200 0_17_LE FT	ACTTGTC ACGCCTA AACGAAC A	for wa rd	27 87 3	27 89 4	22	191 8	60. 7	45. 5	36. 7	No ne	No ne
	SARSC oV_200 0_17_RI GHT	TAGGCAG CTCTCCC TAGCATTG	rev ers e	29 76 9	29 79 0	22	191 8	61. 6	54. 5	45. 3	No ne	No ne

Pool1

POOL2. Average Tm 🌡️ 61.05 °C

	A	B	C	D	E	F	G	H	I	J	K	L
	Name	Sequence	Dir ect ion	Sta rt	En d	Le ng th	Pro du ct Siz e	Tm	%G C	Hai rpi n Tm	Pai r Di me r Tm	Sel f Di me r Tm



	A	B	C	D	E	F	G	H	I	J	K	L
	SARSC oV_200 0_02_L EFT	AGGCCGC TATAACAA TACTAGAT GGA	for wa rd	195 6	19 81	26	192 3	61. 3	42. 3	No ne	No ne	No ne
	SARSC oV_200 0_02_RI GHT	CAGCGAT CTTTTGT TCAACTT GCT	rev ers e	38 55	38 78	24	192 3	60. 8	41. 7	No ne	No ne	No ne
	SARSC oV_200 0_04_L EFT	TCAACAT GCCAATT TAGATTCT TGCA	for wa rd	54 73	54 98	26	192 9	60. 3	34. 6	No ne	8.1	No ne
	SARSC oV_200 0_04_RI GHT	GCTGAAA TCGGGGC CATTTGTA	rev ers e	73 80	74 01	22	192 9	61. 5	50	No ne	8.1	4.4
	SARSC oV_200 0_06_L EFT	GCTGCTG AATGTACA ATTTTAA AGATGC	for wa rd	90 11	90 39	29	20 03	61.1	34. 5	No ne	No ne	No ne
	SARSC oV_200 0_06_RI GHT	AACCAGT GGTGTGT ACCCTTG A	rev ers e	10 99 2	110 13	22	20 03	61. 5	50	47	No ne	17.6
	SARSC oV_200 0_08_L EFT	TCACCTA ATTTAGCA TGGCCTC TT	for wa rd	126 20	12 64 3	24	195 6	60. 1	41. 7	No ne	No ne	3.3
	SARSC oV_200 0_08_RI GHT	CAGGGTC AGCAGCA TACACAA G	rev ers e	14 55 4	14 57 5	22	195 6	61. 5	54. 5	No ne	No ne	No ne
	SARSC oV_200 0_10_LE FT	TGCATAC GTAGACC ATTCTTAT GTTGT	for wa rd	162 91	16 31 7	27	198 5	60. 8	37	32. 8	No ne	0.1
	SARSC oV_200 0_10_RI GHT	GCTTCTT CGCGGGT GATAAAC A	rev ers e	182 54	18 27 5	22	198 5	61. 5	50	No ne	No ne	6.1
	SARSC oV_200 0_12_LE FT	GGACTAC AAAAGAG ATGCTCC AGC	for wa rd	19 87 8	19 90 1	24	192 0	61. 5	50	42. 5	11.2	No ne
	SARSC oV_200	ACCTCTT AGTACCAT	rev ers e	217 75	217 97	23	192 0	60. 5	47. 8	37.1	11.2	12. 7

	A	B	C	D	E	F	G	H	I	J	K	L
	0_12_RIGHT	TGGTCCC A										
	SARSC oV_200 0_14_LEFT	GCTGAAC ATGTCAA CAACTCA TATGA	for wa rd	23 519	23 54 4	26	197 3	60. 1	38. 5	35. 2	No ne	5.1
	SARSC oV_200 0_14_RIGHT	TGCAGTA GCGCGAA CAAAATC T	rev ers e	25 47 0	25 49 1	22	197 3	61. 4	45. 5	46. 2	No ne	12.1
	SARSC oV_200 0_16_LEFT	TCTTATTA CAAATTG GGAGCTT CGCA	for wa rd	27 051	27 07 6	26	199 5	61. 3	38. 5	37. 6	No ne	No ne
	SARSC oV_200 0_16_RIGHT	GCTTCTT AGAAGCC TCAGCAG C	rev ers e	29 02 4	29 04 5	22	199 5	61. 6	54. 5	43. 6	No ne	30. 4

Pool2.

6 PRIMER STOCKS (100 micromolar (μM))

If you have ordered each primer independently and need to generate primer pool stocks:

add 5 μL of each primer from Pool 1 to a 1.5 mL Eppendorf labeled "Pool 1

(100μM)" and each primer from Pool 2 to a 1.5 mL Eppendorf labelled "Pool 2

(100μM)". These are your 100 micromolar (μM) **stocks** of each primer pool.

Note

Primers should be diluted and pooled in the **mastermix cabinet** which should be cleaned with decontamination wipes and UV sterilised before and after use.

7 WORKING PRIMERS (10 micromolar (μM))

Dilute the primer stocks **1:10** in molecular grade water, to generate **10μM primer working stocks**. It is recommend that multiple aliquots of each primer pool are made to in case of degradation or contamination.

Note

Primers need to be used at a **final concentration of** [M] 0.015 micromolar (μM) **per primer.**

8 POOLING OF PRIMERS

Label two [M] 1.5 mL Eppendorf tubes, one as **POOL1** and the another as **POOL2**.

8.1 Add [M] 10 μL of each of the primers of **set 1** [⇒ go to step #7](#) to the Eppendorf tube labelled **POOL1**, final concentration will be [M] 10 micromolar (μM) each

8.2 Add [M] 10 μL of each of the primers of **set 2** [⇒ go to step #7](#) to the Eppendorf tube labelled **POOL2**, final concentration will be [M] 10 micromolar (μM) each

Multiplex PCR

9 In the PCR hood set up the multiplex mastermix reaction as follows in 2 0.2mL PCR tubes. We use the

 KAPA Taq HotStart **Merck MilliporeSigma (Sigma-Aldrich) Catalog #KK1510**

Component	Pool 1	Pool 2	Final Concentration
5X Kapa HotStart Buffer	[M] 2.5 μL	[M] 2.5 μL	[M] 1 X
25 mM MgCl_2	[M] 0.75 μL	[M] 0.75 μL	
[M] 1.5 millimolar (mM)			
10 mM dNTPs	[M] 0.25 μL	[M] 0.25 μL	
[M] 0.2 millimolar (mM) each			
Primer Pool 1 or 2 (10 μM ea)	[M] 0.6 μL	[M] 0.6 μL	
[M] 0.5 micromolar (μM)			
Kapa Taq HotStart Polymerase	[M] 0.1 μL	[M] 0.1 μL	[M] 0.5 U
Nuclease-free water	[M] 5.8 μL	[M] 5.8 μL	

**Final mastermix volume**

10.0 µL

10.0 µL

Note

A PCR mastermix for each pool should be made up in the **mastermix cabinet** and aliquoted into PCR strip tubes. Tubes should be wiped down when entering and leaving the mastermix cabinet.

- 10 In the **extraction and sample addition cabinet** add 2.5 µL cDNA to each tube and mix well by pipetting.

The final volume will be 12.5 µL

Note

The **extraction and sample addition cabinet** should be cleaned with decontamination wipes and UV sterilised before and after use.

- 11 Pulse centrifuge the tubes to collect the contents at the bottom of the tube.

- 12 Set-up the following program on the thermal cycler:

5m 45s

Step	Temperature	Time	Cycles
Heat Activation	98 °C	00:00:30	1
Denaturation	95 °C	00:00:15	25-35
Annealing and Extension	60 °C	00:05:00	25-35
Hold	4 °C	Indefinite	1

Note

Cycle number should be 25 for Ct 18-21 up to a maximum of 35 cycles for Ct 35. We typically use 30 cycles.



Note

It is advisable to check each PCR amplicon by electrophoresis (1% agarose). [MORE INFO](#)

Expected result

Final concentrations of PCR products can range from ~20- 150ng/ul.

Pooling and PCR quantification

- 13 Amplicon quantification to make an equimolar mixture.

Note

At this stage, care should be taken with amplified PCR products. Only open tubes in a designated post-PCR workspace with equipment that is separate from areas where primers and mastermixes are handled.

- 13.1 Put  1 μ L of each pool in a Nanodrop or similar spectrometer and quantify the DNA concentration.
- 13.2 Label a  1.5 mL Eppendorf tube for each sample and make a **equimolar mix** with the two pools.
Calculate to achieve a final concentration of  50 ng/ul
- 14 Quantify DNA using a **Qubit** or other method.
Quantification using Nanodrop is not recommended for a good estimation of the final pool.



Note

Alternatively, these amplicons can also be used for Illumina sequencing, such as found here: [x.doi.org/10.17504/protocols.io.betejeje](https://doi.org/10.17504/protocols.io.betejeje)

We have found that performing an Ampure XP bead clean up at this stage does not improve performance. Therefore, it is not necessary to clean up the PCR reaction at this step.

14.1 Protocol

Protocol

NAME

DNA quantification using the Qubit fluorometer

CREATED BY

Nikki Freed

Preview

- 14.1.1 Prepare a mastermix of Qubit™ working solution for the required number of samples and standards. The Qubit dsDNA kit requires 2 standards for calibration (see note below).

Per sample:

Qubit® dsDNA HS Reagent  1 µL

Qubit® dsDNA HS Buffer  199 µL

**Note**

If you have already performed a calibration on the Qubit machine for the selected assay you can use the previous calibration stored on the machine. We recommend performing a new calibration for every sample batch but a same-day calibration would be fine to use for multiple batches.

To avoid any cross-contamination, we recommend that you remove the total amount of working solution required for your samples and standards from the working solution bottle and then add the required volume to the appropriate tubes instead of pipetting directly from the bottle to each tube.

- 14.1.2 Label the tube lids. Do not label the side of the tube as this could interfere with the sample reading.

Note

Use only thin-wall, clear, 0.5mL PCR tubes. Acceptable tubes include Qubit™ assay tubes (Cat. No. Q32856)

- 14.1.3 Aliquot Qubit™ working solution to each tube:
- standard tubes requires 190µL of Qubit™ working solution
 - sample tubes require anywhere from 180–199µL (depending how much sample you wish to add).

The final volume in each tube must be 200µL once sample/standard has been added.

- 14.1.4 Add 10µL of standard to the appropriate tube.

- 14.1.5 Add 1–20µL of each user sample to the appropriate tube.

Note

If you are adding 1–2µL of sample, use a P-2 pipette for best results.

- 14.1.6 Mix each tube vigorously by vortexing for 3–5 seconds.



14.1.7 Allow all tubes to incubate at room temperature for 2 minutes, then proceed to “Read standards and samples”.

14.1.8 On the Home screen of the Qubit™ 3 Fluorometer, press DNA, then select 1X dsDNA HS as the assay type. The Read standards screen is displayed. Press Read Standards to proceed.

Note

If you have already performed a calibration for the selected assay, the instrument prompts you to choose between reading new standards and running samples using the previous calibration. **If you want to use the previous calibration, skip to step 12.** Otherwise, continue with step 9.

14.1.9 Insert the tube containing Standard #1 into the sample chamber, close the lid, then press Read standard. When the reading is complete (~3 seconds), remove Standard #1.

14.1.10 Insert the tube containing Standard #2 into the sample chamber, close the lid, then press Read standard. When the reading is complete, remove Standard #2.

14.1.11 The instrument displays the results on the Read standard screen. For information on interpreting the calibration results, refer to the Qubit™ Fluorometer User Guide, available for download at thermofisher.com/qubit.

14.1.12 Press Run samples.

14.1.13 On the assay screen, select the sample volume and units:

- Press the + or – buttons on the wheel, or anywhere on the wheel itself, to select the sample volume added to the assay tube (from 1–20µL).
- From the unit dropdown menu, select the units for the output sample concentration (in this case choose ng/µL).

14.1.14 Insert a sample tube into the sample chamber, close the lid, then press Read tube. When the reading is complete (~3 seconds), remove the sample tube.

14.1.15 **The top value (in large font) is the concentration of the original sample and the bottom value is the dilution concentration.** For information on interpreting the sample results, refer to the Qubit™ Fluorometer User Guide.

14.1.16 Repeat step 14 until all samples have been read.



14.1.17 Carefully **record all results** and store run file from the Qubit on a memory stick.

14.1.18 All negative controls should ideally be 'too low' to read on the Qubit machine, but **MUST** be < 1ng per ul. If your negative controls >1ng per ul, considerable contamination has occurred and you must redo previous steps.

Normalisation

15 Label a  0.2 mL PCR tube for each sample.

15.1 Adjust the amount of DNA in the tube to be  100 ng total per sample in  7.5 µL molecular grade water.

Note

For example if your PCR reaction is at least  100 ng/µl add 1ul of the PCR reaction to 6.5ul of molecular grade water. Use 7.5ul of the negative control, even if there is no detectable DNA in the PCR reaction.

Tagmentation

7m

16  Nextera XT DNA Library Preparation Kit **Illumina, Inc. Catalog #FC-131-1024** We have reduced the amount of reagents used per reaction.

16.1 Label a  1.5 mL Eppendorf tube for each sample and add in this order:

Component

Volume

1. Tagment DNA buffer (TD)  3.75 µL

2. cDNA  2 ng/µl  3 µL

Pipette 10 times to mix

16.2 Add **Amplicon Tagment Mix**  0.5 µL and pipette 10 times to mix



16.3 Centrifuge at  280 x g, 20°C, 00:01:00

16.4 Incubate in thermal cycler

 55 °C for  00:10:00

 10 °C indefinitely

10m

16.5 Add  2 µL of Neutralize Tagment Buffer (NT).

Pipette 10 times to mix

16.6 Centrifuge at  280 x g, 20°C, 00:01:00

16.7 Incubate at  Room temperature for  00:05:00  0 °C

5m

16.8 Preserve the samples at  4-8 °C until use.

Indexing

30m

17 La siguiente reacción requiere de la enzima polimerasa 2x Ampigene HS Taq Mix Catalog # ENZ-NUC101-0200 y de los sets de Nextera XT. Cada set contiene 96 combinaciones de TAGs con esto podemos alcanzar 384 muestras usando los Nextera XT Index Kit v2 (Sets A,B,C y D) Catalog # 20027213;20027214;20027215;20027216.

17.1 A cada tubo de reaccion agregar:

Component	Volume
Nuclease-free water	10 µL
2x Ampigene HS Taq Mix.	10 µL
Nextera XT index i5	1 µL
Nextera XT index i7	1 µL
ADN Tagmentado	3 µL

Final volume

25 µL

17.2 Centrifuge at 280 x g, 20°C, 00:01:00

1m

17.3 Set-up the following program on the thermal cycler:

4m 40s

Cover 100 °C

Step	Temperature	Time	Cycles
Heat Activation	72 °C	00:03:00	1
Denaturation	95 °C	00:00:30	1
	95 °C	00:00:10	
Annealing and Extension	55 °C	00:00:30	14
	72 °C	00:00:30	
Hold	4 °C	Indefinite	1

18

Purificación final y pooling

19 Al volumen que se tiene en el tubo agregar 0.8X de perlas magnéticas Ampure XP.

19.1 Mezclar y spin-down.

19.2 Llevar al magneto hasta formar el pellet y desechar el sobrenadante.

19.3 Agregar 150 uL de etanol 80% en posición contraria al pellet. Esperar 30 segundos. Desechar el etanol.

19.4 Retirar el exceso de etanol y secar las perlas por 5 min.

19.5 Resuspender las perlas en 26 uL de Agua libre de nucleasas.



Mix y spin down. Incubar por 5 min a TA.

19.6 Llevar al magneto hasta formar el pellet y transferir 25 uL del sobrenadante a un nuevo tubo previamente etiquetado.

19.7 Cuantificar 2 uL por Qubit HS y analizar los tamaños mediante electroforesis en gel de agarosa al 1.0%.

Secuenciacion de bibliotecas en plataforma Illumina Miniseq

20 A partir de la concentración en ng/uL determinada por Qubit y obtenido el tamaño aproximado del fragmento, llevar cada una de las librerías a una concetracion de 4 nM. Transferir 5 uL de cada librerías a 4 nM aun tubo previamente etiquetado para obtener el pool final.

Seguir el protocolo Library Denaturing and miniseq Sample Loading del kit