

Aug 24, 2020

nCoV-2019 McGill RT Protocol, Lunascript

Forked from a private protocol

DOI

dx.doi.org/10.17504/protocols.io.bjgekjte

Shu-Huang Chen¹, Sarah J Reiling¹, Josh Quick², Ioannis Ragoussis¹

¹McGill University; ²University of Birmingham

McGill Genome Centre



Sarah J Reiling

McGill University



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bjgekjte

Protocol Citation: Shu-Huang Chen, Sarah J Reiling, Josh Quick, Ioannis Ragoussis 2020. nCoV-2019 McGill RT Protocol, Lunascript. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bjgekjte>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 07, 2020

Last Modified: August 24, 2020

Protocol Integer ID: 40166

Keywords: Cov-19 Nanopore Sequencing, Cov-19 Illumina Sequencing,



Abstract

Artic nCoV-2019 McGill modified Lunascript Reverse Transcriptase sequencing protocol.

Materials

MATERIALS



Water, nuclease free







LunaScript RT SuperMix Kit **New England Biolabs Catalog # E3010L**

cDNA preparation

1

Mix the following components in a 0.2 mL 8-strip tube;

Component	Volume
LunaScript RT SuperMix (5x)	 4 µL
Nuclease-free water	 5 µL
Template RNA	 11 µL
Total	 20 µL

Note

Viral RNA input from a clinical sample should be between Ct 18-35. If Ct is between 12-15, then dilute the sample 100-fold in water, if between 15-18 then dilute 10-fold in water. This will reduce the likelihood of PCR-inhibition.

Note

A mastermix should be made up in the **mastermix cabinet**. Tubes should be wiped down when entering and leaving the mastermix cabinet.

2


Gently mix by pipetting and pulse spin the tube to collect liquid at the bottom of the tube.


3

Incubate the reaction as follows:

 25 °C for  00:02:00

 55 °C for  00:20:00

 95 °C for  00:01:00

Place on ice for  00:01:00 or store cDNA at -20C