May 13, 2020 Version 2

Cov-2019 sequencing protocol (RAPID barcoding, 1200bp amplicon) V.2

Version 1 is forked from <u>nCoV-2019 sequencing protocol v2 (Gunlt)</u>

DOI

dx.doi.org/10.17504/protocols.io.bgc8jszw

Nikki Freed¹, Olin Silander¹

¹Massey University

Coronavirus Method De...

Nikki Freed





DOI: dx.doi.org/10.17504/protocols.io.bgc8jszw

Protocol Citation: Nikki Freed, Olin Silander 2020. nCoV-2019 sequencing protocol (RAPID barcoding, 1200bp amplicon). protocols.io <u>https://dx.doi.org/10.17504/protocols.io.bgc8jszw</u>

License: This is an open access protocol distributed under the terms of the **<u>Creative Commons Attribution License</u>**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development We are still developing and optimizing this protocol

Created: May 13, 2020

Last Modified: May 13, 2020

Protocol Integer ID: 36992

] Therefore all Minute than a bound ha
"] "muchelletereleterette
"]

Abstract

To enable faster, easier sequencing of SARS-COV2 genomes with fewer steps than current methods, we use multiplexed 1200 base pair PCR amplicons with the Oxford Nanopore RAPID barcoding kit (RBK004).

This is a modification of the ARTIC amplicon V3 sequencing protocol for MinION for nCoV-2019 developed by Josh Quick, which produces 400 base pair amplicons and uses the Oxford Nanopore Ligation barcoding kit (LSK-109).

We have increased the size of the amplicons to 1200bp and use the RAPID barcode kit (RBK004), which enables requires less time and fewer reagents than the LSK-109 protocol. The amplicons produced in this protocol could also be used for Illumina sequencing.

Primers were all designed using Primal Scheme: <u>http://primal.zibraproject.org</u>/, described here <u>https://www.nature.com/articles/nprot.2017.066</u>.

Primer sequences are here:

https://docs.google.com/spreadsheets/d/1M5I_C56ZC8_2Ycgm9EFieVIVNqxsP7dXAnGoBZy3nDo/edit? usp=sharing

We can ship a small amount of pooled primers to interested labs for further testing, email freednikki@gmail.com or olinsilander@gmail.com

Guidelines

This has so far been testing using only five SARS-CoV2 patient positive samples, with Cq values ranging from 20 to 31. Further testing might be needed to test the method on low viral load samples/high Cq samples.

Materials

STEP MATERIALS

 Primers 25nm, desalted, ideally LabReady formulation from IDT: <u>https://docs.google.com/spreadsheets/d/1M5I_C56ZC8_2Ycgm9EFieVIVNqxsP7dXAnGoBZy3nDo/edit#gid=75</u> <u>5704891</u>

 Extraction kits; Zymo Quick-RNA Viral Kit 	Zymo	R1034
OR ■ i.e. QIAamp Viral RNA Mini	Qiagen	52904
SuperScript IV (50 rxn)dNTP mix (10 mM each)	Thermo Thermo	18090050 R0192
 Random Hexamers (50 μM) OR 	Thermo	N8080127
 Random Primer Mix (60 μM) 	NEB	S1330S
 RNase OUT (125 rxn) 	Thermo	10777019
 Q5 Hot Start HF Polymerase 	NEB	M0493S
Agencourt AMPure XP	Beckman Coulte	<u>er A63880</u>
Rapid Barcoding Kit 1-12	Nanopore	SQK-RBK004
 R9.4.1 flow cell 	Nanopore	FLO-MIN106

Protocol materials

🔀 SQK-RBK004 Rapid Barcoding Kit **Oxford Nanopore Technologies Catalog #**SQK-RBK004

Safety warnings

Please follow standard health and safety guidelines when working with COVID-19 patient samples.

cDN	NA preparation		5m
1	Mix the following components	in an 0.2mL 8-strip tube;	5m
	Component	Volume	
	50µM random hexamers 10mM dNTPs mix (10mM each Template RNA Total	Δ 1 μL Δ 1 μL Δ 11 μL Δ 13 μL	
	Note		
	then dilute the sample 100-fol will reduce the likelihood of PC	sample should be between Ct 18-35. If Ct is between 12-15, d in water, if between 15-18 then dilute 10-fold in water. This CR-inhibition. It is good practice to carry a negative control process from cDNA preparation to sequencing.	
	Note		
		up in the mastermix cabinet and aliquoted into PCR strip down when entering and leaving the mastermix cabinet.	
2	Gently mix by pipetting and pu	Ilse spin the tube to collect liquid at the bottom of the tube.	
3	Incubate the reaction as follow	/S:	6m
	₿ 65 °C for 00:05:00		
	Snap cool in a prechilled metal	l rack or on ice 🕥 00:01:00	
	Note		
	A quick cooling step using a P formation and can decrease v	PCR cooling block or ice helps to inhibit secondary structure ariation in overall coverage.	


4 Add the following to the annealed template RNA :

Component	Vol	ume	
SSIV Buffer		Д	4 μL
100mM DTT		<u></u> ⊿ 1	μL
RNaseOUT RNase Inhibitor	₫ 1	μL	
SSIV Reverse Transcriptase	₫ 1	μL	
Total		Д	20 µL

Note

A mastermix should be made up in the **mastermix cabinet** and added to the denatured RNA in the **extraction and sample addition cabinet**. Tubes should be wiped down when entering and leaving the mastermix cabinet.

- 5 Gently mix by pipetting and pulse spin the tube to collect liquid at the bottom of the tube.
- 6 Incubate the reaction in a preheated PCR machine:



Primer pool preparation

7 If required, resuspend lyophilised primers at a concentration of 100µM each

5m

1h 5m

Note

Primers for this protocol were designed using <u>Primal Scheme</u> and generate overlapping 1200nt amplicons. Primer names and dilutions are listed here:

https://docs.google.com/spreadsheets/d/1M5I_C56ZC8_2Ycgm9EFieVIVNqxsP7dXAnGo BZy3nDo/edit?usp=sharing

We have tested multiplexing 1500 nt and 2000 nt amplicons as well, all work well. These are included in the link. Here we will discuss just the protocol for 1200 nt amplicons as this longer amplicons might be sensitive to RNA degradation.

7.1 Primers used to generate 1200 bp amplicons are here:

Primer Name	Sequence	P o ol	Len gth	Tm	GC %	Start
SARSCoV_1200_1_LE FT	ACCAACCAACTTTCGATCTCTT GT	1	24	60.6 9	41.6 7	30
SARSCoV_1200_1_RIG HT	GGTTGCATTCATTTGGTGACGC	1	22	61.49	50	1205
SARSCoV_1200_3_LE FT	GGCTTGAAGAGAAGTTTAAGGA AGGT	1	26	61.19	42.3 1	2153
SARSCoV_1200_3_RI GHT	GATTGTCCTCACTGCCGTCTTG	1	22	61.5	54.5 5	3257
SARSCoV_1200_5_LE FT	ACCTACTAAAAAGGCTGGTGGC	1	22	60.5 5	50	4167
SARSCoV_1200_5_RI GHT	AGCATCTTGTAGAGCAGGTGGA	1	22	61.16	50	5359
SARSCoV_1200_7_LE FT	ACCTGGTGTATACGTTGTCTTTG G	1	24	60.8	45.8 3	6283
SARSCoV_1200_7_RI GHT	GCTGAAATCGGGGGCCATTTGTA	1	22	61.53	50	7401
SARSCoV_1200_9_LE FT	AGAAGTTACTGGCGATAGTTGTA ATAACT	1	29	60.5 9	34.4 8	8253
SARSCoV_1200_9_RI GHT	TGCTGATATGTCCAAAGCACCA	1	22	60.2 9	45.4 5	9400
SARSCoV_1200_11_LE FT	AGACACCTAAGTATAAGTTTGTT CGCA	1	27	60.74	37.0 4	1034 3
SARSCoV_1200_11_RI GHT	GCCCACATGGAAATGGCTTGAT	1	22	61.8	50	1146 9
SARSCoV_1200_13_L EFT	ACCTCTTACAACAGCAGCCAAA C	1	23	61.55	47.8 3	1245 0
SARSCoV_1200_13_RI GHT	CGTCCTTTTCTTGGAAGCGACA	1	22	61.38	50	1362 1
SARSCoV_1200_15_L EFT	TTTTAAGGAATTACTTGTGTATG CTGCT	1	28	60.0 6	32.1 4	1454 0

SARSCoV_1200_15_RI GHT	ACACACAACAGCATCGTCAGAG	1	22	61.12	50	1573 5
SARSCoV_1200_17_LE FT	TCAAGCTTTTTGCAGCAGAAAC G	1	23	61.28	43.4 8	1662 4
SARSCoV_1200_17_RI GHT	CCAAGCAGGGTTACGTGTAAGG	1	22	61.19	54.5 5	1775 4
SARSCoV_1200_19_L EFT	GGCACATGGCTTTGAGTTGACA	1	22	61.91	50	1859 6
SARSCoV_1200_19_RI GHT	CCTGTTGTCCATCAAAGTGTCC C	1	23	61.62	52.17	1967 8
SARSCoV_1200_21_L EFT	TCTGTAGTTTCTAAGGTTGTCAA AGTGA	1	28	60.5 8	35.7 1	2055 3
SARSCoV_1200_21_RI GHT	GCAGGGGGTAATTGAGTTCTGG	1	22	60.9 5	54.5 5	2164 2
SARSCoV_1200_23_L EFT	ACTTTAGAGTCCAACCAACAGA ATCT	1	26	60.18	38.4 6	22511
SARSCoV_1200_23_R IGHT	TGACTAGCTACACTACGTGCCC	1	22	61.52	54.5 5	2363 1
SARSCoV_1200_25_L EFT	TGCTGCTACTAAAATGTCAGAG TGT	1	25	60.51	40	2463 3
SARSCoV_1200_25_R IGHT	CATTTCCAGCAAAGCCAAAGCC	1	22	61.45	50	2579 0
SARSCoV_1200_27_L EFT	TGGATCACCGGTGGAATTGCTA	1	22	61.75	50	2674 4
SARSCoV_1200_27_RI GHT	TGTTCGTTTAGGCGTGACAAGT	1	22	60.74	45.4 5	2789 4
SARSCoV_1200_29_L EFT	TGAGGGAGCCTTGAATACACCA	1	22	61.1	50	2867 7
SARSCoV_1200_29_R IGHT	TAGGCAGCTCTCCCTAGCATTG	1	22	61.61	54.5 5	2979 0

Primers for Pool 1

Primer Name	Sequence	P o ol	Len gth	Tm	GC %	Start
SARSCoV_1200_2_LE FT	CCATAATCAAGACTATTCAACCA AGGGT	2	28	61.27	39.2 9	1100
SARSCoV_1200_2_RI GHT	ACAGGTGACAATTTGTCCACCG	2	22	61.33	50	2266
SARSCoV_1200_4_LE FT	GGAATTTGGTGCCACTTCTGCT	2	22	61.66	50	3144
SARSCoV_1200_4_RI GHT	CCTGACCCGGGTAAGTGGTTAT	2	22	61.49	54.5 5	4262
SARSCoV_1200_6_LE FT	ACTTCTATTAAATGGGCAGATAA CAACTG	2	29	60.18	34. 48	5257

SARSCoV_1200_6_RI GHT	GATTATCCATTCCCTGCGCGTC	2	22	61.75	54.5 5	6380
SARSCoV_1200_8_LE FT	CAATCATGCAATTGTTTTTCAGC TATTTTG	2	30	60.3 9	30	7298
SARSCoV_1200_8_RI GHT	TGACTTTTTGCTACCTGCGCAT	2	22	61.39	45.4 5	8385
SARSCoV_1200_10_LE FT	TTTACCAGGAGTTTTCTGTGGTG T	2	24	60.3 2	41.6 7	9303
SARSCoV_1200_10_RI GHT	TGGGCCTCATAGCACATTGGTA	2	22	61.5	50	1045 1
SARSCoV_1200_12_LE FT	ATGGTGCTAGGAGAGTGTGGAC	2	22	61.48	54.5 5	11372
SARSCoV_1200_12_RI GHT	GGATTTCCCACAATGCTGATGC	2	22	60.4 8	50	1256 0
SARSCoV_1200_14_LE FT	ACAGGCACTAGTACTGATGTCGT	2	23	61.12	47.8 3	1350 9
SARSCoV_1200_14_RI GHT	GTGCAGCTACTGAAAAGCACGT	2	22	61.94	50	1464 1
SARSCoV_1200_16_LE FT	ACAACACAGACTTTATGAGTGTC TCT	2	26	60.18	38.4 6	1560 8
SARSCoV_1200_16_RI GHT	CTCTGTCAGACAGCACTTCACG	2	22	61.17	54.5 5	1672 0
SARSCoV_1200_18_LE FT	GCACATAAAGACAAATCAGCTCA ATGC	2	27	62.0 3	40.7 4	1762 2
SARSCoV_1200_18_RI GHT	TGTCTGAAGCAGTGGAAAAGCA	2	22	60.6 8	45.4 5	1870 6
SARSCoV_1200_20_L EFT	ACAATTTGATACTTATAACCTCT GGAACAC	2	30	60.15	33.3 3	1957 4
SARSCoV_1200_20_RI GHT	GATTAGGCATAGCAACACCCGG	2	22	61.39	54.5 5	2069 8
SARSCoV_1200_22_L EFT	GTGATGTTCTTGTTAACAACTAA ACGAACA	2	30	61.44	33.3 3	2153 2
SARSCoV_1200_22_RI GHT	AACAGATGCAAATCTGGTGGCG	2	22	62.0 3	50	2261 2
SARSCoV_1200_24_L EFT	GCTGAACATGTCAACAACTCATA TGA	2	26	60.13	38.4 6	2351 8
SARSCoV_1200_24_RI GHT	ATGAGGTGCTGACTGAGGGAAG	2	22	61.74	54.5 5	2473 6
SARSCoV_1200_26_L EFT	GCCTTGAAGCCCCTTTTCTCTA	2	22	60.2 9	50	2569 0
SARSCoV_1200_26_RI GHT	AATGACCACATGGAACGCGTAC	2	22	61.5	50	2685 7
SARSCoV_1200_28_L EFT	TTTGTGCTTTTTAGCCTTTCTGC	2	24	60.14	37.5	2778 4
SARSCoV_1200_28_RI GHT	GTTTGGCCTTGTTGTTGTTGGC	2	22	61.82	50	2900 7

Primers for Pool 2

8 Generate primer pool stocks by adding $\Delta 5 \mu L$ of each odd region primer to a

▲ 1.5 mL Eppendorf labelled "Pool 1 (100µM)" and each even region primer to a

 \blacksquare 1.5 mL Eppendorf labelled "Pool 2 (100 μ M)". The pool is also given in the link

above. These are your 100 μ M stocks of each primer pool.

Note

Primers should be diluted and pooled in the **mastermix cabinet** which should be cleaned with decontamination wipes and UV sterilised before and after use.

9 Dilute this primer pool 1:10 in molecular grade water, to generate 10µM primer stocks. It is recommend that multiple aliquots of each primer pool are made to in case of degradation or contamination.

Note

Primers need to be used at a final concentration of 0.015μ M per primer. In this case (1200 nt amplicons), pool 1 has 30 primers and pool 2 has 28 primers, so the requirement is 1.13μ L for primer pool 1 and 1.05μ L for primer pool 2 (10uM) per 25μ L reaction. However, as these values are relatively close, we round up and down to 1.1ul for both pools, so the pools can be made in a similar fashion. For other schemes, adjust the volume added appropriately.

Multiplex PCR

10 In the mastermix hood set up the multiplex PCR reactions as follows in 0.2mL 8-strip PCR tubes:

Component	Pool 1 Pool 2	
5X Q5 Reaction Buffer	Δ 5 μL	👗 5 μL
10 mM dNTPs	Δ 0.5 μL	Δ 0.5 μL
Q5 Hot Start DNA Polymerase	Δ 0.25 μL	👗 0.25 μL

Primer Pool 1 or 2 (10µM)	-		
	4 1.1	μL 🕹 1.1	ΓμL
Nuclease-free water	프 1	5.9 μL	15.9 μL
Total	<u>д</u>	22.5 μL	22.5 μL
Note			
A PCR mastermix for each aliquoted into PCR strip tul the mastermix cabinet.			
In the extraction and samp mix well by pipetting.	le addition cab	inet add 🖁 2.5	μL cDNA to each tube and
Note			
Note			
The extraction and sample decontamination wipes an			
	d UV sterilised be	efore and after use	
decontamination wipes an	d UV sterilised be	efore and after use	
decontamination wipes an Pulse centrifuge the tubes	d UV sterilised be to collect the co m on the therm	efore and after use	
decontamination wipes an Pulse centrifuge the tubes to Set-up the following progra	d UV sterilised be to collect the co m on the therm	efore and after use ontents at the bott al cycler:	om of the tube.
decontamination wipes an Pulse centrifuge the tubes to Set-up the following progra Step	d UV sterilised be to collect the co m on the therm Temperature	efore and after use ontents at the bott al cycler: Time	om of the tube.
decontamination wipes and Pulse centrifuge the tubes of Set-up the following progra Step Heat Activation	d UV sterilised be to collect the co m on the therm Temperature § 98 °C	efore and after use ontents at the bott al cycler: Time (© 00:00:30	om of the tube. Cycles
decontamination wipes and Pulse centrifuge the tubes of Set-up the following progra Step Heat Activation Denaturation	d UV sterilised be to collect the co m on the therm Temperature § 98 °C § 98 °C	efore and after use ontents at the bott al cycler: Time (© 00:00:30 (© 00:00:15	om of the tube. Cycles 1 25-35
decontamination wipes and Pulse centrifuge the tubes of Set-up the following progra Step Heat Activation Denaturation Annealing and Extension	d UV sterilised be to collect the co m on the therm Temperature § 98 °C § 98 °C § 65 °C	efore and after use ontents at the bott al cycler: Time (© 00:00:30 (© 00:00:15 (© 00:05:00	om of the tube. Cycles 1 25-35 25-35

Expected result Final concentrations of PCR products can range from 20-150ng/ul. Pooling and PCR quantification 14 Label a Δ 1.5 mL Eppendorf tube for each sample and combine the two pools the PCR reaction as follows: Component Volume Pool 1 PCR reaction 🗸 25 uL Pool 2 PCR reaction 🗕 25 μL Total 🗕 50 μL Note At this stage, care should be taken with amplified PCR products. Only open tubes in a designated post-PCR workspace with equipment that is separate from areas where primers and mastermixes are handled. After combining the two pools of amplified DNA, the PCR products can be used for Oxford Nanopore Sequencing, using the RAPID barcode kit RBK004, as described in this protocol (below, Steps 15 onward). Alternatively, these amplicons can be used for Oxford Nanopore Sequencing, following Josh Quick's ligation based protocol (CoV-2019 sequencing protocol v2, dx.doi.org/10.17504/protocols.io.bdp7i5rn, at step 15) using the SQK-LSK109 kit. Alternatively, these amplicons can also be used for Illumina sequencing, such as found here: x.doi.org/10.17504/protocols.io.betejeje We have found that performing an Ampure XP bead clean up at this stage does not improve performance. Therefore, do not clean up the PCR reaction at this step. 14.1 Quantify DNA using a Qubit or other method. Quantification using Nanodrop is not recommended.

Protocol NAME DNA quantification using the Qubit f	
DNA quantification using the Qubit f	
CREATED BY	uorometer
Nikki Freed PI	REVIEW

14.1.1 Prepare a mastermix of Qubit[™] working solution for the required number of samples and standards. The Qubit dsDNA kit requires 2 standards for calibration (see note below).

rei sample.	
Qubit® dsDNA HS Reagent	Δ 1 μL
Qubit [®] dsDNA HS Buffer	Δ 199 μL

Note

Dor complex

If you have already performed a calibration on the Qubit machine for the selected assay you can use the previous calibration stored on the machine. We recommend performing a new calibration for every sample batch but a same-day calibration would be fine to use for multiple batches.

To avoid any cross-contamination, we recommend that you remove the total amount of working solution required for your samples and standards from the working solution bottle and then add the required volume to the appropriate tubes instead of pipetting directly from the bottle to each tube.

14.1.2 Label the tube lids. Do not label the side of the tube as this could interfere with the sample reading.

Note
Use only thin-wall, clear, 0.5mL PCR tubes. Acceptable tubes include Qubit™ assay tubes (Cat. No. Q32856)

- 14.1.3 Aliquot Qubit[™] working solution to each tube:
 - standard tubes requires 190µL of Qubit[™] working solution
 - sample tubes require anywhere from 180–199µL (depending how much sample you wish to add).

The final volume in each tube must be 200µL once sample/standard has been added.

- 14.1.4 Add 10µL of standard to the appropriate tube.
- 14.1.5 Add 1–20 μ L of each user sample to the appropriate tube.

Note

If you are adding $1-2\mu L$ of sample, use a P-2 pipette for best results.

- 14.1.6 Mix each tube vigorously by vortexing for 3–5 seconds.
- 14.1.7 Allow all tubes to incubate at room temperature for 2 minutes, then proceed to "Read standards and samples".
- 14.1.8 On the Home screen of the Qubit[™] 3 Fluorometer, press DNA, then select 1X dsDNA HS as the assay type. The Read standards screen is displayed. Press Read Standards to proceed.

Note

If you have already performed a calibration for the selected assay, the instrument prompts you to choose between reading new standards and running samples using the previous calibration. **If you want to use the previous calibration, skip to step 12**. Otherwise, continue with step 9.

14.1.9 Insert the tube containing Standard #1 into the sample chamber, close the lid, then press Read standard. When the reading is complete (~3 seconds), remove Standard #1.

- 14.1.10 Insert the tube containing Standard #2 into the sample chamber, close the lid, then press Read standard. When the reading is complete, remove Standard #2.
- 14.1.11 The instrument displays the results on the Read standard screen. For information on interpreting the calibration results, refer to the Qubit[™] Fluorometer User Guide, available for download at thermofisher.com/qubit.
- 14.1.12 Press Run samples.
- 14.1.13 On the assay screen, select the sample volume and units:
 - Press the + or buttons on the wheel, or anywhere on the wheel itself, to select the sample volume added to the assay tube (from 1–20µL).
 - From the unit dropdown menu, select the units for the output sample concentration (in this case choose ng/μL).
- 14.1.14 Insert a sample tube into the sample chamber, close the lid, then press Read tube. When the reading is complete (~3 seconds), remove the sample tube.
- 14.1.15 The top value (in large font) is the concentration of the original sample and the bottom value is the dilution concentration. For information on interpreting the sample results, refer to the Qubit[™] Fluorometer User Guide.
- 14.1.16 Repeat step 14 until all samples have been read.
- 14.1.17 Carefully **record all results** and store run file from the Qubit on a memory stick.
- 14.1.18 All negative controls should ideally be 'too low' to read on the Qubit machine, but MUST be < 1ng per ul. If your negative controls >1ng per ul, considerable contamination has occurred and you must redo previous steps.

Normalisation

- 15 Label a 40.2 mL PCR tube for each sample.
- 15.1 Adjust the amount of DNA in the tube to be <u>Δ 100 ng</u> total per sample in <u>Δ 7.5 μL</u> molecular grade water. For example if your PCR reaction is at 100ng/ul, add 1ul of the PCR reaction to 6.5ul of molecular grade water. Input to the Rapid Barcoding kit will vary depending on the amplicon length but we have determined 50-200 ng works for efficient

barcoding of this amplicon length. Use 7.5ul of the negative control, even if there is no detectable DNA in the PCR reaction.

Rapid barocoding using the SQK RBK004

16 Mulitple samples can be run on the same flow cell by barcoding. Up to 12 samples at a time can be run. Amplicons from each sample will be individually barcoded in the following steps. These follow the RBK004 protocol from Oxford Nanopore. It is highly recommended to use their protocol for the following steps.

SQK-RBK004 Rapid Barcoding Kit **Oxford Nanopore Technologies Catalog** #SQK-RBK004

16.1 Add $\boxed{_}$ 7.5 μ L of each diluted PCR reaction from step 15 to the labeled PCR tube. Set up the following reaction for each sample:

Component	Volume	
DNA amplicons from step 15 (100ng total)	Δ 7.5 μL	
Fragmentation Mix RB01-12 (one for each sample, included in kit)	Δ 2.5 μL	
Total		

🗕 10 μL

- 16.2 Mix gently by flicking the tube, and spin down.
- 16.3 Incubate the reaction in a PCR machine:
 - \$30 °C
 for
 \$00:01:00

 \$80 °C
 for
 \$00:01:00

 \$4 °C
 for
 \$00:00:30
- 16.4 Pool all barcoded samples, noting the total volume.

17 Ampure Bead Cleanup. Use a 1:1 ratio of sample to beads.

5m

5m

15m

	Protocol
	Amplicon clean-up using SPRI beads
	CREATED BY Nikki Freed PREVIEW
17.1	Vortex SPRI beads thoroughly to ensure they are well resuspended, the solution should be a homogenous brown colour.
	X Agencourt AMPure XP Beckman Coulter Catalog #A63880
17.2	Add an equal volume (1:1) of SPRI beads to the sample tube and mix gently by either flicking or pipetting. For example add $\boxed{4}$ 50 µL room temperature SPRI beads to a
17.3	$\boxed{4}$ 50 µL reaction. Pulse centrifuge to collect all liquid at the bottom of the tube.
17.5	
17.4	Incubate for 👀 00:05:00 at room temperature.
17.5	Place on magnetic rack and incubate for $\bigcirc 00:02:00$ or until the beads have pelleted and the supernatant is completely clear.
17.6	Carefully remove and discard the supernatant, being careful not to touch the bead pellet
17.7	Add $\underline{\square}$ 200 μ L of freshly prepared room-temperature [M] 80 % volume ethanol to the pellet.

- 17.8 Keeping the magnetic rack on the benchtop, rotate the bead-containing tube by 180°.Wait for the beads to migrate towards the magnet and re-form a pellet. Remove the ethanol using a pipette and discard.
- 17.9 and repeat ethanol wash.
- 17.10 Pulse centrifuge to collect all liquid at the bottom of the tube and carefully remove as much residual ethanol as possible using a P10 pipette.
- ^{17.11} With the tube lid open incubate for 00:01:00 or until the pellet loses it's shine (if the pellet dries completely it will crack and become difficult to resuspend).
- 17.12 Remove the tube from the magnetic rack. Resuspend pellet in ▲ 10 µL molecular grade water or Elution buffer, mix gently by flicking and incubate for 00:02:00.
 W Elution Buffer (EB) Qiagen Catalog #19086
- 17.13 Place on magnet and transfer sample to a clean 1.5mL Eppendorf tube ensuring no beads are transferred into this tube.
- 17.14 Quantify Δ 1 μL product using the Quantus Fluorometer using the ONE dsDNA assay.
 QuantiFluor(R) ONE dsDNA System, 100rxn Promega Catalog #E4871

Equipment	
Quantus	NAME
Fluorometer	TYPE
Promega	BRAND
E6150	SKU
https://www.promega.co.uk/products/microplate-readers-fluorometers- luminometers/fluorometers/quantus-fluorometer	LINK

Add $\underline{\square} 1 \mu \underline{\square}$ of RAP (from the RBK004 kit) to $\underline{\square} 10 \mu \underline{\square}$ cleaned, barcoded DNA from step 17. Mix gently by flicking the tube, and spin down.

19 Incubate the reaction for 🚫 00:05:00 at room temperature.

20 The prepared library is used for loading into the MinION flow cell according to Oxford Nanopore Rapid Barcoding (RBK004) protocol. Store the library on ice until ready to load.

MinION sequencing

21 Start the sequencing run using MinKNOW.

1m

5m

10m

	Protoco			
		NAME		
		Starting a MinION sec	luencing run using Minl	KNOW
	CREATED			-\
	Nikki Fr	∋ed	PREVIE	
21.1	If required detected.	d plug the MinION into the o	computer and wait for the I	VinION and
21.2	Choose fl	ow cell 'FLO-MIN106' from	the drop-down menu.	

- 21.3 Then select the flowcell so a tick appears.
- 21.4 Click the 'New Experiment' button in the bottom left of the screen.
- 21.5 On the New experiment popup screen, select the running parameters for your experiment from the individual tabs:

Experiment: Name the run in the experiment field, leave the sample field blank.

Kit: Selection: Select RBK004

Run Options: Set the run length to 6 hours (you can stop the run once sufficient data has been collected as determined using RAMPART).

Basecalling: Select 'fast basecalling'.

Output: The number of files that MinKNOW will write to a single folder. By default this is set to 4000 but can be reduced to make RAMPART update more frequently.

Click 'Start run'.

.

- 21.6 Monitor the progress of the run using the MinKNOW interface.
- 22 Depending on the variation in coverage of each amplicon, generally, you will need approx 10,000 to 20,000 reads or 10-20Mb **per sample** to confidently assemble and call variants. This can typically be achieved on a minION flow cell in under two hours when runnning 12 samples. Shorter, if running fewer samples.
- 23 The primer scheme .bed and .tsv files necessary for the ARTIC variant calling pipeline are <u>here</u>