Oct 16, 2020

NEB Instant Sticky-end Ligase Master Mix-CHEM 584

Forked from Cloning with NEB Instant Sticky-end Ligase Master Mix (M0370)

DOI

dx.doi.org/10.17504/protocols.io.bninmcde

Ken Christensen¹, Harold Bien²

¹Brigham Young University; ²Stony Brook University

Ken Christensen Brigham Young University





DOI: dx.doi.org/10.17504/protocols.io.bninmcde

External link: <u>https://www.neb.com/protocols/2012/08/27/protocol-transfer-master-mix-to-ice-prior-to-reaction-set-up-</u> <u>mix-tube-by-finger-flicking-before-u</u>

Protocol Citation: Ken Christensen, Harold Bien 2020. NEB Instant Sticky-end Ligase Master Mix-CHEM 584. protocols.io https://dx.doi.org/10.17504/protocols.io.bninmcde

License: This is an open access protocol distributed under the terms of the **<u>Creative Commons Attribution License</u>**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development We are still developing and optimizing this protocol

Created: October 16, 2020

Last Modified: October 16, 2020

Protocol Integer ID: 43310

Abstract

Taken directly from NEB's website, this is the protocol for using their optimized sticky-end ligation and transformation. For more information, see https://www.neb.com/products/m0370-instant-sticky-end-ligase-master-mix.

Instant Sticky-end Ligase Master Mix is a ready-to-use 2X solution of T4 DNA ligase and a proprietary ligation enhancer in an optimized reaction buffer. Specifically formulated to rapidly ligate cohesive-end (2–4 bp) substrates and improve transformation, the mix simplifies reaction set-up and provides an optimized ratio of enzyme and buffer components. No thawing of the master mix is required as it maintains a liquid state during storage at -20°C* and no incubation time is necessary to achieve ligation efficiencies sufficient for successful sub-cloning of sticky-end substrates. Just add the master mix to DNA with compatible ends, mix and transform; thereby reducing the valuable time needed for routine ligations. Ligations for subcloning can be carried out in small volumes with low concentrations, allowing users to conserve precious DNA samples, and be used directly to transform many strains of chemically competent *E. coli* without dilution.

* Freezers vary in their actual internal temperature. Our testing demonstrates that the master mix is liquid at -20°C. Freeze-thaw testing at -70°C has confirmed that the performance after 20 freeze/thaw cycles is close to the original mix.


Product Source

Purified from an *E. coli* strain containing a recombinant gene encoding T4 DNA Ligase.

Reaction Volume Definition

1X Instant Sticky-end Ligase Master Mix with DNA substrates in a 10 μ l reaction volume. A 10 μ l reaction contains 1800 cohesive end units of T4 DNA Ligase.

Attachments



Guidelines

Once, completed, ligation products can be used immediately for transformation or stored at -20C until ready to use. In-house testing has demonstrated that maximal transformation efficiency is achieved using between 20–100 ng of vector (sticky) and a corresponding 3-fold molar excess of the insert to be ligated into the vector.

Chemically competent strains of *E. coli* (commercially available or prepared by user) can be transformed by ligation products prepared using the Instant Sticky-end Ligase Master Mix. Electrocompetent cells are not compatible. Users of competent cells from other vendors may need to dilute ligation reactions 4-fold, prior to transformation, in order to achieve maximum transformation efficiency. Not all cells from other vendors will benefit from this additional step. The following protocol is recommended by NEB. Other protocols can be used but the volume of ligation reaction used should not exceed 5 µl reaction per 50 µl cells.

Transformation efficiencies around 2×10^6 cfu/µg are typically achieved for recombinant cohesive end substrates (vector + insert), using cells with a 7×10^8 calculated efficiency with uncut DNA. This corresponds to several hundred colonies on a plate when 100 µl of a 1 ml outgrowth is plated at a 1:5 dilution. As with all ligation and transformation protocols, many factors affect the calculated transformation efficiency, including purity and integrity of DNA ends, competence of the cells being transformed, media choices, incubation temperatures and times and biological effects (intact ORF in high-copy vector, toxic genes, etc.).

Materials

MATERIALS

X Instant Sticky-end Ligase Master Mix - 50 rxns New England Biolabs Catalog #M0370S STEP MATERIALS

🔀 Instant Sticky-end Ligase Master Mix - 50 rxns New England Biolabs Catalog #M0370S

🔀 Instant Sticky-end Ligase Master Mix - 50 rxns New England Biolabs Catalog #M0370S

Protocol materials

X Instant Sticky-end Ligase Master Mix - 50 rxns New England Biolabs Catalog #M0370S

🔀 Instant Sticky-end Ligase Master Mix - 50 rxns New England Biolabs Catalog #M0370S

X Instant Sticky-end Ligase Master Mix - 50 rxns New England Biolabs Catalog #M0370S

🔀 Instant Sticky-end Ligase Master Mix - 50 rxns New England Biolabs Catalog #M0370S

🔀 Instant Sticky-end Ligase Master Mix - 50 rxns New England Biolabs Catalog #M0370S

Safety warnings

1 Do not heat inactivate.

Heat inactivation dramatically reduces transformation efficiency.

Preparation

1 Place master mix tube on ice and flick a few times to mix.

Instant Sticky-end Ligase Master Mix - 50 rxns New England Biolabs Catalog #M0370S

Ligation

2 Combine 20–100 ng of vector with a 3-fold molar excess of insert and q.s. to 5 μ l with dH₂O.

🗕 5 μL

3 Add 5 μl of Instant Sticky-end Ligase Master Mix, mix thoroughly by pipetting up and down 7-10 times, and place on ice. The sample is now ready to be used for transformation.

🗕 5 μL

Instant Sticky-end Ligase Master Mix - 50 rxns New England Biolabs Catalog #M0370S

Transformation into Mix & Go Competent Cells

- 4 Thaw a 50 μl aliquot of Mix & Go competent cells on ice.
- 5 Add 5 μl of the ligation reaction to the cells and mix by finger-flicking. Do not vortex the tube.

🗛 5 μL

6 Incubate the tube on ice for 5 minutes. Do not mix.



7 For anitbiotics other than ampicillin (e.g., kanamycin), add 200 μl of warm SOC to the tube. For ampicillin selective plasmids, no SOC addition or outgrowth is needed so you can proceed to STEP 9 immediately.

🗕 200 μL

Incubate for one hour at 37°C with rotation or shaking (200–250 rpm).
01:00:00

5m

5m

9 Spread the entire tube of cells (e.g., 50 or 250 μl) onto an appropriate antibiotic selection plate and incubate overnight at 37°C.