

Jul 22, 2019

Version 1

# NEBNext® Ultra™ II End Repair/dA-Tailing Module (NEB #E7546) V.1



DOI

[dx.doi.org/10.17504/protocols.io.4nngvde](https://dx.doi.org/10.17504/protocols.io.4nngvde)

Menna Teffera<sup>1</sup>, New England Biolabs<sup>1</sup>

<sup>1</sup>New England Biolabs

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: [info@neb.com](mailto:info@neb.com)



Isabel Gautreau

New England Biolabs

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4nngvde>

External link: <https://www.neb.com/-/media/catalog/datacards-or-manuals/manuale7546.pdf?rev=221a72d9f6624cfda4fc2fe760691f02>

**Protocol Citation:** Menna Teffera, New England Biolabs 2019. NEBNext® Ultra™ II End Repair/dA-Tailing Module (NEB #E7546).  
**protocols.io** <https://dx.doi.org/10.17504/protocols.io.4nngvde>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** June 24, 2019

**Last Modified:** June 22, 2023

**Protocol Integer ID:** 25006

**Keywords:** nebnex ultra ii dna library prep, nebnex ultra ii ligation module, ultra ii dna library prep kit for illumina, ultra ii dna library construction workflow, ultra ii dna library prep kit, use with the nebnex, ultra ii ligation module, use with nebnex singleplex, nebnex singleplex, nebnex, multiplex oligos for illumina, repaired dna, ultra ii end repair, fragmented dna, multiplex oligo, nanopore, neb, dna, ultra ii q5 master mix

## Abstract

The NEBNext Ultra II End Repair/dA-Tailing Module is optimized to convert 500 pg-1 µg of fragmented DNA to repaired DNA having 5′ phosphorylated, 3′ dA-tailed ends.

This module is part of the Ultra™ II workflow, and is optimized for use with the NEBNext™ Ultra II Ligation Module (NEB #E7595), for Illumina®-compatible library construction.

This module is also compatible with some Oxford Nanopore MinION™ workflows.

This module is designed for use with NEBNext Singleplex or Multiplex Oligos for Illumina (NEB #E7350, #E7335, #E7500, #E7600 or #E7535), NEBNext Ultra II Ligation Module (NEB #E7595), and NEBNext Ultra II Q5 Master Mix (NEB #M0544).

Kits that include reagents for every step in the Ultra II DNA library construction workflow are also available (NEBNext Ultra II DNA Library Prep Kit for Illumina (NEB [#E7645](#)) and NEBNext Ultra II DNA Library Prep with Sample Purification Beads (NEB [#E7103](#)).

## Guidelines

**Safe Stop Point:** This is a point where you can safely stop the protocol and store the samples prior to proceeding to the next step in the protocol.


**Caution:** Signifies a step in the protocol that has two paths leading to the same point.

**Color:** A color listed before or after a reagent name indicates the cap color of the reagent to be added.



## Materials

### MATERIALS

 NEBNext Ultra II End Prep Reaction Buffer **New England Biolabs Catalog #E7647**

 NEBNext Ultra II End Prep Enzyme Mix **New England Biolabs Catalog #E7646**

## Troubleshooting

### Before start

**Starting Material:** 500 pg–1 µg fragmented DNA. We recommend that DNA be sheared in 1X TE. If the DNA volume post shearing is less than 50 µl, add 1X TE to a final volume of 50 µl. Alternatively, 10 mM Tris-HCl, pH 8.0 or 0.1X TE can be used.

## NEBNext End Prep

- 1 Mix the following contents in a sterile nuclease-free tube:

	<b>Component</b>	<b>Volume</b>
	(green) NEB Next Ultra II End Prep Enzyme Mix	3 $\mu$ l
	(green) NEB Next Ultra II End Prep Reaction Buffer	7 $\mu$ l
	Frag men ted DNA	50 $\mu$ l
	<b>Total Volume</b>	<b>60 <math>\mu</math>l</b>

- 2 Set a 100  $\mu$ l or 200  $\mu$ l pipette to 50  $\mu$ l and then gently pipette the entire volume up and down at least 10 times to mix thoroughly. Perform a quick spin to collect all liquid from the sides of the tube.

### Note

Note: It is important to mix well. The presence of a small amount of bubbles will not interfere with performance.



- 3 Place in a thermocycler, with the heated lid set to  $\geq 75^{\circ}\text{C}$ , and run the following program:

1h

 00:30:00 at  20 °C

 00:30:00 at  65 °C

Hold at  4 °C

#### Note

Safe Stop Point: If necessary, samples can be stored at  $-20^{\circ}\text{C}$ ; however, a slight loss in yield (~20%) may be observed. We recommend continuing with adaptor ligation before stopping.

- 4 Proceed directly to NEBNext Ultra II Ligation Module NEB #E7595.