



Jul 22, 2019

Version 1

# NEBNext Ultra II Ligation Module (NEB # E7595) for NEBNext Ultra II End Repair/dA Tailing Module (NEB #E7546) V.1



DOI

[dx.doi.org/10.17504/protocols.io.4nqgvdw](https://dx.doi.org/10.17504/protocols.io.4nqgvdw)

New England Biolabs<sup>1</sup>, Menna Teffera<sup>1</sup>

<sup>1</sup>New England Biolabs



Isabel Gautreau

New England Biolabs

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4nqgvdw>

External link: <https://www.neb.com/-/media/catalog/datacards-or-manuals/manuale7595.pdf?rev=6d5ae698cf394a27864e788e6056e8d8>

**Protocol Citation:** New England Biolabs, Menna Teffera 2019. NEBNext Ultra II Ligation Module (NEB # E7595) for NEBNext Ultra II End Repair/dA Tailing Module (NEB #E7546). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4nqgvdw>



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** June 24, 2019

**Last Modified:** June 22, 2023

**Protocol Integer ID:** 25008

**Keywords:** nebnext ultra ii ligation module, use with the nebnext, ultra ii ligation module, nebnext, ultra ii fs dna module, ultra ii end repair, neb, e7546, tailing module, e7595, e7810, dna, module

## Abstract

This module is part of the Ultra™ II workflow, and is optimized for use with the NEBNext® Ultra II End Repair/dA-Tailing Module (NEB #**E7546**), for Illumina®-compatible library construction.

The NEBNext Ultra II Ligation Module is optimized for use with the NEBNext Ultra II End Repair/dA-Tailing Module (NEB #**E7546**) or the NEBNext Ultra II FS DNA Module (NEB #**E7810**).

## Guidelines

**Safe Stop Point:** This is a point where you can safely stop the protocol and store the samples prior to proceeding to the next step in the protocol.

**Caution:** Signifies a step in the protocol that has two paths leading to the same point.

**Color:** A color listed before or after a reagent name indicates the cap color of the reagent to be added.

## Adaptor Dilution Guidelines

The appropriate adaptor dilution for your sample input and type may need to be optimized experimentally. The dilutions provided here are a general starting point.

**Table 1.1: Adaptor Dilution**

	Input	Adaptor Dilution (Volume of adaptor : Total volume)	Working Adaptor Concentration
	1 µg–101 ng	No Dilution	15 µM
	100 ng–5 ng	10-Fold (1:10)	1.5 µM
	less than 5 ng	25-Fold (1:25)	0.6 µM



## Materials

### MATERIALS

- ✕ Ligation Enhancer **New England Biolabs Catalog #E7374** in Kits E7370 or E7445
- ✕ NEBNext Adaptor for Illumina **New England Biolabs Catalog #E7337** in Kits E7335, E7500, E771
- ✕ NEBNext Ultra II Ligation Master Mix **New England Biolabs Catalog #E7648**
- ✕ USER Enzyme (Multiplex Oligos for Illumina) **New England Biolabs Catalog #E7338**

### STEP MATERIALS

- ✕ USER Enzyme (Multiplex Oligos for Illumina) **New England Biolabs Catalog #E7338**
- ✕ NEBNext Ultra II Ligation Master Mix **New England Biolabs Catalog #E7648**
- ✕ NEBNext Ligation Enhancer **New England Biolabs Catalog #E7374**
- ✕ NEBNext Adaptor for Illumina **New England Biolabs**

Materials that may be needed that are not included in this kit:

Tris-HCL Buffer

10 mM NaCl

## Protocol materials

- ✕ USER Enzyme (Multiplex Oligos for Illumina) **New England Biolabs Catalog #E7338**
- ✕ NEBNext Ultra II Ligation Master Mix **New England Biolabs Catalog #E7648**
- ✕ NEBNext Ligation Enhancer **New England Biolabs Catalog #E7374**
- ✕ NEBNext Adaptor for Illumina **New England Biolabs**
- ✕ Ligation Enhancer **New England Biolabs Catalog #E7374** in Kits E7370 or E7445
- ✕ NEBNext Adaptor for Illumina **New England Biolabs Catalog #E7337** in Kits E7335, E7500, E771
- ✕ NEBNext Ultra II Ligation Master Mix **New England Biolabs Catalog #E7648**
- ✕ USER Enzyme (Multiplex Oligos for Illumina) **New England Biolabs Catalog #E7338**
- ✕ NEBNext Ultra II Ligation Master Mix **New England Biolabs Catalog #E7648**
- ✕ NEBNext Ligation Enhancer **New England Biolabs Catalog #E7374**
- ✕ NEBNext Adaptor for Illumina **New England Biolabs**
- ✕ USER Enzyme (Multiplex Oligos for Illumina) **New England Biolabs Catalog #E7338**

## Troubleshooting



## Before start

### Starting Material:

500 pg–1 µg fragmented DNA that has been end repaired and dA-Tailed using the NEBNext End Repair/dA-Tailing Module (**NEB #E7546**).

#### Note

**Caution:** If DNA input is  $\leq 100$  ng, dilute the NEBNext Adaptor for Illumina in 10 mM Tris-HCl or 10 mM Tris-HCl with 10 mM NaCl as indicated in Table 1.1.

## Ligation/End Prep

- 1 Add the following products directly to the End Prep Reaction Mixture:

Component	Volume
End Prep Reaction Mixture	60 $\mu$ l
(red) NEB Next Ultra II Ligation Master Mix*	30 $\mu$ l
(red) NEB Next Ligation Enhancer	1 $\mu$ l
(red) NEB Next Adaptor for Illumina*	2.5 $\mu$ l
<b>Total Volume</b>	<b>93.5 <math>\mu</math>l</b>


\* Mix the Ultra II Ligation Master Mix by pipetting up and down several times prior to adding to the reaction.

\*\* The NEBNext adaptor is provided in NEBNext Singleplex (NEB #E7350) or Multiplex (NEB #E7335, #E7500, #E7710, #E7730, #E7600, #E7535, and #E6609) Oligos for Illumina.

 NEBNext Ultra II Ligation Master Mix **New England Biolabs Catalog #E7648**



 NEBNext Ligation Enhancer **New England Biolabs Catalog #E7374**

 NEBNext Adaptor for Illumina **New England Biolabs**



#### Note

Note: The Ligation Master Mix and Ligation Enhancer can be mixed ahead of time and is stable for at least 8 hours @ 4°C. We do not recommend premixing the Ligation Master Mix, Ligation Enhancer and adaptor prior to use in the Adaptor Ligation Step.

- 2 Set a 100 µl or 200 µl pipette to 80 µl and then pipette the entire volume up and down at least 10 times to mix thoroughly. Perform a quick spin to collect all liquid from the sides of the tube.


#### Note

**Caution:** The NEBNext Ultra II Ligation Master Mix is very viscous. Care should be taken to ensure adequate mixing of the ligation reaction, as incomplete mixing will result in reduced ligation efficiency. The presence of a small amount of bubbles will not interfere with performance.

- 3 Incubate at  20 °C for  00:15:00 in a thermocycler **with the heated lid off**.

15m

- 4 Add 3 µl of (red) USER™ Enzyme to the ligation mixture from Step 3.

 USER Enzyme (Multiplex Oligos for Illumina) **New England Biolabs Catalog #E7338**

#### Note

Note: Steps 4 and 5 are only required for use with NEBNext Adaptors. USER enzyme can be found in the NEBNext Singleplex (NEB #E7350) or Multiplex (NEB #E7335, #E7500, #E7710, #E7730, #E7600 and #E6609) Oligos for Illumina.



- 5 Mix well and incubate at 37 °C for 00:15:00 with the heated lid set to ≥ 47 °C .

15m

- 6 DNA is now ready for size selection or cleanup.

**Note**

Note: Please see NEB #E7645 manual for recommended size selection/cleanup and PCR amplification protocols.

**Note**

**Safe Stop Point:** Samples can be stored overnight at –20°C.