



Jul 23, 2019 Version 1

🌐 NEBNext Ultra II Ligation Module (NEB # E7595) for NEBNext Ultra II End Repair/dA Tailing Module (NEB #E7546) V.1



DOI

dx.doi.org/10.17504/protocols.io.4nqgvdw

New England Biolabs¹, Menna Teffera¹

¹New England Biolabs



Isabel Gautreau

New England Biolabs

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.4nqgvdw

External link: <https://www.neb.com/-/media/catalog/datacards-or-manuals/manuale7595.pdf?rev=6d5ae698cf394a27864e788e6056e8d8>

Protocol Citation: New England Biolabs, Menna Teffera 2019. NEBNext Ultra II Ligation Module (NEB # E7595) for NEBNext Ultra II End Repair/dA Tailing Module (NEB #E7546). [protocols.io https://dx.doi.org/10.17504/protocols.io.4nqgvdw](https://dx.doi.org/10.17504/protocols.io.4nqgvdw)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: June 25, 2019

Last Modified: June 23, 2023

Protocol Integer ID: 25008

Abstract

This module is part of the Ultra™ II workflow, and is optimized for use with the NEBNext® Ultra II End Repair/dA-Tailing Module (NEB #[E7546](#)), for Illumina®-compatible library construction.

The NEBNext Ultra II Ligation Module is optimized for use with the NEBNext Ultra II End Repair/dA-Tailing Module (NEB #[E7546](#)) or the NEBNext Ultra II FS DNA Module (NEB #[E7810](#)).

Guidelines

Safe Stop Point: This is a point where you can safely stop the protocol and store the samples prior to proceeding to the next step in the protocol.

Caution: Signifies a step in the protocol that has two paths leading to the same point.

Color: A color listed before or after a reagent name indicates the cap color of the reagent to be added.

Adaptor Dilution Guidelines




The appropriate adaptor dilution for your sample input and type may need to be optimized experimentally. The dilutions provided here are a general starting point.

Table 1.1: Adaptor Dilution





Input	Adaptor Dilution (Volume of adaptor or: Total volume)	Working Adaptor Concentration
1 µg–101 ng	No Dilution	15 µM
100 ng–5 ng	10-Fold (1:10)	1.5 µM
less than 5 ng	25-Fold (1:25)	0.6 µM

Materials

MATERIALS

-  Ligation Enhancer **New England Biolabs Catalog #E7374** in Kits E7370 or E7445
-  NEBNext Adaptor for Illumina **New England Biolabs Catalog #E7337** in Kits E7335, E7500, E771
-  NEBNext Ultra II Ligation Master Mix **New England Biolabs Catalog #E7648**
-  USER Enzyme (Multiplex Oligos for Illumina) **New England Biolabs Catalog #E7338**





STEP MATERIALS

-  USER Enzyme (Multiplex Oligos for Illumina) **New England Biolabs Catalog #E7338**
-  NEBNext Ultra II Ligation Master Mix **New England Biolabs Catalog #E7648**
-  NEBNext Ligation Enhancer **New England Biolabs Catalog #E7374**
-  NEBNext Adaptor for Illumina **New England Biolabs**

Materials that may be needed that are not included in this kit:

Tris-HCL Buffer
10 mM NaCl

Protocol materials

-  USER Enzyme (Multiplex Oligos for Illumina) **New England Biolabs Catalog #E7338**
-  NEBNext Ultra II Ligation Master Mix **New England Biolabs Catalog #E7648**
-  NEBNext Ligation Enhancer **New England Biolabs Catalog #E7374**
-  NEBNext Adaptor for Illumina **New England Biolabs**
-  Ligation Enhancer **New England Biolabs Catalog #E7374** in Kits E7370 or E7445
-  NEBNext Adaptor for Illumina **New England Biolabs Catalog #E7337** in Kits E7335, E7500, E771
-  NEBNext Ultra II Ligation Master Mix **New England Biolabs Catalog #E7648**
-  USER Enzyme (Multiplex Oligos for Illumina) **New England Biolabs Catalog #E7338**
-  NEBNext Ultra II Ligation Master Mix **New England Biolabs Catalog #E7648**
-  NEBNext Ligation Enhancer **New England Biolabs Catalog #E7374**
-  NEBNext Adaptor for Illumina **New England Biolabs**
-  USER Enzyme (Multiplex Oligos for Illumina) **New England Biolabs Catalog #E7338**

Before start

Starting Material:

500 pg–1 µg fragmented DNA that has been end repaired and dA-Tailed using the NEBNext End Repair/dA-Tailing Module (**NEB #E7546**).

Note

Caution: If DNA input is ≤ 100 ng, dilute the NEBNext Adaptor for Illumina in 10 mM Tris-HCl or 10 mM Tris-HCl with 10 mM NaCl as indicated in Table 1.1.

Ligation/End Prep

- 1 Add the following products directly to the End Prep Reaction Mixture:


Component	Volume
End Prep Reaction Mixture	60 µl
(red) NEBNext Ultra II Ligation Master Mix*	30 µl
(red) NEBNext Ligation Enhancer	1 µl
(red) NEBNext Adaptor for Illumina**	2.5 µl
Total Volume	93.5 µl

* Mix the Ultra II Ligation Master Mix by pipetting up and down several times prior to adding to the reaction.

** The NEBNext adaptor is provided in NEBNext Singleplex (NEB #E7350) or Multiplex (NEB #E7335, #E7500, #E7710, #E7730, #E7600, #E7535, and #E6609) Oligos for Illumina.

 NEBNext Ultra II Ligation Master Mix **New England Biolabs Catalog #E7648**

 NEBNext Ligation Enhancer **New England Biolabs Catalog #E7374**

 NEBNext Adaptor for Illumina **New England Biolabs**



**Note**

Note: The Ligation Master Mix and Ligation Enhancer can be mixed ahead of time and is stable for at least 8 hours @ 4°C. We do not recommend premixing the Ligation Master Mix, Ligation Enhancer and adaptor prior to use in the Adaptor Ligation Step.

- 2 Set a 100 µl or 200 µl pipette to 80 µl and then pipette the entire volume up and down at least 10 times to mix thoroughly. Perform a quick spin to collect all liquid from the sides of the tube.

Note

Caution: The NEBNext Ultra II Ligation Master Mix is very viscous. Care should be taken to ensure adequate mixing of the ligation reaction, as incomplete mixing will result in reduced ligation efficiency. The presence of a small amount of bubbles will not interfere with performance.

- 3 Incubate at  20 °C for  00:15:00 in a thermocycler **with the heated lid off**.

15m




- 4 Add 3 µl of (red) USER™ Enzyme to the ligation mixture from Step 3.



USER Enzyme (Multiplex Oligos for Illumina) **New England**
Biolabs Catalog #E7338

Note

Note: Steps 4 and 5 are only required for use with NEBNext Adaptors. USER enzyme can be found in the NEBNext Singleplex (NEB #E7350) or Multiplex (NEB #E7335, #E7500, #E7710, #E7730, #E7600 and #E6609) Oligos for Illumina.

- 5 Mix well and incubate at  37 °C for  00:15:00 with the heated lid set to ≥  47 °C .

15m

- 6 DNA is now ready for size selection or cleanup.



Note

Note: Please see NEB #E7645 manual for recommended size selection/cleanup and PCR amplification protocols.

Note

Safe Stop Point: Samples can be stored overnight at -20°C .