

Sep 30, 2019

Version 2

Nerve Conduction Velocity Tests V.2

DOI

dx.doi.org/10.17504/protocols.io.7tghnjw



Eva Feldman¹

¹University of Michigan - Ann Arbor

Diabetic Complications Consortium

Tech. support email: rmcindoe@augusta.edu



Lili Liang

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7tghnjw>

External link: <https://www.diacomp.org/shared/document.aspx?id=35&docType=Protocol>



Protocol Citation: Eva Feldman 2019. Nerve Conduction Velocity Tests. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.7tghnjw>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 30, 2019

Last Modified: September 30, 2019

Protocol Integer ID: 28232

Keywords: nerve conduction velocity (NCV), diabetic neuropathy, nerve conduction velocity tests summary, nerve conduction velocity, nerve study, dermal temperature probe, body temperature, electrode, diabetic complication

Abstract

Summary:

To confirm the presence of diabetic neuropathy, nerve conduction velocity (NCV) studies are performed. The animals are anesthetized with 30/2.5 mg/kg ketamine/xylazine to prevent discomfort. Body temperature is monitored with a dermal temperature probe and maintained at 32° C with a warming lamp during NCV. Body temperature is maintained at 37°C after NCV using a warming pad to ease animal stress from anesthetic. The nerve studies last less than 30 min per rat or mouse. The electrodes are cleaned with 70% alcohol between animals to maintain pathogen-free status.

Diabetic Complication:



Neuropathy



Materials

Reagents and Materials:

Equipment:

◆ **Nicolet VikingQuest Portable System with Nerve Conduction Studies VikingQuest software run on Windows NT**

◆ HP laser printer

◆ Nicolet 12mm .4mm diameter disposable platinum EEG subdermal needles

◆ Niclot disposable ground

◆ Heating lamp

◆ **Infrared thermometer**

◆ Heating pad

◆ Flexible tape measure

◆ 8"x 8" Styrofoam

◆ Ketaset 100mg/ml (3 parts) and Rompun 20 mg/ml (1 part) --- stock solution should be diluted 1/10 for mice and as is for rats.

Troubleshooting

Settings:

- 1
 - ◆ Motor tests
 - o Duration .02 ms
 - o Range 25 mA
 - o low frequency filter 1 Hz
 - o High frequency filter 10 kHz
 - o Sensitivity 1 mV
 - o Time 2 ms/div
 - ◆ Sensory test
 - o Duration .02 ms
 - o Range 25 mA
 - o low frequency filter 1 Hz
 - o High frequency filter 10 kHz
 - o Sensitivity 50 μ V
 - o Time 2 ms/div

Procedure:

- 2 To confirm the presence of diabetic neuropathy, nerve conduction velocity (NCV) studies are performed. The animals will be anesthetized with 30-100/2.5-10 mg/kg ketamine/xylazine to prevent discomfort. Body temperature is monitored with a dermal temperature probe and maintained at 37° C with a warming lamp during NCV. Skin temperature is maintained at 34° C during NCV using an infrared thermometer. Body temperature is maintained at 37°C after NCV using a warming pad to ease animal stress from anesthetic. The nerve studies will last less than 30 min per rat or mouse. The electrodes are cleaned with 70% alcohol between animals to maintain pathogen-free status.

Sciatic-tibial motor NCV is determined by stimulating distally at the sciatic notch and distally at the ankle via bipolar electrodes with supramaximal stimulation. The conduction velocity is calculated using two the points of stimulation along the nerve and measuring the resultant onset latency and distance.

Sensory NCV is determined by stimulating the sural nerve distally at the ankle via bipolar electrodes with supramaximal stimulation and recording at the fourth and fifth digit. The conduction velocity is calculated using the onset latency and distance.



Tail motor latency is determined by stimulating distally along the tail at a recorded distance of 3 cm. The onset latency is used for the latency measurement.

Tail sensory NCV is determined by stimulating proximally along the tail at a recorded distance of 3 cm. The conduction velocity is calculated using the onset latency and distance.