

Feb 26, 2019

# Neural aggregate formation

 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.x84fryw](https://dx.doi.org/10.17504/protocols.io.x84fryw)

Celeste M M. Karch<sup>1</sup>, Rita Martinez<sup>1</sup>, Jacob Marsh<sup>1</sup>


<sup>1</sup>Washington University in St Louis

Neurodegeneration Method Development Community  
Tech. support email: [ndcn-help@chanzuckerberg.com](mailto:ndcn-help@chanzuckerberg.com)



**Celeste M M. Karch**

Washington University in St Louis



OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.x84fryw](https://dx.doi.org/10.17504/protocols.io.x84fryw)

**Protocol Citation:** Celeste M M. Karch, Rita Martinez, Jacob Marsh 2019. Neural aggregate formation. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.x84fryw>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 17, 2019

**Last Modified:** February 26, 2019

**Protocol Integer ID:** 20476



## Attachments



IPSC CORTICAL

DIFFER...

179KB

## Guidelines

This protocol is part of the IPSC CORTICAL DIFFERENTIATION collection.

This method should be performed using sterile technique.













## Materials

Please refer to the attached full manuscript for required materials.

## Safety warnings




- ⚠ Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.



- 1 Harvest iPSCs for neural aggregate formation when iPSCs have reached 75- 85% confluency. Aspirate medium and rinse with  2 mL of DPBS.
- 2 Add  1 mL of Accutase. Incubate at  37 °C for  00:05:00 . Gently tap plate to dislodge cells.
- 3 Dilute Accutase with  4 mL of DMEM/F12 medium and collect cell suspension in 15ml conical tube.
- 4 Centrifuge cells at 750 rpm for  00:03:00 . Then carefully aspirate medium from iPSC pellet.
- 5 Add  3 mL of neural induction medium to iPSC pellet. Using a hemacytometer, count iPSCs. Adjust volume of iPSC suspension to 450-650,000 cells/mL using neural induction medium supplemented Rock inhibitor (10  $\mu$ M final).
- 6 Add  100  $\mu$ L of iPSC suspension per well to a v-bottom 96-well plate.
- 7 Centrifuge plate at 750 rpm for  00:03:00 to sediment iPSC into spheres.
- 8 Incubate cells at  37 °C , 5% CO<sub>2</sub> and 95% humidified chamber for  24:00:00 . After 24 hrs, carefully remove all medium from well and replace with  100  $\mu$ L per well of Neural Induction Medium.

#### Note

Do not disturb or break apart spheres. The spheres are very delicate at this stage.

- 9 Incubate neurospheres in 96 well plate for  96:00:00 . Perform half volume medium changes daily (removed  50  $\mu$ L and replace with  50  $\mu$ L ).