

Feb 27, 2019

# Neural progenitor banking

 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.x89frz6](https://dx.doi.org/10.17504/protocols.io.x89frz6)

Celeste M M. Karch<sup>1</sup>, Rita Martinez<sup>1</sup>, Jacob Marsh<sup>1</sup>


<sup>1</sup>Washington University in St Louis

Neurodegeneration Method Development Community  
Tech. support email: [ndcn-help@chanzuckerberg.com](mailto:ndcn-help@chanzuckerberg.com)



**Celeste M M. Karch**

Washington University in St Louis



OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.x89frz6](https://dx.doi.org/10.17504/protocols.io.x89frz6)

**Protocol Citation:** Celeste M M. Karch, Rita Martinez, Jacob Marsh 2019. Neural progenitor banking. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.x89frz6>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 18, 2019

**Last Modified:** February 27, 2019

**Protocol Integer ID:** 20481



## Attachments



IPSC CORTICAL

DIFFER...

179KB

## Guidelines


This protocol is part of the IPSC CORTICAL DIFFERENTIATION collection.

This method should be performed using sterile technique.

## Materials

Please refer to the attached full manuscript for required materials.





## Safety warnings

 Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.



- 1 Upon reaching at least 85% confluency, harvest neural progenitor cells as described in protocol below.

**Neural progenitor expansion protocol**

- 2 Perform a cell count in  3 mL of NIM using a hemacytometer.
- 3 Add equal volumes of NIM and 2x neural freezing medium to the NPC cell suspension for a final  $1 \times 10^6$  cells/mL.
- 4 Gently mix solution and distribute  1 mL into sterile cryovials. Store cryovials in Styrofoam containers at  -80 °C for  48:00:00 and then transfer to liquid nitrogen for long-term storage.