



Jan 26, 2026

Neural Progenitor Cell (NPC) differentiation

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvje1nwgk5/v1>

Lexi Bounds¹, Alejandro Barrera¹, Maria ter Weele¹, Siyan Liu¹, Euphy Wu², Shengyu Li¹, Revathy Venukuttan¹, Ruhi Rai¹, Wancen Mu², Nahid Iglesias¹, Paola Giusti-Rodriguez³, Timothy E. Reddy¹, Yun Li², Raluca Gordan⁴, Andrew Allen¹, Michael Love², Patrick Sullivan², Greg Crawford¹, Charles Gersbach¹

¹Duke University; ²University of North Carolina; ³University of Florida; ⁴University of Massachusetts

Gersbach Lab



Andrea R Daniel

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvje1nwgk5/v1>

Protocol Citation: Lexi Bounds, Alejandro Barrera, Maria ter Weele, Siyan Liu, Euphy Wu, Shengyu Li, Revathy Venukuttan, Ruhi Rai, Wancen Mu, Nahid Iglesias, Paola Giusti-Rodriguez, Timothy E. Reddy, Yun Li, Raluca Gordan, Andrew Allen, Michael Love, Patrick Sullivan, Greg Crawford, Charles Gersbach 2026. Neural Progenitor Cell (NPC) differentiation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n2bvje1nwgk5/v1>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 30, 2025

Last Modified: January 26, 2026

Protocol Integer ID: 231163

Keywords: iPSC, NPC, Neural progenitor cells, differentiations, differentiation into neural progenitor cell, neural progenitor cell, npc, cell, differentiation this protocol, ipsc culture, differentiation

Funders Acknowledgements:

NIH

Grant ID: HG012053

Disclaimer

This protocol was adapted from the work of Lexi Bounds in the Gersbach lab at Duke University.

Abstract

This protocol describes methods for WTC11 iPSC culture and differentiation into neural progenitor cells.



Materials

WTC11 iPSCs from the Coriell Institute for Biomedical Research (GM25256)

mTeSR1 (StemCell Technologies #85850)

Matrigel (Corning #354230)

Versene (ThermoFisher #15040066)

KnockOut™ Serum Replacement (ThermoFisher #10828028)

Y-27632 (Dihydrochloride) (StemCell Technologies (#72302)

6-well culture plates (Corning #3471)

SB431542 (StemCell Technologies #72234): Resuspend 10 mg in 2.06 mL DMSO.

LDN193189 (StemCell Technologies #72147): Resuspend 10 mg in 2.09 mL of DMSO.

STEMdiff™ Neural Rosette Selection Reagent (StemCell Technologies #05832)

Penicillin-Streptomycin (10,000 units/mL) (ThermoFisher #15140122)

B27 Supplement w/o Vitamin A (50X) (ThermoFisher #12587010)

N-2 Supplement (100X) (ThermoFisher #12587010)

Animal-Free Recombinant Human EGF (Peprotech #AF-100-15): Resuspend 100 ug in 1 mL ultra-pure H₂O.

Animal-Free Recombinant Human FGF-basic (Peprotech #AF-100-18B): o Resuspend 50 ug in 500 uL ultra-pure H₂O.

- Neural basal medium (NBM): 500 mL DMEM/F-12 + 5 mL N2 + 10 mL B27 + 5 mL P/S
- Neural induction medium (NIM): 50 mL NBM + 0.5 uL LDN + 50 uL SB431542
- Neural expansion medium (NEM): 50 mL NBM + 10 uL bFGF + 10 uL EGF

iPSC culture

- 1 All cell lines were cultured at 37C and 5% CO₂.
- 2 WTC11 iPSCs were obtained from the Coriell Institute for Biomedical Research (GM25256)
- 3 Cells are propagated in mTeSR1 (StemCell Technologies #85850) on Matrigel (Corning #354230)-coated tissue culture plates and passaged as colonies using Versene (ThermoFisher #15040066).
- 4 Cells are frozen in freezing solution comprised of 60% mTeSR1 media, 10% DMSO, and 30% KnockOut™ Serum Replacement (ThermoFisher #10828028), with 10uM ROCK Inhibitor (RI) Y-27632 (Dihydrochloride) (StemCell Technologies (#72302)).

NPC Differentiation

- 5 WTC11 iPSCs are differentiated into NPCs using embryoid-body (EB) formation and small molecules as previously described with minor modifications (Refs 1, 2)
- 6 WTC11 iPSCs are seeded in low-adhesion 6-well culture plates (Corning #3471) in NBM containing 10 uM SB431542 (StemCell Technologies #72234) and 0.1 uM LDN193189 (StemCell Technologies #72147) (NIM) with 10 uM Y-27632.
- 7 24 hours later (Day 1), 2 mL of media is removed and replaced with 2 mL NIM.
- 8 On Days 2-4, 3 mL of media is removed and replaced with 3 mL NIM.
- 9 On Day 5, EBs are transferred onto Matrigel-coated 6-well plates containing 1 mL NIM per well.
- 10 24 hours later (Day 6), all media is removed and replaced with 3 mL NIM.
- 11 On Days 7-12, 3 mL NIM was changed daily.



12 On Day 13, neural rosettes are harvested using STEMdiff™ Neural Rosette Selection Reagent according to the manufacturer's protocol (StemCell Technologies #05832) and plated onto Matrigel-coated plates in NEM ('Passage 0'):

- 12.1
1. Aspirate media.
 2. Wash w/ 1 mL DMEM/F-12.
 3. Aspirate.
 4. Add 1 mL neural rosette selection reagent.
 5. Incubate 1 hour at 37C. (may need longer incubation, check at 30min, 45min)
 6. Aspirate selection reagent.
 7. Using P1000, wash each well 2x w/ 1 mL DMEM/F-12 and collect cells.
 - a. Try not to break up the rosettes.
 8. Centrifuge 200-300g for 5min.
 9. Aspirate supernatant.
 10. Resuspend cells in 2 mL NEM.
 11. Plate to one well of 6-well.
- Repeat for all wells.

13 NPCs are propagated for five passages in NEM before characterization as described below. All experiments are performed with NPCs of passage number 7-15.

NPC culture

- 14 WTC11 NPCs are cultured in NEM. Change media daily with 2 mL NEM. Passage cells w/ Accutase when 80-90% confluent (include RI in media when passaging, like iPSCs) Seed cells at 200-500K cells per well (cell line/differentiate dependent). Collect cells for RNA at P2, P5. Freeze cells when possible in NBM + KSR + DMSO + RI.



Protocol references

1. Generation, Expansion, and Differentiation of Human Pluripotent Stem Cell(hPSC) Derived Neural Progenitor Cells (NPCs). Brafman DA. *Methods Mol Biol.* 2015;1212:87-102. doi: 10.1007/7651_2014_90.PMID: 25063499
2. Inhibition of STEP(61) ameliorates deficits in mouse and hiPSC-based schizophrenia models. Xu J, Hartley BJ, Kurup P, Phillips A, Topol A, Xu M, Ononenyi C, Foscue E, Ho SM, Baguley TD, Carty N, Barros CS, Müller U, Gupta S, Gochman P, Rapoport J, Ellman JA, Pittenger C, Aronow B, Nairn AC, Nestor MW, Lombroso PJ, Brennand KJ. *Mol Psychiatry.* 2018 Feb;23(2):271-281. doi: 10.1038/mp.2016.163. Epub 2016 Oct 18.PMID: 27752082
3. Transcriptional signatures of schizophrenia in hiPSC-derived NPCs and neurons are concordant with post-mortem adult brains. Hoffman GE, Hartley BJ, Flaherty E, Ladrán I, Gochman P, Ruderfer DM, Stahl EA, Rapoport J, Sklar P, Brennand KJ. *Nat Commun.* 2017 Dec 20;8(1):2225. doi: 10.1038/s41467-017-02330-5. PMCID: PMC5738408