

Feb 26, 2019

# Neural progenitor expansion

 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.x88frzw](https://dx.doi.org/10.17504/protocols.io.x88frzw)

Celeste M M. Karch<sup>1</sup>, Rita Martinez<sup>1</sup>, Jacob Marsh<sup>1</sup>


<sup>1</sup>Washington University in St Louis

Neurodegeneration Method Development Community  
Tech. support email: [ndcn-help@chanzuckerberg.com](mailto:ndcn-help@chanzuckerberg.com)



**Celeste M M. Karch**

Washington University in St Louis



OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.x88frzw](https://dx.doi.org/10.17504/protocols.io.x88frzw)

**Protocol Citation:** Celeste M M. Karch, Rita Martinez, Jacob Marsh 2019. Neural progenitor expansion. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.x88frzw>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 17, 2019

**Last Modified:** February 26, 2019

**Protocol Integer ID:** 20480



## Attachments



IPSC CORTICAL

DIFFER...

179KB

## Guidelines

This protocol is part of the IPSC CORTICAL DIFFERENTIATION collection.

This method should be performed using sterile technique.

## Materials

Please refer to the attached full manuscript for required materials.

## Safety warnings

- ⚠ Please refer to the SDS (Safety Data Sheet) for information about hazards, and to obtain advice on safety precautions.

- 1 Using the pre-coated PLO-laminin plate, aspirate laminin from 2-3 wells. Remove supernatant from 15 mL conical tube containing 2 wells of neural rosette clusters. See protocol below.

### Protocol




NAME

## Neural Rosette Formation and Selection





CREATED BY





Celeste M M. Karch

**PREVIEW**

- 1.1 On Day 5 of neural aggregate formation, remove media (by pipetting) and carefully wash spheres with  100  $\mu$ L of pre-warmed DMEM/F12. Repeat 2 times.







### Note

Do not break apart spheres. Neural spheres are very delicate at this stage. An alternative approach is to remove  50  $\mu$ L of spent media and wash with  50  $\mu$ L DMEM/F12 . Add  50  $\mu$ L fresh neural induction media . Transfer  100  $\mu$ L of spheres and media to the new PLO/laminin-coated well. This approach will transfer more dead cells into the new well.

- 1.2 Remove the last wash and add  50  $\mu$ L of neural induction media to each well.
    - 1.3 Aspirate laminin from one well of the pre-coated plate. Using 200  $\mu$ L sterile tips, carefully pipet up spheres from wells using  100  $\mu$ L volume and transfer thirty-two spheres per well. Repeat above steps for the remaining wells. Incubate cells in  37  $^{\circ}$ C , 5% CO<sub>2</sub> and 95% humidified chamber and distribute evenly by making a "T" motion.
    - 1.4 After  24:00:00 , examine attached aggregates. Remove medium and replace with 2mLs/well fresh neural induction medium daily.






**Note**

If some aggregates have not attached, carefully pipet out all medium and replace with 1 ml/well fresh neural induction medium. Once 90-100% of aggregates attach, exchange medium daily with 2 mls/well neural induction medium.









- 1.5 Monitor spheres daily under microscope for formation of neural rosette structures. Neural rosettes are ready to harvest when spheres have completely flattened and clusters are clearly visible (3-7 days after plating, line dependent).
- 1.6 Harvest neural rosettes by aspirating spent medium. Add  1 mL of pre-warmed DMEM/F12 to each well to remove unattached cells (repeat if necessary).
- 1.7 Add  1 mL of Neural Rosette Selection reagent to each well and incubate for up to  01:00:00 at  37 °C (check cells at  00:20:00 . Cells are typically collected after 30-45 min incubation. Look for rosette structure to be rounding up without the disturbance of other surrounding cells).
- 1.8 Carefully remove Neural Rosette Selection reagent with a pipet, being careful not to disturb rosette clusters. Add  1 mL DMEM/F12 to each well, then using a p1000 detach rosette clusters by rinsing over them.
- 1.9 Transfer rosette material from 1 well into a 15 mL conical tube for cryopreservation of neural rosettes and from 2 wells into a separate 15 mL conical tube for neural progenitor expansion. Do not triturate clusters.

**Note**

To maintain a pure culture, it is best to leave some rosettes behind rather than collect all of the rosettes and additional cells.


- 1.10 Centrifuge rosette clusters at 750 rpm for  00:03:00 .
- 2 Add appropriate amount of neural induction media (NIM) to achieve a final volume of  2 mL per well. Pipet clusters up and down breaking them into 1/4 or 1/5 the size of the rosette cluster.
- 3 Add  2 mL of neural rosette cluster suspension to each well and incubate at  37 °C , 5% CO<sub>2</sub> and 95% humidified chamber for  24:00:00 .



- 4 Examine adherent cells under microscope, aspirate off medium and replace with  2 mL of fresh neural induction medium daily for 3-5 days or until they reach ~80% confluent.
- 5 Make  12 mL of DMEM/F12 supplemented with 10% FBS by adding  2 mL of FBS to  10 mL of DMEM/F12 to inhibit trypsin activity (termed "Complete Media").
- 6 Remove medium and rinse cells with 1 mL/well of DMEM/F12 .
- 7 Harvest cells by adding 1 mL/well of 0.05% trypsin and incubate at  37 °C for  00:03:00 .
- 8 Inhibit trypsin activity by adding  4 mL of DMEM/F12 supplemented with 10% FBS to each well.
- 9 Collect cells in 15mL conical tube. Centrifuge at 750 rpm for  00:03:00 .


**Note**

To minimize bubbles and increase cell yield, add 1mL Complete Media to 15mL conical tube. After trypsinization, add 3mL DMEM/F12 to collect cells and transfer to 15mL conical tube.

- 10 Aspirate supernatant and resuspend in  12 mL of NIM, pipet at least 3 times to break up large clumps.

**Note**

Small clumps and single cells are acceptable to passage.

- 11 Passage cells to pre-coated PLO-lamin plates (see above) by adding 2 mLs/well of suspended cells for a total of 6 wells.
- 12 Feed with  2 mL of fresh NIM daily until they reach 85-95% confluency.



- 13 **IMPORTANT:** NPCs are best maintained by passaging at a 1:3 dilution. Passaging cells that are more dilute will result in slow growth. Cells will typically become confluent after 2-4 days when plated as a 1:3 dilution. If NPCs are observed to slow in their growth, passage cells and replate without dilution.