



Oct 01, 2021

NGM Agar



In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.zy7f7zn

Adrien Assie¹, Buck Samuel¹

¹Baylor College of Medicine

Adrien Assie: I am not the author of this protocol. This is a standard protocol for C. elegans maintenance available on wormbook.



Adrien Assie

Baylor College of Medicine

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.zy7f7zn

External link: http://www.wormbook.org/chapters/www_strainmaintain/strainmaintain.html

Protocol Citation: Adrien Assie, Buck Samuel 2021. NGM Agar. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.zy7f7zn>

Manuscript citation:

Brenner, S. (1974). Genetics 77, 71.

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 11, 2019

Last Modified: October 01, 2021

Protocol Integer ID: 22271

Keywords: C. elegans, Growing media, Caenorhabditis , nematode

















Disclaimer

This is a standard protocol for *C. elegans* maintenance available on wormbook.

Abstract

NGM agar recipe from wormbook



- 1 Start with  975 mL water
- 2  3.0 g NaCl
- 3  2.5 g Peptone
- 4  17 g Agar
- 5 Autoclave with stir bar
- 6 Cool to  55 °C (faster in water bath)
- 7 Then add the following while stirring on heat plate (DO NOT OVERHEAT)
- 8  0.5 mL of  1 Molarity (M) CaCl₂ (sterile)
- 9  1 mL of 5 mg/mL Cholesterol (dissolved in ethanol)
- 10  1 mL of  1 Molarity (M) MgSO₄ (sterile)
- 11  25 mL of  1 Molarity (M) Potassium Phosphate Buffer, pH 6.0 (sterile)
- 12 Pour into petri dishes using sterile technique
 - 60 mm dishes = About  7 mL per plate
 - 100 mm dishes = About  20 mL per plate

