

Jul 20, 2020

Nickel-NTA Protein Purification

DOI

dx.doi.org/10.17504/protocols.io.bihhkb36

Nay Chi Khin¹, Diep Ganguly², Suyan Yee¹, Kai Xun Chan¹

¹The Australian National University; ²University of Pennsylvania

Pogson Group



Diep R Ganguly

The Australian National University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bihhkb36>

Protocol Citation: Nay Chi Khin, Diep Ganguly, Suyan Yee, Kai Xun Chan 2020. Nickel-NTA Protein Purification. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bihhkb36>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: July 13, 2020

Last Modified: July 20, 2020

Protocol Integer ID: 39177

Keywords: Ni-NTA, Protein purification, nta protein purification protocol, tagged recombinant protein, recombinant protein, purification protocol, using ni, protein,

Abstract

Protocol to purify 6xHis-tagged recombinant proteins expressed in *E. coli* using Ni-NTA pull-down.

Guidelines

Try to perform all steps at 4 °C as best as possible and keep protein samples on ice. For all spin-concentrate and buffer exchange steps, use slow speed (e.g. 3500-5000 rcf) to prevent precipitation.

Materials

MATERIALS

☒ Sodium dihydrogen phosphate **P212121**

☒ Disodium hydrogen phosphate

☒ Triton-X100

☒ 1.5 mL Eppendorf tubes

☒ Glycerol **Bio Basic Inc. Catalog #GB0232.SIZE.500ml**

☒ Poly-Prep Chromatography Columns **Bio-Rad Laboratories Catalog #731-1550**

☒ DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #43816-10ML**

☒ Sodium chloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S3014**

☒ Imidazole **Merck MilliporeSigma (Sigma-Aldrich) Catalog #I5513**

☒ Bovine Thrombin Protein **Thermo Fisher Catalog #1996000JL**

☒ HisPur® Ni-NTA Resin **Thermo Fisher Catalog #88221**

☒ Pierce® HRV 3C Protease Solution Kit (2 units/μL) **Thermo Fisher Catalog #88946**

☒ TEV protease **New England Biolabs Catalog #P8112S**

☒ Tris Hydrochloride (Tris-HCl) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #RES3098T-B7**

☒ Amicon® Ultra-4 Centrifugal Filter Unit **Merck Millipore (EMD Millipore) Catalog #UFC8003**

☒ Amicon® Ultra-15 Centrifugal Filter Unit **Merck Millipore (EMD Millipore) Catalog #UFC9003**

Troubleshooting

Before start

<https://www.protocols.io/private/3797D176057F0C37477116D259C75682>

Prepare buffers

1h

1 Buffer recipes

1h

Resuspension buffer

- 50 mM Tris-HCl, pH 8.0
- 2 mM EDTA

1X Binding buffer

- 50 mM sodium phosphate buffer, pH 8.0 (ensure >1 pH units away from pI of expressed protein)
- 500 mM sodium chloride
- 0.5% Triton X-100
- 10% glycerol
- 10 mM imidazole

N.b. This should reflect the cell lysis buffer (e.g.

<https://www.protocols.io/private/3797D176057F0C37477116D259C75682>).

1X Wash buffer

- 50 mM sodium phosphate buffer, pH 8.0
- 500 mM sodium chloride
- 0.5% Triton X-100
- 10% glycerol
- 20 mM imidazole

1X Elute buffer

- 50 mM sodium phosphate buffer, pH 8.0
- 500 mM sodium chloride
- 0.5% Triton X-100
- 10% glycerol
- 250 mM imidazole

N.b. DTT cleaves/interferes with Ni^{2+} -binding, use low DTT concentrations (<2 mM), in the above buffers, if needed.

Digestion buffer

- 50 mM Tris-HCl, pH 8.0
- 150 mM NaCl

For *storage*, supplement *digestion buffer* with:

- 1 mM DTT
- 20% glycerol



Ni-NTA Affinity Purification

35m

- 2 Keep solubilised protein samples (from completion of <https://www.protocols.io/private/3797D176057F0C37477116D259C75682>) on ice to thaw completely. 15m
- 3 Centrifuge samples (14,000 rcf for 10-20 minutes) and syringe filter (0.45 µm) supernatant into clean microfuge tube. Keep on ice. 30m
- 4 Add an appropriate amount of Ni-NTA resin slurry (50% slurry in 20% ethanol) into gravity-flow column. For 15 mL bacterial cultures resuspended in 2 mL lysis buffer, equating to ~1-4 mg protein (1-2 mg/mL), 500 µL resin slurry is added.

Top-up column with 5x volume of *1X binding buffer* (e.g. 500 µL slurry = 250 µL resin = ~2 mL 1X binding buffer).

Allow binding buffer to drop into a waste bottle to equilibrate the resin.
- 5 Transfer syringe filtered soluble protein (i.e. supernatant from step 3) into equilibrated resin, cap column tightly at both ends (including parafilm to prevent leakage) and incubate with shaking for 1-2 hours @ 4°C (max 10 mL per column). 2h
- 6 Keep column up-right and allow resin to settle. Once settled, open cap and capture flow-through (FT).
- 7 Wash resin three-times with *1X Wash buffer* (2x resin volume, i.e. 1 mL) and collect separately (i.e. W1, W2, W3).
- 8 Elute four-times with *1X Elute buffer* (1X resin volume, i.e. 500 µL), collect separately (i.e. E1, E2, E3, E4).
- 9 Run SDS-PAGE of pre-FT (20 µL), FT (20 µL), W1-W3 (20 µL), and E1-E4 (8 µL) to quality-check purification.
- 10 Store all samples @ 4 °C for the short-term. Keep columns containing resin for second Ni-NTA step.

Protease digestion

16h

- 11 Pool elutes containing protein fragment of expected size.



- 12 Perform buffer exchange via centrifuge filtration (or dialysis) into *Digestion buffer* to reduce [imidazole] < 1mM.

For example, if *1X Elution buffer* contains 250 mM imidazole, you will need 4x centrifugation at a dilution of 1:4 per spin (4 °C @ 3,500 rcf for 20 minutes).
- 13 Add appropriate protease (e.g. TEV, Thrombin, HRV 3C) and incubate overnight (or > 16hrs) at 4 °C. Ensure to add any required cofactors (e.g. citrate for TEV).
- 14 Equilibrate Ni-NTA resin (from previously used and stored columns) with 10-20 mL *Digestion buffer* and allow resin to settle. Let buffer flow-through.
- 15 Add digested protein to resin and collect flow-through (dig-FT). Can perform multiple elutes with *Digestion buffer* to obtain as much protein as possible.
- 16 Run elutes on SDS-PAGE gel to test for cleavage.
- 17 Pool and concentrate the pure fraction(s) using buffer exchange.
- 18 Prepare aliquots of purified protein supplemented with 10-20% glycerol and 1mM DTT. Store at -20 °C.

16h