


Jan 12, 2024

Nigrostriatal organotypic cultures to study neuromelanin accumulation in dopaminergic circuits

DOI

dx.doi.org/10.17504/protocols.io.kxygx381wg8j/v1



Mario Fernandez-Ballester¹, Federico N. Soria¹

¹Achucarro Basque Center for Neuroscience



Federico N. Soria

Achucarro Basque Center for Neuroscience

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx381wg8j/v1>

Protocol Citation: Mario Fernandez-Ballester, Federico N. Soria 2024. Nigrostriatal organotypic cultures to study neuromelanin accumulation in dopaminergic circuits. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygx381wg8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: January 08, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 93056

Keywords: ASAPCRN, organotypic culture, dopaminergic neurons, neuromelanin, organotypic cultures from parasagittal brain slice, neuromelanin accumulation in dopaminergic circuit, neuromelanin accumulation, nigrostriatal organotypic culture, survival of dopaminergic neuron, integrity of the nigrostriatal pathway, dopaminergic neuron, nigrostriatal pathway, organotypic culture, basic cytoarchitecture of the brain, parasagittal brain slice, dopaminergic degeneration, dopaminergic circuit, neuron, extracellular matrix, brain, cell

Funders Acknowledgements:

Aligning Science Across Parkinson

Grant ID: ASAP-020505

Abstract

In this protocol we describe the preparation and maintenance of rat organotypic cultures from parasagittal brain slices. We use a 13° slicing angle and a customized cocktail of small molecules and growth factors, to maximize the integrity of the nigrostriatal pathway and ensure the survival of dopaminergic neurons, respectively. The slices maintain the basic cytoarchitecture of the brain, including glia and the extracellular matrix (not including vessels nor immune system). The cultures can be used to study dopaminergic degeneration, cell-cell and cell-matrix interactions and are particularly suitable for AAV-mediated overexpression of transgenes, time-lapse live imaging and longitudinal studies.



Materials

Pro-dopaminergic Organotypic Medium

- Neurobasal-A (Gibco)
- 1X B27 (Gibco)
- 1X Glutamax (Gibco)
- 50 U/mL Penicillin-Streptomycin (ThermoFisher)
- 1.25 ug/mL Amphotericin B (ThermoFisher)
- *200 μ M Ascorbic Acid (Sigma)*
- *100 μ M Bucladesine (cAMP, MedChemExpress)*
- *2.5 μ M DAPT (Stemcell)*
- *2 ng/mL BDNF (Stemcell)*
- *2 ng/mL GDNF (Stemcell)*
- *1 ng/mL TGF-beta (MedChemExpress)*

The factors in italics have to be added fresh at every medium change (3 times a week)

Other materials and reagents

- Millipore Cell Culture Inserts (0.4 μ m pore size, 30 mm diameter) (Merck-Millipore #PICM0RG50)
- HBSS without Ca^{2+} and Mg^{2+}
- Low-melting point agarose (CondaLab #8050)
- AraC (Cytosine-beta-D-arabinofuranoside, Sigma-Aldrich #C1768)
- p60 petri dish
- p100 petri dish
- p35 petri dish
- Dissection tools
- Surgical blades
- Spatula

McIlwain Tissue Chopper

Equipment

McIlwain Tissue Chopper	NAME
Cavey Laboratory Engineering	BRAND
51350	SKU




Troubleshooting

Culture preparation

- 1 Prepare a 6-well culture plate with **Millipore 0.4 μ m culture inserts** (see Materials) placed on top of 1 mL of **Pro-dopaminergic organotypic medium** (see Materials) per well. Incubate it at  37 °C and 5% CO₂ to warm up.

Note

It is important to minimize air bubbles between the medium and the membrane insert, as this will prevent access of nutrients to the slice.


- 2 Prepare a p100 petri dish with 25 mL of Hank's Balanced Salt solution without Ca²⁺ and Mg²⁺ (ThermoFisher) supplemented with D-Glucose 36 mM (Gibco) and a p60 petri dish with 5 mL of the same dissection solution, and keep them  On ice .

- 3 Desinfect all surfaces and dissection tools. Prepare the **McIlwain Tissue Chopper** (see Materials) with a fresh blade and clean holder plates.

- 4 Euthanize P10-P12 rat pups by quick decapitation using large scissors. Separate the head skin with a scalpel and expose the skull. Cut the skull through the middle line with fine scissors. Use forceps to pull apart the skull and expose the brain.

Note

We prefer P10-P12 (contrary to the typical P5-P7 ages used for cortico-striatal slices) to maximize both the volume of the substantia nigra and dopaminergic innervation of the striatum (Oo & Burke, 1997).


- 5 Remove the brain with a spatula and transfer to a filter paper to remove moist. Place the brain carefully into a p35 petri dish with a warm 4% solution of LM-Agarose (see Materials). Temperature of the agarose solution should be no more than  37 °C .

- 6 Let the agarose gellify  On ice for a few minutes.

- 7 As soon as the agarose is stiff, crop with a blade a block of agarose containing the brain. Paste the block with cyanocrilate glue to the holder plate so the ventral part is in contact

with the plate.

- 8 Position the plate in the chopper with the brain midline oriented 13° from the blade (see Image).



13° angle sagittal sectioning

Obtain sagittal slices of $350\ \mu\text{m}$ width at 13° angle. When the blade arrives to the midline, rotate the holder plate so the blade is aligned with the rightmost line and section again.

Expected result

This method of sectioning will render sagittal slices containing the full nigrostriatal pathway.


- 9 Separate the agarose block from the plate with a scalpel, avoiding damage to the ventral part. Pour the slices into the p100 petri dish with HBSS (from step 2).
- 10 Under the dissection microscope, remove the agarose and separate the slices with a couple of spatulas (preferentially plastic-ones, which can be custom crafted cutting the back of a Pasteur pipette). Only the 2nd or 3rd slice lateral from the midline contain the full nigrostriatal pathway. Transfer the selected slices to the p60 petri dish for to wash them (you can use an inverted glass pipette with a suction rubber bulb to minimize damaging of the slices).

**Note**

Cortex and hippocampus can be removed with a scalpel to reduce the amount of tissue to be cultured.

Expected result

From each brain, 2-4 slices (1-2 per hemisphere) with the full nigrostriatal pathway can be obtained.

- 11 Transfer the slice to the preheated culture inserts with medium using the inverted glass pipette with a suction rubber bulb. Each insert can contain only one nigrostriatal slice. Remove excess HBSS with a pipette and incubate the slices at  37 °C and 5% CO₂.

Expected result

To maximize slice survivability, no more than 15-20 min should pass from euthanasia to placing the slices into the incubator.

Maintenance

- 12 At DIV1 (24h after culture preparation), replace the feeding medium with 1 mL of the same medium but adding antimetabolic 4.4 uM AraC (Sigma-Aldrich).

Note

AraC will prevent the overproliferation of astrocytes, minimizing the formation of a glial scar on top of the slice.

- 13 At DIV3, replace the AraC-containing medium with fresh medium and add 1 uL of AAV9-CMV-hTyr (10^{12} vg/mL) over the substantia nigra. Do not touch the tissue.

Note

To improve spatial precision (e.g. when infecting a selected region from the slice), AAV infection can be done with a micromanipulator and a glass pipette.



- 14 For culture maintenance, change medium 3-times per week (small molecules and growth factors have to be added fresh in each medium change). Slices can survive in culture more than 30 days.

Expected result

Slices will not flatten as much as cortico-striatal slices.

Neuromelanin observation

- 15 At 10 DIV, neuromelanin spots will appear in the substantia nigra as brown spots visible with the naked eye.

Expected result

The AAV will infect any cell (this can be visualized with an anti-human tyrosinase antibody), but **neuromelanin will be produced only in catecholaminergic-producing regions and restricted to neurons**. This is why it is important to maximize the survivability of dopaminergic (TH-positive) neurons.

- 16 At 15 DIV, the hTyr expression is maximal and neuromelanin accumulation can be seen all over the substantia nigra and connected regions (the nigrostriatal pathway becomes brown). Since the slices do not flatten and cell density is higher, any fluorescence imaging should be done in *confocal* or *two-photon* setups where optical sectioning can be achieved. Neuromelanin should be ideally imaged with a color camera, to ensure that the accumulate is brown (neuromelanin) and not black/gray (cell death).

Protocol references

Oo, T. F., & Burke, R. E. (1997). The time course of developmental cell death in phenotypically defined dopaminergic neurons of the substantia nigra. *Brain research. Developmental brain research*, 98(2), 191–196. [https://doi.org/10.1016/s0165-3806\(96\)00173-3](https://doi.org/10.1016/s0165-3806(96)00173-3)