



Feb 10, 2026

Nile red staining protocol to quantify PHB in *Saccharomyces cerevisiae* cultivations

 [PLOS Computational Biology](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8xwxrv2w/v1>

Anna Ylinen¹, Tuula Tenkanen¹

¹VTT Technical Research Centre of Finland Ltd



Tuula Tenkanen

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8xwxrv2w/v1>

External link: <https://doi.org/10.1371/journal.pcbi.1014087>

Protocol Citation: Anna Ylinen, Tuula Tenkanen 2026. Nile red staining protocol to quantify PHB in *Saccharomyces cerevisiae* cultivations. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5jyl8xwxrv2w/v1>

**Manuscript citation:**

Tenkanen T, Ylinen A, Jouhten P, Penttilä M, Castillo S (2026) PHA synthase variant design using a conditional variational autoencoder. PLOS Computational Biology 22(3). doi: [10.1371/journal.pcbi.1014087](https://doi.org/10.1371/journal.pcbi.1014087)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 09, 2026

Last Modified: February 10, 2026

Protocol Integer ID: 242904

Keywords: PHB, PHA, Nile red, phb in *saccharomyces cerevisiae*, nile red staining protocol, nile red staining, other phas as nile, pha granule, surface of pha granule, *saccharomyces cerevisiae*, other pha, phb, red staining, producing microbe, nile, other intracellular structures such as lipid droplet

Funders Acknowledgements:

Nessling Foundation

Grant ID: 201700259, 201800005

Centre for Young Synbio Scientists (CYSS) funded by Jenny and Antti Wihuri Foundation

Abstract

Nile red staining can be used to screen PHB producing microbes as it binds on the surface of PHA granules. This protocol has been modified from Spiekermann *et al.* (1999) (DOI 10.1007/s002030050681) and Zuriani *et al.* (2013) (DOI 10.1007/s12257-012-0607-z). We recommend three or four technical replicates for each sample. As a negative control you will need a strain that does not produce PHB or other PHAs as Nile red do not only bind to the surface of PHA granules, but also to other intracellular structures such as lipid droplets.

Before start

Prepare Nile red solution with a concentration of 30 mg/ml in DMSO.



- 1 Dilute all cultivation samples to same OD with water for example in a 1.5 ml Eppendorf tube. OD
2-5 is preferred but range 0.5-10 works as well. You will need 100 µl of diluted sample/ measurement. Therefore, 400 µl + a little extra is needed for 4 replicate measurements.
- 2 Pipet 100 µl of each cell sample (all samples now at same OD) to a black 96 well plate.
- 3 Pipet 20 µl of the Nile red solution (30 mg/ml in DMSO) to each well and mix by pipetting. Incubate for 10 minutes in RT and measure the fluorescence with a fluorescence plate reader.

Protocol references

Spiekermann, P., Rehm, B.H.A., Kalscheuer, R., Baumeister, D., Steinbüchel, A., 1999. A sensitive, viable-colony staining method using Nile red for direct screening of bacteria that accumulate polyhydroxyalkanoic acids and other lipid storage compounds. *Archives of Microbiology* 171, 73–80. <https://doi.org/10.1007/s002030050681>

Zuriani, R., Vigneswari, S., Azizan, M.N.M., Majid, M.I.A., Amirul, A.A., 2013. A high throughput Nile red fluorescence method for rapid quantification of intracellular bacterial polyhydroxyalkanoates. *Biotechnology and Bioprocess Engineering* 18, 472–478. <https://doi.org/10.1007/s12257-012-0607-z>