



Feb 03, 2025

Nissl Staining with Thionin

DOI

<https://dx.doi.org/10.17504/protocols.io.261ger8edl47/v1>

Tracy Larson¹

¹University of Virginia



Tracy Larson

Department of Biology, University of Virginia

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ger8edl47/v1>

Protocol Citation: Tracy Larson 2025. Nissl Staining with Thionin. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261ger8edl47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 02, 2025



Last Modified: February 03, 2025

Protocol Integer ID: 119454

Keywords: nissl stain with thionin, nissl stain, selective staining, cell bodies of neuron, clear identification of neuron, staining process, components of nissl body, neuronal cytoplasm, studying neuronal morphology, neuronal morphology, neuron, cytoarchitecture of various brain region, microscope, clusters of rough endoplasmic reticulum, cell body, thionin solution, thionin, distribution within the tissue, primary dye, rough endoplasmic reticulum, cell density, histological technique, nissl body, various brain region, tissue section, excess stain, cytoarchitecture, tissue, purple structure, contrast, brain, dye

Abstract

A "Nissl stain with thionin" is a histological technique used to visualize the cell bodies of neurons in tissue sections, where the primary dye, thionin, binds to the negatively charged components of Nissl bodies (clusters of rough endoplasmic reticulum) within the neuronal cytoplasm. This selective staining allows for clear identification of neurons and their distribution within the tissue. The method is particularly useful for studying neuronal morphology, cell density, and the cytoarchitecture of various brain regions. During the staining process, tissue sections are exposed to a thionin solution, which stains the Nissl bodies, and the excess stain is often removed through a differentiating solution to enhance contrast. Under the microscope, neurons appear as distinct blue or purple structures, making it easy to distinguish them from non-neuronal tissue.



Materials

4% PFA

Paraformaldehyde 4.0g

1x PBS 100ml

**heat to get PFA into solution, then cool to 4C before use*

***don't boil!!*

10xPBS, phosphate buffered saline

To make 1 L of 10X stock:

NaCl 80.0 g

KCl 2.0 g

Na₂HPO₄ 14.4 g

KH₂PO₄ 2.4 g

nano H₂O to 1.0 L

**dilute to 1X working concentration*

Thionin Stock

Stock 1.3% thionin:

13 gm Thionin

1000 ml distilled H₂O

**Stir and heat gently for 1 hour. Filter the solution after the dye has dissolved. Store in a stoppered bottle. Use high purity thionin, such as Sigma T3387.*

Thionin Buffer

1 M Acetic Acid (HAc):

58.5 ml glacial acetic acid

dilute to 1 liter with distilled water

2. 1 M Sodium Hydroxide (NaOH):

50 gm sodium hydroxide pellets

dissolve to 1 liter with distilled water

**Mix buffer reagents together and adjust pH before mixing with the thionin stock in the formulas below.*

Thionin Stain

80.0 ml 1M HAc

14.4 ml 1M NaOH

305.6 ml Thionin stock

Tissue Fixation

- 1 Fix tissue in 4% PFA for 15 min
- 2 Rinse in 1x PBS at least once for 5 min

Staining

- 3 Rinse in dH₂O for 3-5 min
- 4 Stain in Thionin stain for 2-7 min, depending on number of times stain has previously been used
 - * Typical schedule: 1-3 uses = 2-3 min; 4-6 uses = 3-4 min; 6-10 uses = 5-7 min
 - **Re-make solution after 10 uses
- 5 Dehydrate:
 - 70% EtOH with a few drops of Acetic Acid, 15-30 sec, dip twice
 - 70% EtOH, 15-30 sec, dip twice
 - 95% EtOH, 30 sec to several minutes
 - *This step differentiates the nuclei stain (purple) from the cytoplasm stain (blue)
 - **Watch for desired "darkness" of stain (ideally medium blue with a slight hue of purple)
- 6 100% EtOH, 30 sec with several swishes back and forth
 - *Important to completely dehydrate/remove all water from the sections before moving to xylenes

Mounting Coverslips

- 7 Clear in xylenes for 3-5 min
 - *Look for clearing of "fuzzy" or "drippy" (this is moisture trapped in the tissue) appearance before mounting
- 8 Drop a dime-sized drop of DPX onto each slide and gently place a 24mmx60mm coverslip on top