

Sep 11, 2019

Nodule cleaning



In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.s4xegxn



Ming-Dao Chia

OPEN ACCESS



DOI: dx.doi.org/10.17504/protocols.io.s4xegxn

Protocol Citation: Ming-Dao Chia: Nodule cleaning. protocols.io https://dx.doi.org/10.17504/protocols.io.s4xegxn

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

These are protocols used in the Honours thesis of the same title: https://doi.org/10.25911/5d6cf6fb99437

Created: August 30, 2018

Last Modified: September 11, 2019

Protocol Integer ID: 15223

Abstract

Adapted for high throughput of small samples for sequencing.

Based on doi:10.1371/journal.pone.0047677



- 1 Cleanly cut nodules from plant roots and drop into a 2ml Eppendorf tube per plant. 2 Add 1ml 100% EtOH to each tube
- 3 Incubate for 5 seconds
- 4 Remove as much liquid as possible without removing loose nodules
- 5 Add 1ml 6% sodium hypochlorite
- 6 Remove as much liquid as possible
- 7 Add 2ml sterile water
- 8 Remove as much liquid as possible **≣)** go to step #7 Repeat once
- 9 Close tube, move to sequencing/culturing