


Oct 03, 2024

🌐 Non-destructive microplastic isolation from water and sediment samples

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l2q264l1y/v1>



Jonas Stadfeld¹, Sneha Suresh^{1,2}, Srijak Bhatnagar^{1,2}

¹University of Calgary; ²Athabasca University



Srijak Bhatnagar

Athabasca University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2q264l1y/v1>

Protocol Citation: Jonas Stadfeld, Sneha Suresh, Srijak Bhatnagar 2024. Non-destructive microplastic isolation from water and sediment samples. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l2q264l1y/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 23, 2024

Last Modified: October 03, 2024

Protocol Integer ID: 100443

Keywords: Microplastics, Plasticsphere, Microplastics Microbiome, Nile Red, Density Separation, comprehensive study of the microplastic microbiome, microplastic microbiome, sediment samples microplastics in aquatic ecosystem, sediment samples microplastic, microplastic, microbial, microbial community, microbe, diverse microbial community, metagenomic, associated microbiota, unique habitats for diverse microbial community, sediment sample, aquatic ecosystem, rna gene, nile red dye, integrity of the plastic, red dye

Abstract

Microplastics in aquatic ecosystems serve as unique habitats for diverse microbial communities, collectively referred to as the plasticsphere. Investigating these microbes using both culture-dependent and culture-independent techniques (e.g., 16S rRNA gene sequencing, metagenomics) requires a method that preserves the integrity of the plastic and associated microbiota. This protocol outlines a non-destructive isolation technique for microplastics from water and sediment samples using Nile Red dye, thus allowing for a comprehensive study of the microplastic microbiome and its functional potential without compromising the plastic matrix or the microbial communities it harbors.

Image Attribution

Jonas Stadfeld and Sneha Suresh

Troubleshooting



Nile Red Solution (2 µg/mL)

- 1 Add  50 mg of Nile Red to a Falcon Tube.
- 2 Using a serological pipette, carefully add  50 mL of molecular-grade acetone to the Falcon Tube.
- 3 Screw the Falcon tube cap and vortex vigorously to dissolve the Nile Red dye ( 00:00:30). 30s
- 4 Pick up the dissolved Nile red dye into a 50 mL syringe by pulling the plunger.
- 5 Filter sterilization: Attach a sterile 0.22 µm Luer-lock syringe filter to the syringe and push the Nile red into a new sterilized  50 mL falcon tube.
- 6 Label the tube as concentrated Stock solution ([M] 1 mg/mL), cover it in aluminum foil, and store at  4 °C
- 7 In a sterilized 500 mL Nalgene Glass media bottle (orange cap), add  499 mL of Milli Q water.
Tip: Add 500 mL using a measuring cylinder and remove 1 mL using a pipette.
- 8 Add  1 mL of filter-sterilized Nile red stock solution ([M] 1 mg/mL). Close the lid and shake gently.
Cover it in Aluminum foil and label.
- 9 Store in the dark at  4 °C for long-term storage.

Microplastic Isolation - Water

- 10 NOTE: Use only metal and glass equipment. Tubing should be silicone.
- 11 Filter at least 40L of water through a 20µm Cellulose filter (47 mm).

Note: Depending upon the murkiness, the filter may clog before 40L is filtered. In that case, multiple filters can be combined to form ONE replicate 40L sample.

The Filter can be stored at  -20 °C if used for DNA extraction.

12 Place the filter in a petri dish and add  1 mL of Nile Red solution ([MI] 2 µg/mL).

13 Incubate the filter in dark at  4 °C for  24:00:00

1d



14 Place the filter in a fluorescence microscope under a 529 nm wavelength filter



15 Using sterile forceps, pick the red fluorescing microplastics and transfer them to a 1.5 mL microcentrifuge tube containing  100 µL sterile 1X PBS.



Note: Nile Red can stain roots and other biological materials. In background noise, microplastics fluoresce strongly, lacking compartments and multicellular structure.

16 Optional step: Rinse the microplastics with  100 µL 1X PBS to remove environmental contaminants by pelleting the plastics with centrifugation (1000 xg for  00:01:00), discarding the supernatant and repeating the step twice with fresh sterile 1X PBS.

1m



Microplastic Isolation - Sediment

12h

17 Freeze dry the  0.5 kg  Sediment Sample at 80 °C to remove all water.



Note: The duration of freeze-drying varies based on the water content of the sample. Expect 24 hours to three days.

18 While the freeze dryer runs, prepare  900 mL of Calcium Chloride Solution (final density of 1.4 Kg/L) for each sample.

18.1 Weigh the  400 g of CaCl₂ and transfer it to a sterile 1L Nalgene media bottle.

18.2 Add a sterile stir bar and  1 L of MilliQ water and stir until dissolved.

- 19 Optional Step: weigh the freeze dried sediment to compare microplastics across various samples *
- 20 Using as sterile spatula, transfer the freeze-dried sediments to a 1 L sterile beaker.
- 21 Add  300 mL of the CaCl₂ solution and mix vigorously using the spatula.
- 22 Cover with aluminum foil and let the sample stand for  12:00:00 . 12h
- 23 Decant the CaCl₂ solution into a clean, sterile beaker and store the decanted solution at  4 °C .  
- 24 Repeat Steps 19-21 two more times with  12:00:00 and  24:00:00 of 1d 12h
sedimentation time for Step 21.
- 25 Filter all the decanted CaCl₂ solution through a 20 μm cellulose filter.
- 26 Use Steps 11 -15 (Microplastics Isolation - Water) to isolated sediment microplastics.