

Oct 31, 2024

NPC1 inhibitor treatments and immunoblotting of whole-cell lysates from cell culture systems

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.261ger657l47/v1

Felix Kraus¹, Harper JW^{2,3}

¹Department of Cell Biology, Blavatnik Institute, Harvard Medical School, 240 Longwood Ave, Boston MA 02115, USA;

²Harvard Medical School;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA



Felix Kraus

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ger657l47/v1>

Protocol Citation: Felix Kraus, Harper JW 2024. NPC1 inhibitor treatments and immunoblotting of whole-cell lysates from cell culture systems. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261ger657l47/v1>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 31, 2024

Last Modified: October 31, 2024

Protocol Integer ID: 111354

Keywords: ASAPCRN, inhibition of the lysosomal cholesterol transporter npc1, lysosomal cholesterol transporter npc1, autophagy protein, abundance of autophagy protein, cell lysate, cell lysates from cell culture system, npc1 inhibitor treatment, hela cell culture system, cell culture system, u18666a inhibitor, protein, immunoblotting

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-025160

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This is a protocol for assessing abundance of autophagy proteins after inhibition of the lysosomal cholesterol transporter NPC1 via the U18666A inhibitor by Western blotting from whole-cell lysates derived from HeLa cell culture systems.

Troubleshooting



Safety warnings

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.



Cell Culture and Treatments

5m

- 1 Seed HeLa TMEM192-3xHA cells of desired genotypes into 6-well plates.
- 2 Upon reading desired confluency (~50 - 70%, depending on the duration of inhibitor treatment), treat cells according to the following four treatments:
 1. Fed
 2. NPC1 inhibitor U18666A (2 μ M)
- 3 Treat cells for 1 to 3 days (depending on the experiment).
Change media and inhibitor daily.
- 4 Aspirate growth media, wash twice with 1xPBS and harvest cells in ice-cold PBS on ice by scraping cells from the wells.

SDS-PAGE and immunoblotting

5m

- 5 For whole-cell lysates, prepare samples following standard protocols, and final samples should be in LDS buffer with DTT or similar. Incubate samples at 80 °C for 00:05:00 .
- 6 Load samples into a NuPAGE Novex [®] 4-12% Bis-Tris Midi Protein Gels and separate by electrophoresis in 1xTris/Glycine/SDS buffer.
- 7 Transfer proteins to PVDF or nitrocellulose membranes by standard wet transfer in 20% methanol Tis/Glycine buffer.
- 8 Block membrane in blocking buffer (5% non-fat dry milk or 3% BSA in TBST) at Room temperature for 01:00:00 .
- 9 Incubate membrane in primary antibody solution (blocking solution plus primary antibody at 1:500-1:1,000, depending on the primary antibody) at 4 °C for 12:00:00 to 16:00:00 .

5m

1h

1d 4h





Primary antibodies:

NCOA4

FTH1



SQSTM1 / p62
LC3B
panGARBAP
Actin

- 10 Wash membrane six times with TBST for  00:05:00 each wash. 5m
- 11 Incubate membrane in secondary antibody solution (blocking solution plus secondary antibody conjugated to HRP at 1:5,000-1:10,000) at  Room temperature for  01:00:00 . 1h
- 12 Wash membrane four times with TBST for  00:05:00 each wash. 5m
- 13 Apply Western Lightning Plus Chemiluminescence substrate (Revvity) to membrane and acquire blot images using a ChemiDoc MP imager.
- 14 Process raw image files with Image Lab software (Bio-Rad).